



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY
APPLICATION #: OPM-0056-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: [X] New [] Renewal [] Update to Pre-CBC 2013 OPA Number:

Manufacturer Information

Manufacturer: TRUMPF Medical Systems, Inc.

Manufacturer's Technical Representative: Jeff Saunders

Mailing Address: 1046 LeGrand Boulevard, Charleston, SC 29492

Telephone: 843-534-0606 Email: Jeffrey.Saunders@us.trumpf.com


Product Information

Product Name: TRULIGHT 1000

Product Type: Wall Mounted Examination Light

Product Model Number: None

General Description: Wall Mounted Examination Light



Applicant Information

Applicant Company Name: TRUMPF Medical Systems, Inc.

Contact Person: Jeff Saunders

Mailing Address: 1046 LeGrand Boulevard, Charleston, SC 29492

Telephone: 843-534-0606 Email: Jeffrey.Saunders@us.trumpf.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

Signature of Applicant: [Signature] Date: 9/30/2013

Title: Manager, Customer Support Company Name: Trumpf Medical Systems

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

Registered Design Professional Preparing Engineering Recommendations

Company Name: CYS Structural Engineers, Inc.

Name: Dieter T. Siebald California License Number: SE #4346

Mailing Address: 2495 Natomas Park Drive, Suite 650, Sacramento, CA 95833

Telephone: 916-920-2020 Email: dieters@cyseng.com

OSHPD Special Seismic Certification Preapproval (OSP)

- Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required)
- Special Seismic Certification is no preapproved

Certification Method(s)

- Testing in accordance with: ICC-ES AC156 FM 1950-10
- Other* (Please Specify): _____


*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing.

- Analysis
- Experience Data
- Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

List of Attachments Supporting the Manufacturer's Certification

- Test Report Drawings Calculations Manufacturer's Catalog
- Other(s) (Please Specify): _____

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY

Signature:  Date: February 12, 2014

Print Name: Jeffrey Y. Kikumoto

Title: Senior Structural Engineer

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"






TRULIGHT 1000
WALL MOUNTED EQUIPMENT SUPPORT & ATTACHMENT
TABLE OF CONTENTS
OPM-0056-13

	PAGE
GENERAL NOTES	2
ABBREVIATIONS	4
DESIGN CRITERIA	5
TRULIGHT 1000	
ELEVATION	6
PLAN	7
WALL BEARING BRACKET	8
WALL SUPPORT	
STUD WALL SUPPORT	9
STUD CONNECTION AT TOP	10
TUBE STEEL POST SUPPORT	11
TUBE STEEL POST TOP & BOTTOM CONNECTIONS	12
MOUNTING PLATE ATTACHMENT	13
UNDERSIDE OF METAL DECK DETAIL	14

- NOTES:**
1. THESE DRAWINGS ARE PREPARED FOR TRUMPF MEDICAL SYSTEMS, INC., CHARLESTON, SOUTH CAROLINA.
 2. THE CONTRACTOR AND INSPECTOR OF RECORD SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OSHPD PRE-APPROVAL PROGRAMS WEBSITE.
 3. THIS PRE-APPROVAL COVERS THE SUPPORTS & ATTACHEMENTS OF THE UNIT TO THE SUPPORTING STRUCTURE. THE UNIT AND WALL BEARING HARDWARE ARE SUPPLIED BY TRUMPF. ALL THREAD RODS, MOUNTING PLATE AND REQUIRED WALL SUPPORT ELEMENTS SHALL BE SUPPLIED AND INSTALLED BY THE CONTRACTOR.



SHEET TITLE: TABLE OF CONTENTS

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 1 of 14
--	--------------------------------------	---


L:\Jobs\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:57am Login:camachom DimScale:1 LTScale:6



GENERAL NOTES:

1. THIS OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.
2. IT IS THE RESPONSIBILITY OF THE STRUCTURAL ENGINEER OF RECORD FOR A SITE SPECIFIC PROJECT TO VERIFY:
 - A. THE ADEQUACY OF THE NEW OR EXISTING STRUCTURE TO RESIST THE FORCES AND WEIGHT SPECIFIED FOR EACH EQUIPMENT IN ADDITION TO ALL OTHER LOADS. PROVIDE AND DESIGN SUPPLEMENTARY MEMBERS AS REQUIRED.
 - B. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS.
 - C. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. THE SPACING SHOWN IN THE TEST LOADS TABLE ON PAGE 3 IS THE REQUIRED MINIMUM SPACING OF THE 1/2" DIAMETER ANCHOR BOLTS. THE REQUIRED SPACING FROM ANCHORS OF OTHER DIAMETERS AND EMBEDMENTS MAY VARY AND SHALL BE EVALUATED BY THE SEOR.
 - D. THAT THE INSTALLATION IS IN CONFORMANCE WITH THE CBC 2013 AND WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL.
 - E. THAT THE ACTUAL EQUIPMENT'S WEIGHT, CENTER OF GRAVITY (CG) LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS, AND THE MATERIAL AND GAUGE OF THE EQUIPMENT WHERE ATTACHMENTS ARE MADE, AGREE WITH THE INFORMATION SHOWN ON THE PRE-APPROVAL DOCUMENTS.
- 3A. EXPANSION ANCHORS INSTALLED IN NORMAL WEIGHT OR SAND-LIGHTWEIGHT CONCRETE SHALL BE CARBON STEEL HILTI KB-TZ EXPANSION ANCHORS COMPLYING WITH ESR-1917 REISSUED MAY 1, 2013.
- B. INSTALLATION: INSTALL THE EXPANSION ANCHORS IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC EVALUATION REPORT FOR THE SPECIFIC ANCHOR AND THE PARAMETERS GIVEN IN THE TABLE ON PAGE 3.
- C. TESTING:
 - JOB TESTING: FOR VERIFYING SATISFACTORY INSTALLATION WORKMANSHIP, PERFORM JOB SITE TESTING IN ACCORDANCE WITH THE TENSION LOAD TABLE PROVIDED IN THIS DOCUMENT. TEST 50% OF THE INSTALLED ANCHORS. THE TEST LOAD MAY BE APPLIED BY ANY METHOD THAT WILL EFFECTIVELY MEASURE THE TENSION IN THE ANCHOR SUCH AS DIRECT PULL WITH A HYDRAULIC JACK OR CALIBRATED SPRING LOADING DEVICES OR CALIBRATED TORQUE WRENCH METHOD. ALL TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE INSPECTOR OF RECORD. IF ANY ANCHOR FAILS THE TEST, TEST ALL ANCHORS. THE TEST SHALL BE PERFORMED 24 HOURS OR MORE AFTER INSTALLATION. TESTING MAY BE DONE PRIOR TO EQUIPMENT INSTALLATION. ALSO REFER TO CBC 1913A.7 "FIELD TESTS FOR POST-INSTALLED ANCHORS IN CONCRETE".
 - FAILURE/ACCEPTANCE CRITERIA: THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:
 - HYDRAULIC RAM METHOD: APPLY AND HOLD TEST LOAD FOR A MINIMUM OF 15 SECONDS. THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE APPLICABLE TEST LOAD WHERE WASHERS ARE USED. FOR WEDGE TYPE ANCHORS, SUCH AS HILTI KB-TZ, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER UNDER THE NUT BECOMES LOOSE.
 - TORQUE WRENCH METHOD: THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: ONE-HALF (1/2) TURN OF THE NUT.

OPM-0056-13



SHEET TITLE: GENERAL NOTES


<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	<p>TEL (916) 920-2020 www.cyseng.com</p>	<p>Job No: 13084 Date: 2-11-2014 Page: 2 of 14</p>
--	--	--

L:\Jobs\13084 Trumph - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:57am Login:camachom DimScale:1 LScale:6



GENERAL NOTES CONTINUED:

3D. TEST VALUES: APPLY TEST LOADS TO ANCHORS WITHOUT REMOVING THE NUT IF POSSIBLE, SEE TABLE BLW.




MECHANICAL ANCHOR

ANCHOR DIA (INCH) da	INSTALLATION EMBED (INCH) hnom	EFFECTIVE EMBED (INCH) hef	HOLE DEPTH (INCH) ho	MIN CONC THICKNESS (INCH) h	MIN CONC EDGE DISTANCE (INCH)	MIN AB SPACING (INCH)	TEST LOAD		CONDITION OF ANCHORAGE
							TENSION LOAD (LBS)	TORQUE (FT-LBS)	
3/8	25/16	2	25/8	4*	9	4	-	25	CASE 1

* CONC THICKNESS MAY BE 3/4" AT CONC OVER MTL DECK. SEE TABLE 6 OF ESR-1917.


4. A SINGLE CASE OF ANCHORAGE IS SPECIFIED AND PRESENTED IN THIS PRE-APPROVAL:



DATE: 02/12/2014

CASE 1: ANCHORAGE DTLS LOCATED AT ANY LEVEL OF A BLDG ($z/h \leq 1.0$), IT IS ASSUMED THAT THE FLRS ARE BUILT OF A MIN 3/4" SAND-LWC OR NWC TOPPING OVER MTL DECK OR 4" NWC SLAB ($f'_c = 3000$ PSI, MIN).

5. THIS PRE-APPROVAL MAY BE USED AT ANY GEOGRAPHICAL LOCATION IN THE STATE OF CALIFORNIA. WHERE S_{ps} IS LESS THAN OR EQ TO 2.50.



SHEET TITLE: GENERAL NOTES (CONTINUED)



<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	Job No: 13084
	Date: 2-11-2014
	Page: 3 of 14

L:\Jobs\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:57am Login:camachom DimScale:1 LScale:6



ABBREVIATIONS:

@	AT	KSI	KIPS PER SQUARE INCH
AB	ANCHOR BOLT	LBS	POUNDS
ABV	ABOVE	LRFD	LOAD & RESISTANCE FACTOR DESIGN
ALT	ALTERNATE	LWC	LIGHT WEIGHT CONCRETE
ASD	ALLOWABLE STRENGTH DESIGN	MAX	MAXIMUM
ASTM	AMERICAN SOCIETY FOR TESTING & MATERIALS	MFR	MANUFACTURER
BLDG	BUILDING	MIN	MINIMUM
BLKG	BLKG	MTL	METAL
BLW	BELOW	NWC	NORMAL WEIGHT CONCRETE
BOTT	BOTTOM	OPM	OSPHD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION
BRG	BEARING		
BTW	BETWEEN	OSHPD	OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT
CBC	CALIFORNIA BUILDING CODE		
CG	CENTER OF GRAVITY	PG(S)	PAGE(S)
CLR	CLEAR	PL	PLATE
CLSE	CALIFORNIA LICENSED STRUCTURAL ENGINEER	PSI	POUNDS PER SQUARE INCH
CL	CENTERLINE	RDP	REGISTERED DESIGN PROFESSIONAL
CONC	CONCRETE	REQ	REQUIRED
CONN	CONNECTION	SAD	SEE ARCHITECTURAL DRAWINGS
DBL	DOUBLE	SEOR	STRUCTURAL ENGINEER OF RECORD
DTL(S)	DETAIL(S)	SIM	SIMILAR
DIA (φ)	DIAMETER	SLWC	SAND LIGHT WEIGHT CONCRETE
DWG	DRAWING	SMS	SHEET METAL SCREWS
(E)	EXISTING CONDITION	STD	STANDARD
EA	EACH	STL	STEEL
ELEV	ELEVATION	THRD	THREAD OR THREADED
EQ	EQUAL	TOC	TOP OF CONCRETE
EQUIP	EQUIPMENT	TYP	TYPICAL
ES	EACH SIDE	V	ANCHORAGE SHEAR REACTION DUE TO SEISMIC FORCE
FF	FINISH FLOOR	VERT	VERTICAL
FLR	FLOOR	W _p	COMPONENT OPERATING WEIGHT
F _p	HORIZONTAL SEISMIC FORCE PER ASCE 7-10 SEISMIC FORCE REQUIREMENTS	W	WITH
FT (')	FOOT/FEET		
F _{pv}	VERTICAL SEISMIC FORCE PER ASCE 7-10 SEISMIC FORCE REQUIREMENTS		
F _y	SPECIFIED YIELD STRENGTH OF REINFORCING, PSI OR SPECIFIED MINIMUM YIELD STRESS OF STEEL, KSI		
GA	GAUGE		
GWB	GYPSUM WALLBOARD		
HSS	HOLLOW STRUCTURAL SECTION		
HT	HEIGHT		
ICC	INTERNATIONAL CODE COUNCIL		
IN (")	INCH		



SHEET TITLE: ABBREVIATIONS

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 4 of 14
--	--------------------------------------	---

L:\Jobs\13\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6





DESIGN CRITERIA

ANCHORAGE DESIGN FOR EXAM LIGHT IS PER 2013 CBC AT LRFD LEVEL FORCES.

- $W_p = 20$ LBS
- $z/h = 1.0$ UPPER FLRS ABV THE BASE
- $S_{DS} = 2.50$
- $I_p = 1.5$
- $a_p = 2.5$
- $R_p = 2.5$
- MIN $F_p = 0.30 S_{DS} I_p W_p = 1.125 W_p$
- MAX $F_p = 1.60 S_{DS} I_p W_p = 6.000 W_p$
- $F_p = 0.4 a_p S_{DS} I_p W_p (1+2 z/h)/R_p = 4.5 W_p = 90$ LBS
- $F_{pv} = \pm 0.20 S_{DS} W_p = \pm 0.50 W_p = \pm 10$ LBS

ANCHORAGE DESIGN FOR WALL IS PER 2013 CBC AT LRFD LEVEL FORCES.

- $z/h = 1.0$
- $S_{DS} = 2.50$
- $I_p = 1.5$
- $a_p = 1.0$
- $R_p = 2.5$
- $\phi = 2.50$ FOR CONC ANCHORS ONLY
- $F_p = 0.4 a_p S_{DS} I_p W_p (1+2 z/h)/R_p = 1.8 W_p$
- $F_{pv} = \pm 0.20 S_{DS} W_p = \pm 0.50 W_p$






SHEET TITLE: DESIGN CRITERIA

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 5 of 14
--	--------------------------------------	---

L:\Jobs\13\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000






SHEET TITLE: TRULIGHT 1000 ELEVATION

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 6 of 14
--	--------------------------------------	---

L:\Jobs\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000




SHEET TITLE: TRULIGHT 1000 PLAN


 CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 7 of 14
---	--------------------------------------	---

L:\Jobs\13\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000




PLAN VIEW



SIDE VIEW

NOTE:
WALL BRG BRACKET IS SUPPLIED BY TRUMPF
AND MOUNTING FL. IS SUPPLIED BY CONTRACTOR




SHEET TITLE: TRULIGHT 1000 WALL BEARING BRACKET DETAIL

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	Job No: 13084
	Date: 2-11-2014
	Page: 8 of 14

L:\Jobs\13\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000

TRUMPF



NOTE:

THE WALL STUDS & CONNECTIONS TO SUPPORTING STRUCTURE HAVE BEEN DESIGNED FOR A MAX ALLOWABLE INTERSTORY DRIFT OF $D_{pl} = (0.015)(14' \times 12' / ') (1.5) = 3.78"$



SHEET TITLE: STUD WALL SUPPORT



CYS STRUCTURAL ENGINEERS, INC.


2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com


Job No:	13084
Date:	2-11-2014
Page:	9 of 14

L:\Jobs\13\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camacthom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000



PERPENDICULAR TO FLUTES



PARALLEL TO FLUTES




SHEET TITLE: STUD CONNECTION AT TOP

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 10 of 14
--	--------------------------------------	--

L:\Jobs\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom Dimscale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000



2- 3/8" Ø EXPANSION ANCHORS.
SEE PGS 2, 3 & 12 FOR
INSTALLATION REQUIREMENTS

8"

PL 3/8x3 1/2x0'-10"
(ASTM A36, Fy = 36 KSI MIN)

NOTE:
THE HSS & CONNECTIONS TO SUPPORTING STRUCTURE HAS
BEEN DESIGNED FOR A MAX ALLOWABLE INTERSTORY DRIFT
OF $D_{pl} = (0.015)(14' \times 12' / 12) / (1.5) = 3.78"$




SHEET TITLE: TUBE STEEL POST SUPPORT

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	Job No: 13084
	Date: 2-11-2014
	Page: 11 of 14


L:\Jobs\13\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000



(E) NWC OR SLWC (3000 PSI MIN)
OVER STL-W DECK (20 GA MIN)

OPM-0056-13
TOP CONNECTION
BY: Jeffrey Y. Kikumoto
DATE: 02/12/2014



BOTTOM CONNECTION




SHEET TITLE: TUBE STEEL POST TOP & BOTTOM CONNECTIONS

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	Job No: 13084
	Date: 2-11-2014
	Page: 12 of 14


L:\Jobs\13084 Trumph - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6

TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000




PLAN AT HSS BLKG

PLAN AT HSS POST



ELEV AT HSS BLKG

SECTION AT HSS POST



SHEET TITLE: MOUNTING PLATE ATTACHMENT

<p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13084 Date: 2-11-2014 Page: 13 of 14
--	--------------------------------------	--

L:\Jobs\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom Dimscale:1 LTScale:6



TRUMPF MEDICAL SYSTEMS
TRULIGHT 1000



ALLOWABLE VERT LOAD (LBS) T_{ASD}	ANCHOR DIA (INCH) d_a	EFFECTIVE EMBED (INCH) h_{ef}	HOLE DEPTH (INCH) h_o	MIN SPACING (INCH) S_{min}	MIN EDGE DISTANCE (INCH) C_{min}	TORQUE (FT-LBS)
610	1/2	2 1/4	2 5/8	6 3/4	7	40

NOTES:

1. MIN SPACING REQUIREMENT IS ALONG FLUTE LENGTH ONLY.
2. ALLOWABLE LOADS TAKEN FROM OSHPD OPD-0001-13.



SHEET TITLE: UNDERSIDE OF METAL DECK



CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	13084
Date:	2-11-2014
Page:	14 of 14

L:\Jobs\13084 Trumf - Trulight 1000 Wall Mounted OPM-0056-13\STRU\S1.dwg Time:Feb28,2014-09:58am Login:camachom DimScale:1 LTScale:6