



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT  
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL  
OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY  
APPLICATION #: OPM-0065-13

OSHPD Preapproval of Manufacturer's Certification (OPM)


Type:  New  Renewal  Update to Pre-CBC 2013 OPA Number: OPA-2928-10

Manufacturer Information

Manufacturer: Siemens Healthcare Diagnostics  
Manufacturer's Technical Representative: Fernando Jusino  
Mailing Address: 511 Benedict Ave., Tarrytown, NY 10591  
Telephone: 914-524-2573 Email: Fernando.jusino@siemens.com

Product Information


Product Name: Advia Centaur XP & XPT  
Product Type: Medical Analyzer  
Product Model Number: XP & XPT  
General Description: Immunoassay Diagnostic System




Applicant Information

Applicant Company Name: Siemens Healthcare Diagnostics  
Contact Person: Fernando Jusino  
Mailing Address: 511 Benedict Ave., Tarrytown, NY 10591  
Telephone: 914-524-2573 Email: Fernando.jusino@siemens.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

Signature of Applicant:  Date: 11/22/2013  
Title: Senior Product Engineer Company Name: Siemens Healthcare Diagnostics

\*Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs\*





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT  
FACILITIES DEVELOPMENT DIVISION**

**Registered Design Professional Preparing Engineering Recommendations**

Company Name: CYS Structural Engineers, Inc.

Name: Dieter T. Siebald California License Number: S4346

Mailing Address: 2495 Natomas Park Drive, Suite 650, Sacramento, CA 95833

Telephone: 916-920-2020 Email: dieters@cyseng.com

**OSHPD Special Seismic Certification Preapproval (OSP)**

- Special Seismic Certification is preapproved under (Separate application for OSP is required)
- Special Seismic Certification is no preapproved

**Certification Method(s)**

- Testing in accordance with:  ICC-ES AC156  FM 1950-10
- Other\* (Please Specify): \_\_\_\_\_


\*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing.

- Analysis
- Experience Data
- Combination of Testing, Analysis, and/or Experience Data (Please Specify): \_\_\_\_\_

**List of Attachments Supporting the Manufacturer's Certification**

- Test Report  Drawings  Calculations  Manufacturer's Catalog
- Other(s) (Please Specify): \_\_\_\_\_

**OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY**


Signature:  Date: February 12, 2014

Print Name: Jeffrey Y. Kikumoto

Title: Senior Structural Engineer

Condition of Approval (if applicable): \_\_\_\_\_



"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"



ADVIA CENTAUR XP & XPT  
TABLE OF CONTENTS  
OPM-0065-13

	PAGE
GENERAL NOTES .....	2
ABBREVIATIONS & DESIGN CRITERIA .....	5
ADVIA CENTAUR XP PLAN DIMENSIONS .....	6
ADVIA CENTAUR XP ELEV DIMENSIONS .....	7
ADVIA CENTAUR XPT PLAN DIMENSIONS .....	8
ADVIA CENTAUR XPT ELEV DIMENSIONS .....	9
ANCHORAGE BRACKET DETAIL .....	10
ANCHORAGE DETAIL TO EQUIPMENT FRAME .....	11
ANCHORAGE DETAIL TO CONC FILL OVER METAL DECK .....	12
ANCHORAGE DETAIL TO CONC SLAB .....	15

- NOTES:**
1. THESE DRAWINGS ARE PREPARED FOR SIEMENS HEALTHCARE DIAGNOSTICS, INC., TARRYTOWN, NEW YORK.
  2. THE CONTRACTOR AND INSPECTOR OF RECORD SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OSHPD PRE-APPROVAL PROGRAMS WEBSITE.
  3. THIS PRE-APPROVAL COVERS THE SUPPORTS & ATTACHMENTS OF THE UNIT TO THE SUPPORTING STRUCTURE. THE UNIT, ANCHORAGE BRACKETS & CONNECTION HARDWARE TO UNIT IS SUPPLIED & INSTALLED BY SIEMENS. THROUGH BOLTS, STRUT HARDWARE & ATTACHMENTS AT SOFFIT UNDER METAL DECK & EXPANSION BOLTS SHOWN ON PAGES 12 TO 16 SHALL BE SUPPLIED & INSTALLED BY THE CONTRACTOR.



SHEET TITLE: TABLE OF CONTENTS



**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833


TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	1 of 16

L:\Jobs\13\106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom DimScale:1 L1Scale:6

**GENERAL NOTES:**

1. THIS OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.
2. IT IS THE RESPONSIBILITY OF THE STRUCTURAL ENGINEER OF RECORD FOR A SITE SPECIFIC PROJECT TO VERIFY:
  - A. THE ADEQUACY OF THE NEW OR EXISTING STRUCTURE TO RESIST THE FORCES AND WEIGHT SPECIFIED FOR EACH EQUIPMENT IN ADDITION TO ALL OTHER LOADS. PROVIDE AND DESIGN SUPPLEMENTARY MEMBERS AS REQUIRED.
  - B. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS.
  - C. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. THE SPACING SHOWN IN THE TEST LOADS TABLE ON PAGE 3 IS THE REQUIRED MINIMUM SPACING OF THE 1/2" AND 5/8" DIAMETER ANCHOR BOLTS. THE REQUIRED SPACING FROM ANCHORS OF OTHER DIAMETERS AND EMBEDMENTS MAY VARY AND SHALL BE EVALUATED BY THE SEOR.
  - D. THAT THE INSTALLATION IS IN CONFORMANCE WITH THE CBC 2013 AND WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL.
  - E. THAT THE ACTUAL EQUIPMENT'S WEIGHT, CENTER OF GRAVITY (CG) LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS, AND THE MATERIAL AND GAGE OF THE EQUIPMENT WHERE ATTACHMENTS ARE MADE, AGREE WITH THE INFORMATION SHOWN ON THE PRE-APPROVAL DOCUMENTS.
- 3A. EXPANSION ANCHORS INSTALLED IN NORMAL WEIGHT OR SAND-LIGHTWEIGHT CONCRETE SHALL BE CARBON STEEL HILTI KB-TZ EXPANSION ANCHORS COMPLYING WITH ESR-1917 REISSUED MAY 1, 2013.
- B. INSTALLATION: INSTALL THE EXPANSION ANCHORS IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC EVALUATION REPORT FOR THE SPECIFIC ANCHOR AND THE PARAMETERS GIVEN IN THE TABLE ON PAGE 3.
- C. JOB TESTING: FOR VERIFYING SATISFACTORY INSTALLATION WORKMANSHIP, PERFORM JOB SITE TESTING IN ACCORDANCE WITH THE TENSION LOAD TABLE PROVIDED IN THIS DOCUMENT. TEST 50% OF THE INSTALLED ANCHORS. THE TEST LOAD MAY BE APPLIED BY ANY METHOD THAT WILL EFFECTIVELY MEASURE THE TENSION IN THE ANCHOR SUCH AS DIRECT PULL WITH A HYDRAULIC JACK OR CALIBRATED SPRING LOADING DEVICES OR CALIBRATED TORQUE WRENCH METHOD. ALL TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE INSPECTOR OF RECORD. IF ANY ANCHOR FAILS THE TEST, TEST ALL ANCHORS. THE TEST SHALL BE PERFORMED 24 HOURS OR MORE AFTER INSTALLATION. TESTING MAY BE DONE PRIOR TO EQUIPMENT INSTALLATION. ALSO REFER TO CBC 1913A.7 "FIELD TESTS FOR POST-INSTALLED ANCHORS IN CONCRETE".
- D. FAILURE/ACCEPTANCE CRITERIA: THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:
  - HYDRAULIC RAM METHOD: APPLY AND HOLD TEST LOAD FOR A MINIMUM OF 15 SECONDS. THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE APPLICABLE TEST LOAD WHERE WASHERS ARE USED. FOR WEDGE TYPE ANCHORS, SUCH AS HILTI KB-TZ, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER UNDER THE NUT BECOMES LOOSE
  - TORQUE WRENCH METHOD: THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: ONE-HALF (1/2) TURN OF THE NUT.




SHEET TITLE: GNERAL NOTES

<p><b>CYS STRUCTURAL ENGINEERS, INC.</b> 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13106 Date: 2/12/2014 Page: 2 of 16
--	--------------------------------------	---

L:\Jobs\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:cmachom DimScale:1 LScale:6

**GENERAL NOTES CONTINUED:**

3E. TEST VALUES: APPLY TEST LOADS TO ANCHORS WITHOUT REMOVING THE NUT IF POSSIBLE, SEE TABLE BELOW.



ANCHOR DIA (INCH) da	INSTALLATION EMBED (INCH) hnom	EFFECTIVE EMBED (INCH) hef	HOLE DEPTH (INCH) ho	MIN CONC THICKNESS (INCH) h <sub>min</sub>	MIN CONC EDGE DISTANCE (INCH)	MIN AB SPACING (INCH)	TEST LOAD		CONDITION OF ANCHORAGE
							TENSION LOAD (LBS)	TORQUE (FT-LBS)	
1/2	2 3/8	2	2 5/8	4	12	4	1338	40	CASE 3
1/2	3 5/8	3 1/4	4	6	12	4	2344	40	CASE 2A
5/8	3 9/16	3 3/8	3 3/4	6	12	4	2236	60	CASE 2B

4. BOLTS THROUGH CONCRETE ON METAL DECK:

- A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER SNUG TIGHT (THE SNUG TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS NOTED OTHERWISE.
- B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16")
- C. THROUGH BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION & TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.

5. SCREW ANCHORS TO BOTTOM OF CONCRETE FILL OVER METAL DECK:

- A. HILTI KH-EZ (ICC ESR-3027) TENSION TEST LOAD FOR CASE 1.

ANCHOR DIA (INCH) da	INSTALLATION EMBED (INCH) hnom	EFFECTIVE EMBED (INCH) hef	HOLE DEPTH (INCH) ho	MIN CONC THICKNESS (INCH) h <sub>min</sub>	MIN CONC EDGE DISTANCE (INCH)	MIN AB SPACING (INCH)	TENSION TEST LOAD (LBS)
1/4	1 5/8	1.18	2	3/4	1 1/4*	10*	400

\* SEE PAGE 13 OF 16 IN THIS OPM AND FOOTNOTE 2, TABLE 2 IN ESR-3027



SHEET TITLE: GENERAL NOTES (CONTINUED)



**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com


Job No:	13106
Date:	2/12/2014
Page:	3 of 16

L:\Jobs\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom DimScale:1 L1Scale:6



**GENERAL NOTES CONTINUED:**

6. FOUR (4) CASES OF ANCHORAGE ARE SPECIFIED AND PRESENTED IN THIS PRE-APPROVAL:





**CASE 1:** ANCHORAGE DETAILS LOCATED AT UPPER FLOORS ABOVE THE BASE OF A BUILDING ( $z/h <= 1.0$ ), IT IS ASSUMED THAT THE FLOORS ARE BUILT OF A MINIMUM 3/4" SAND-LIGHTWEIGHT CONCRETE TOPPING OVER METAL DECK ( $f'c = 3000$  PSI, MINIMUM).

**CASE 2A:** ANCHORAGE DETAILS LOCATED AT OR BELOW THE BASE OF A BUILDING ( $z/h=0$ ). THE FLOORS ARE ASSUMED TO BE BUILT OF A MINIMUM 6" NORMAL-WEIGHT CONCRETE SLAB ( $f'c = 3000$  PSI, MINIMUM).

**CASE 2B:** ANCHORAGE DETAILS LOCATED AT OR BELOW THE BASE OF THE BUILDING ( $z/h=0$ ). THE FLOORS ARE ASSUMED TO BE BUILT OF A MINIMUM 5" NORMAL-WEIGHT CONCRETE SLAB ( $f'c = 3000$  PSI, MINIMUM).

**CASE 3:** ANCHORAGE DETAILS LOCATED AT OR BELOW THE BASE OF THE BUILDING ( $z/h=0$ ). THE FLOORS ARE ASSUMED TO BE BUILT OF A MINIMUM 4" NORMAL-WEIGHT CONCRETE SLAB ( $f'c = 3000$  PSI, MINIMUM). FOR THIS CASE THE MAXIMUM  $S_{DS}$  IS LIMITED TO 1.80.

7. THIS PRE-APPROVAL MAY BE USED AT ANY GEOGRAPHICAL LOCATION IN THE STATE OF CALIFORNIA. WHERE  $S_{DS}$  IS LESS THAN OR EQUAL TO 2.50, EXCEPT FOR CASE 3 ANCHORAGE WHERE  $S_{DS}$  MUST BE LESS THAN OR EQUAL TO 1.80.



SHEET TITLE: GENERAL NOTES (CONTINUED)



**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833


TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	4 of 16

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimscaler:1 LScale:6

**ABBREVIATIONS:**

⊙	AT	IN (")	INCH
ABV	ABOVE	KSI	KIPS PER SQUARE INCH
ASTM	AMERICAN SOCIETY FOR TESTING & MATERIALS	LBS	POUNDS
BLDG	BUILDING	LRFD	LOAD AND RESISTANCE FACTOR DESIGN
BLW	BELOW	LWC	LIGHT WEIGHT CONCRETE
CBC	CALIFORNIA BUILDING CODE	MAX	MAXIMUM
CG	CENTER OF GRAVITY	MFR	MANUFACTURER
⊕	CENTERLINE	MIN	MINIMUM
CONC	CONCRETE	MTL	METAL
DBL	DOUBLE	NO. (#)	NUMBER OR POUNDS
DIA (∅)	DIAMETER	NWC	NORMAL WEIGHT CONCRETE
(E)	EXISTING CONDITION	OSHPD	OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT
ELEV	ELEVATION	PSI	POUNDS PER SQUARE INCH
EQUIP	EQUIPMENT	STL	STEEL
f'c	MINIMUM ULTIMATE COMPRESSIVE STRENGTH OF CONCRETE	Tu	ANCHORAGE TENSION REACTION DUE TO SEISMIC FORCE AT LRFD
FLR	FLOOR	THRD	THREAD OR THREADED
FT (')	FOOT/FEET	TYP	TYPICAL
Fy	SPECIFIED YIELD STRENGTH OF REINFORCING, PS OF STEEL, KSI	Vu	ANCHORAGE SHEAR REACTION DUE TO SEISMIC FORCE AT LRFD
GA	GAUGE	W/	WITH
		W/O	WITHOUT
		Wp	OPERATING WEIGHT
		WT	WEIGHT



**DESIGN CRITERIA**

ANCHORAGE DESIGN FOR FLOOR MOUNTED LABORATORY EQUIPMENT IS PER 2013 CBC AT LRFD LEVEL FORCES


$\alpha_p = 1.0$        $R_p = 2.5$        $I_p = 1.5$        $\Omega_0 = 2.5$  (CONCRETE ANCHORS)

$W_p$  AS NOTED ON EQUIPMENT DRAWING SHOWN ON PAGES 6 & 8

FOR CASE 1 – UPPER FLOORS ABOVE THE BASE,  $z/h \leq 1.0$   
 $S_{ps} = 2.50$        $F_p = 1.80 W_p$        $F_v = 0.50 W_p$

FOR CASES 2A & 2B – SLAB AT OR BELOW BASE,  $z/h = 0$   
 $S_{ps} = 2.50$        $F_p = 1.125 W_p$        $F_v = 0.50 W_p$


FOR CASE 3 – SLAB AT OR BELOW BASE,  $z/h = 0$   
 $S_{ps} = 1.80$        $F_p = 0.810 W_p$        $F_v = 0.360 W_p$



SHEET TITLE: ABBREVIATIONS & DESIGN CRITERIA

	<b>CYS STRUCTURAL ENGINEERS, INC.</b> 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833	TEL (916) 920-2020 www.cyseng.com	Job No: 13106 Date: 2/12/2014 Page: 5 of 16
--	---	--------------------------------------	---

L:\Jobs\13\106 Siemens - Advia Centaur OPM-0065-13\STRU\51.dwg Time:Feb12,2014-10:32am Login:cmachom Dimscaler:1 LScale:6



SHEET TITLE: ADVIA CENTAUR XP  
PLAN DIMENSIONS



**CYS STRUCTURAL ENGINEERS, INC.**



2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	6 of 16

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimscaler:1 L1Scale:6





SHEET TITLE: ADVIA CENTAUR XP  
ELEVATION DIMENSIONS




**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	7 of 16

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimscaler:1 LScale:6



SHEET TITLE: ADVIA CENTAUR XPT  
PLAN DIMENSIONS





**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	8 of 16



L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:cmachom Dimscaler:1 LScale:6





SHEET TITLE: ADVIA CENTAUR XPT  
ELEVATION DIMENSIONS

<p><b>CYS STRUCTURAL ENGINEERS, INC.</b> 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13106 Date: 2/12/2014 Page: 9 of 16
--	--------------------------------------	---

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimscaler:1 LScale:6




**NOTE:**  
FOR CASE 2B ANCHORAGE HOLES MUST BE REAMED OUT BY THE CONTRACTOR TO 1 1/16"  $\phi$  TO ACCOMMODATE THE 5/8"  $\phi$  ANCHORS. CONTRACTOR IS RESPONSIBLE FOR REPLACING BRACKET IF IT IS DAMAGED DURING REAMING.



SHEET TITLE: ANCHORAGE BRACKET DETAIL

 <b>CYS STRUCTURAL ENGINEERS, INC.</b> 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833	TEL (916) 920-2020 www.cyseng.com	Job No: 13106 Date: 2/12/2014 Page: 10 of 16
---	--------------------------------------	--

L:\Jobs\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimscaler:1 L1Scale:6



SHEET TITLE: ANCHORAGE DETAIL  
TO EQUIPMENT FRAME

 <p><b>CYS STRUCTURAL ENGINEERS, INC.</b> 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020	Job No: 13106
	www.cyseng.com	Date: 2/12/2014
		Page: 11 of 16


L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimscaler:1 LScale:6



MAX ANCHOR FORCES  
AT LRFD AT EA AB

	Tu	Vu
CASE 1 $z/h \leq 1.0$ W/O $\Omega_o$	914#	309#
CASE 1 $z/h \leq 1.0$ W/ $\Omega_o$	914#	773#



( $\Omega_o=2.5$ ) OVERSTRENGTH FACTOR  
MUST BE APPLIED TO SHEAR  
FORCE ONLY



STRUT  $\bar{P}$  PER PGS 13 & 14

12" MIN  
EDGE DIST

DATE: 2/12/14



SHEET TITLE: ANCHORAGE DETAIL  
TO CONCRETE FILL OVER METAL DECK





**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	12 of 16

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom DimScale:1 LScale:6



SHEET TITLE: ANCHORAGE DETAIL  
TO CONCRETE FILL OVER METAL DECK


<p><b>CYS STRUCTURAL ENGINEERS, INC.</b> 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 13106 Date: 2/12/2014 Page: 13 of 16
--	--------------------------------------	--

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:cmachom Dimscaler:1 LScale:6

SIEMENS HEALTHCARE DIAGNOSTICS  
ADVIA CENTAUR XP & XPT



EXTEND STRUT LENGTH TO  
NEXT ADJ LOW FLUTE IF AB'S  
ARE LESS THAN 2" FROM  
STRUT  $\bar{r}$  SUPPORT SCREWS  
LENGTH SHALL  
ENGAGE  
2 LOW FLUTES MIN



BOTT OF (E) MTL DECK  
(LOW FLUTES SHADED  
FOR CLARITY)

BRACKET ABV

6.00" TYP

$\frac{5}{16}$ " STRUT  $\bar{r}$ , TYP. FIELD  
VERIFY DECK ORIENTATION  
FOR COORD W/ AB LAYOUT

STRUT  $\bar{r}$  SUPPORT SCREWS:  
HILTI KH-EZ  $\frac{1}{4}$ " $\phi$ x $2\frac{5}{8}$ " CONC SCREWS  
(ICC ESR-3027) 1 EA END OF  $\bar{r}$ . USE  
18 FT-LBS MAX INSTALLATION TORQUE  
USING A CALIBRATED TORQUE WRENCH,  
TYP

AB PER  
CASE 1


1.00" MAX

EITHER SIDE OF  
FLUTE  $\bar{c}$ , TYP

**PLAN VIEW**

ANCHORS PERPENDICULAR TO FLUTES

EXTEND STRUT LENGTH TO  
NEXT ADJ LOW FLUTE IF AB'S  
LENGTH SHALL ENGAGE ARE LESS THAN 2" FROM  
2 LOW FLUTES MIN STRUT  $\bar{r}$  SUPPORT SCREWS



$\frac{5}{16}$ " STRUT  $\bar{r}$ . FIELD  
VERIFY DECK ORIENTATION  
FOR COORD W/ AB LAYOUT


BOTT OF MTL DECK  
(LOW FLUTES SHADED  
FOR CLARITY)

STRUT  $\bar{r}$  SUPPORT SCREWS:  
HILTI KH-EZ  $\frac{1}{4}$ " $\phi$ x $2\frac{5}{8}$ " CONC SCREWS  
(ICC ESR-3027) 1 EA END OF  $\bar{r}$

AB PER CASE 1

**PLAN VIEW**

ANCHORS PARALLEL TO FLUTES



SHEET TITLE: ANCHORAGE DETAIL  
TO CONCRETE FILL OVER METAL DECK



**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com


Job No:	13106
Date:	2/12/2014
Page:	14 of 16

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom DimScale:1 L1Scale:6

MAX ANCHOR FORCES  
AT LRFD AT EA AB

	Tu	Vu
CASE 2A z/h = 0	1434#	483#
CASE 2B z/h = 0	1445#	483#

INCLUDES OVERSTRENGTH  
FACTOR ( $\Omega_0$ )



SHEET TITLE: ANCHORAGE DETAIL  
TO CONCRETE SLAB



**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com


Job No:	13106
Date:	2/12/2014
Page:	15 of 16

L:\Jobs\13\106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom Dimstyle:1 L1Scale:6

MAX ANCHOR FORCES  
AT LRFD AT EA AB


	Tu	Vu
CASE 3 z/h = 0	1018#	349#

INCLUDES OVERSTRENGTH  
FACTOR ( $\Omega_o$ )



**NOTE:**

THIS ANCHORAGE DTL CAN ONLY BE USED AT  
GEOGRAPHICAL LOCATIONS WHERE  $S_{ps}$  IS LESS  
THAN OR EQUAL TO 1.80.



SHEET TITLE: ANCHORAGE DETAIL  
TO CONCRETE SLAB



**CYS STRUCTURAL ENGINEERS, INC.**

2495 NATOMAS PARK DRIVE, SUITE 650  
SACRAMENTO, CA 95833

TEL (916) 920-2020  
www.cyseng.com

Job No:	13106
Date:	2/12/2014
Page:	16 of 16

L:\Jobs\13\13106 Siemens - Advia Centaur OPM-0065-13\STRU\S1.dwg Time:Feb12,2014-10:32am Login:camachom DimScale:1 L1Scale:6