

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

OFFICE USE ONLY APPLICATION FOR OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) **APPLICATION #:** OPM-0151-13 **OSHPD Preapproval of Manufacturer's Certification (OPM)** Type: ⊠ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number: **Manufacturer Information** Manufacturer: Beckman Coulter K.K. Manufacturer's Technical Representative: Hirouki Onishi Global Business Product Development, 454-32 Higashino Nagaizumi-cho Sunto-gun, Shizuoka 411-Mailing Address: 0931 Japan OK Email: honishi@beckman.com Telephone: +81-55-999-1288 **Product Information** Product Name: AU5800 Series Analyzers Product Type: Other mechanical and electrical components Product Model Number: AU5811+ISE, AU5821+ISE, AU5831+ISE and AU5841+ISE General Description: Automated Chemistry Analyzer **Applicant Information** Applicant Company Name: EASE Co. Contact Person: Jonathan Roberson, S.E. Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709 Telephone: (909) 606-7622 Email: J.Roberson@EASECo.com I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013. Signature of Applicant: 10/31/14 Date:

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs'





EASE Co.

Title: Principal Engineer

Company Name:



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations									
Company Name: EASE Co.									
Name: Jonathan Roberson, S.E. California License Number: S4197									
Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709									
Telephone: 909-606-7667 Email: J.Roberson@EASECo.com									
OSHPD Special Seismic Certification Preapproval (OSP)									
 □ Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required) □ Special Seismic Certification is not preapproved 									
Certification Method(s)									
☐ Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-10 ☐ Other* (Please Specify):									
*Use of test criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing. Analysis Experience Data DATE: 03/12/2015 Combination of Testing, Analysis, and/or Experience Data (Please Specify):									
List of Attachments Supporting the Manufacturer's Certification									
 ☐ Test Report ☐ Drawings ☐ Calculations ☐ Manufacturer's Catalog ☐ Other(s) (Please Specify): 									
OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY									
Signature: Date: 03/12/2015 Print Name: William Staehlin Title: SSE Condition of Approval (if applicable):									
Condition of Approval (if applicable):									

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"



os

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY OSH-FD-700 (REV 1/24/13)

Page 2 of 2



5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development PREAPPROVAL OF MANUFACTURER'S CERTIFICATION OPM-0151-13

THIS PREAPPROVAL CONFORMS TO THE 2013 CALIFORNIA BUILDING CODE

MANUFACTURER: EQUIPMENT NAME:

BECKMAN COULTER

AU5800 SERIES ANALYZERS

Sheet: 1 of 17

Date: 3/10/15

GENERAL NOTES

- THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2013 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2013 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2013 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 0.90, 1.2, 1.6, 1.9 & 2.2. SEE DETAIL FOR APPLICABILITY
- 4. FORCES PER ASCE 7-10 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE SDS = 0.90, a_p = 1.0, l_p = 1.5, R_p = 1.5, $z/h \le 1$ AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_o WHERE SDS = 1.20, a_p = 1.0, l_p = 1.5, R_p = 1.5, $z/h \le 1$ AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_o WHERE SDS = 1.60, a_p = 1.0, l_p = 1.5, R_p = 1.5, R
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB ON GRADE DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION BELOW GRADE. (i.e. z/h = 0.0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2013 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.

www.EquipmentAnchorage.com

BECKMAN COULTER

DES. J. ROBERSON 11-1410 JOB NO. 3/10/15 DATE


AU5800 SERIES ANALYZERS

10. EXPANSION ANCHORS:


A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	2"	N/A	N/A	See Sheet 16 of 17	25 FT-LB	1186 lb
5/8"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	3-1/8"	12"	24"	5"	60 FT-LB	3135 lb
5/8"	Normal Weight	3000	Hilti HIT HY 200	ESR-3187	4"	12"	24"	6"	60 FT-LB	4540 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 24" AWAY MINIMUM (i.e. - CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE **EDGE DISTANCES.**
- C. TESTING OF EXPANSION ANCHORS PER 2013 CBC, 1913A.7: TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL 3 INSPECTOR AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO OSHPD
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:
 - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
 - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE
 - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
 - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
 - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL



www.EquipmentAnchorage.com

BECKMAN COULTER


DES. J. ROBERSON

11-1410

AU5800 SERIES ANALYZERS

3/10/15 DATE

JOB NO.



*NOTE: 1-UNIT EQUALS 1 ANALYZER. 2-UNIT EQUALS 2 ANALYZERS. ETC... ALL CONFIG INCLUDE ONE ISE AND ONE SAMPLER UNIT IN ADDITION TO THE ANALYZERS Tu AND Vu VALUES INCLUDE Ω₀ FACTOR

NOTES:

1. FORCES ARE DETERMINED PER 2013 CALIFORNIA BUILDING CODE AND ASCE 7-10 STRENGTH DESIGN IS USED. (ap = 1.0, lp = 1.5, Rp = 1.5, Ω_0 = 1.5, z/h = 0)

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN, THIS PREAPPROVAL ENCOMPASS ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE. SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.

SEE GENERAL NOTES: SHEETS 1 AND 2



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

OF

BECKMAN COULTER

MODEL AU5811 + ISE


DES. J. ROBERSON

JOB NO. 11-1410

DATE 3/10/15

4

17 SHEETS





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER


MODEL AU5821 + ISE

DES. J. ROBERSON
JOB NO. 11-1410

DATE 3/10/15

5

of 17 sheets





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER


MODEL AU5831 + ISE

DES. J. ROBERSON
JOB NO. 11-1410

DATE 3/10/15

6

of 17 sheets



www.EquipmentAnchorage.com


BECKMAN COULTER

MODEL AU5831 + ISE

DES. J. ROBERSON

11-1410 JOB NO.

3/10/15 DATE





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

MODEL AU5841 + ISE


DES. J. ROBERSON

JOB NO. 11-1410

DATE 3/10/15

8

OF 17 SHEETS



www.EquipmentAnchorage.com

BECKMAN COULTER


MODEL AU5841 + ISE

DES. J. ROBERSON

11-1410 JOB NO.

3/10/15 DATE

SHEETS



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

DES. J. ROBERSON

11-1410


10

AU5800 SERIES ANALYZERS

DATE 3/10/15

JOB NO.

17 SHEETS



*NOTE: 1-UNIT EQUALS 1 ANALYZER, 2-UNIT EQUALS 2 ANALYZERS, ETC.. ALL CONFIGINCLUDE ONE ISE AND ONE SAMPLER UNIT IN ADDITION TO THE ANALYZERS Tu AND Vu VALUES DO NOT INCLUDE Ω_0 FACTOR

4

2004

1606

2497

2375

120

0.90

NOTES:

1. FORCES ARE DETERMINED PER 2013 CALIFORNIA BUILDING CODE AND ASCE 7-10 STRENGTH DESIGN IS USED. (ap = 1.0, lp = 1.5, Rp = 1.5, Ω_0 = 1.5, z/h \leq 1)

3-UNIT

4-UNIT

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASS ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

AU5831 + ISE

AU5841 + ISE

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.

4. SEE GENERAL NOTES: SHEETS 1 AND 2



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

MODEL AU5811 + ISE

DES. J. ROBERSON

JOB NO. 11-1410

DATE 3/10/15


SHEET 1

17 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

MAX Sps < 2.20

CONCRETE SLAB ON METAL DECK





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

MÓDEL AU5821 + ISE

DES. J. ROBERSON

3/10/15

JOB NO. 11-1410

DATE

12


SHEET


= 17 _{SHEETS}

SEISMIC SUPPORTS & ATTACHMENTS

MAX Sps ≤ 1.60

CONCRETE SLAB ON METAL DECK





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

MÓDEL AU5831 + ISE

DES. J. ROBERSON


JOB NO. 11-1410

DATE

3/10/15

13

17 SHEETS



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

MÓDEL AU5841 + ISE


DES. J. ROBERSON

JOB NO. 11-1410

DATE 3/10/15

14

17 SHEETS



www.EquipmentAnchorage.com

BECKMAN COULTER

AU5800 SERIES ANALYZERS

DES. J. ROBERSON


11-1410 JOB NO.

3/10/15 DATE

SHEET


SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK



ISE TO ANALYZER CONNECTION ISE TO SAMPLER CONNECTION





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

DES. J. ROBERSON

3/10/15

JOB NO. 11-1410

DATE

16

of 17 sheets

AU5800 SERIES ANALYZERS

SEISMIC SUPPORTS & ATTACHMENTS CONCRETE DETAILS 40" MIN EDGE DISTANCE JONCRETE EDGE N.W. OR SAND L.W. CONC. (3000 PSI MIN.) Tυ FLUTE EQUIP 12" TYP . ⊿. I" MAX OFFSET Vu STRUT MIN. 20 GA (TYP) STEEL W-DECK 4.5" FLUTE I" MIN USE 3/8" HILTI KB-TZ TYP MIN Ý. EXPANSION ANCHORS (MIN. EMBED. (hef) = 2") FLUTE DETAIL HEX NUT TOP & BOT OF FLANGE

(2 ANCHORS MIN PER STRUT)

L3 X 3 X I/4" X I'-5" MIN (A36) AT EACH ANCHOR OF STRUT, PROVIDE TAPPED HOLE BY: WTHROUGH STRUT FLANGE.

OPM (TYP) AT CONDITIONS WHERE NUT

CANNOT BE PROVIDED AT TOP SIDE

MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

BECKMAN COULTER

BRACKET DETAILS

DES. J. ROBERSON

JOB NO. 11-1410

DATE 3/10/15

17

17 SHEETS

