



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY
APPLICATION #: OPM-0152-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: [X] New [] Renewal [] Update to Pre-CBC 2013 OPA Number:

Manufacturer Information

Manufacturer: SECHRIST INDUSTRIES, INC.

Manufacturer's Technical Representative: Deepak Talati

Mailing Address: 4225 E. La Palma Avenue, Anaheim, CA 92807

Telephone: 714-579-8362 Email: DTalati@SechristUSA.com

Product Information

Product Name: Sechrist Monoplace Hyperbaric Chambers

Product Type: Other Mechanical or Electrical Component

Product Model Number: 3200, 3300H(S), 3600H(S), & 4100H(S)

General Description: Floor mounted Hyperbaric chamber with Brackets

Applicant Information

Applicant Company Name: Sechrist Industries, Inc

Contact Person: Deepak Talati

Mailing Address: 4225 E. La Palma Avenue, Anaheim, CA 92807


Telephone: 714-579-8362 Email: DTalati@SechristUSA.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

Signature of Applicant: Deepak Talati Date: 10/29/2014

Title: VP Operations And R&D Company Name: Sechrist Industries, Inc

Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dvnamic Needs





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

Registered Design Professional Preparing Engineering Recommendations

Company Name: Makeitright Inc.

Name: Joseph La Brie California License Number: S3566

Mailing Address: 55 E. Huntington Drive Ste 277, Arcadia, CA 91006

Telephone: (626) 445-0366 Email: labrie@makeitright.net

OSHPD Special Seismic Certification Preapproval (OSP)

- Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required)
- Special Seismic Certification is not preapproved

Certification Method(s)

- Testing in accordance with: ICC-ES AC156 FM 1950-10
- Other* (Please Specify): Hyperbaric Chamber does not require Special Seismic Certification. Anchorage is for anchorage to concrete slabs and decks.

*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing.

- Analysis
- Experience Data
- Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

List of Attachments Supporting the Manufacturer's Certification

- Test Report
- Drawings
- Calculations
- Manufacturer's Catalog
- Other(s) (Please Specify): _____

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY


Signature: *William Staehlin* Date: 02/19/2015

Print Name: William Staehlin

Title: SSE

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





HYPERBARIC CHAMBER

OPM-0152-13

MODEL Numbers:

3200, 3300H(S), 3600H(S) & 4100H(S)


DATE: 02/19/2015


OPM-0152-13

TABLE OF CONTENTS

SECHRIST INDUSTRIES INC. – HYPERBARIC CHAMBER
MODEL NOS. 3200, 3300H(S), 3600H(S), 4100H(S)

TC-1	TABLE OF CONTENTS
GN-1	GENERAL NOTES
GN-2	GENERAL NOTES
1.0	HYPERBARIC CHAMBER DIMENSION
2.0	SEISMIC TIEDOWN BRACKET & SPACER DETAIL
2.1	SEISMIC TIEDOWN BASE DETAIL
3.0	HYPERBARIC CHAMBER ELEVATION – SLAB ON GRADE
3.1	BRACKET ANCHORAGE DETAIL – SLAB ON GRADE
4.0	HYPERBARIC CHAMBER ELEVATION – UPPER LEVEL
4.1	BRACKET ANCHORAGE DETAIL – UPPER LEVEL
5.0	MIN. CONC. OVER METAL DECK & ANCHORAGE DISTANCE MINIMUM REQUIREMENTS
5.1	MIN. CONC. OVER METAL DECK & ANCHORAGE DISTANCE MINIMUM REQUIREMENTS



SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		 <small>hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947</small>
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: TABLE OF CONTENTS	DATE 02/16/15	PAGE TC-1		
02/19/2015		OPM-00152-13: Reviewed for Code Compliance by William Staehlin		Page 4 of 15



GENERAL NOTES

- 1.) THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2013 CALIFORNIA BUILDING CODE (CBC). THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2013 CALIFORNIA BUILDING CODE (CBC).
- 2.) SITE VERIFICATION IS REQUIRED. THE CONTRACTOR SHALL VERIFY ALL EXISTING CONDITIONS IN THE AREA(S) OF WORK PRIOR TO THE BEGINNING OF THE PROJECT. DO NOT SCALE THE DRAWINGS; ALL DIMENSIONS MUST BE VERIFIED IN THE FIELD. THE STRUCTURAL ENGINEER SHALL BE NOTIFIED, IMMEDIATELY, IF ANY DISCREPANCIES ARE FOUND.
- 3.) DESIGN CRITERIA:
 - a.) SLAB ON GRADE (SUPPORTS & ATTACHMENTS FOR MODEL NOS. 3200, 3300H(S), 3600H(S) & 4100H(S))
 $Sds = 2.2 \text{ MAX.}; l = 1.5 ; z/h = 0$ (GROUND OR SLAB ON GRADE)
UPPER FLOOR (SUPPORTS & ATTACHMENTS FOR MODEL NOS. 3200, 3300H(S), & 3600H(S))
 $Sds = 2.5 \text{ MAX.}; l = 1.5 ; z/h \leq 1$ (UPPER FLOOR)
UPPER FLOOR (SUPPORT & ATTACHMENTS FOR MODEL NO. 4100H(S))
 $Sds = 1.8 \text{ MAX.}; l = 1.5 ; z/h \leq 1$ (UPPER FLOOR)
 - b.) PER ASCE 7-10 INCLUDING SUPPLEMENT NO. 1 & TABLE 13.6-1 :
 $a_p = 1 ; R_p = 1.5 ; \Omega_o = 1.5$ (APPLY Ω_o FACTOR FOR ANCHORAGE TO CONCRETE)
- 4.) CENTER OF GRAVITY (C.G.) WEIGHT IS A MAXIMUM.
- 5.) STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL :
 - a.) CHECK THAT THE EXISTING/NEW STRUCTURE WILL BE ABLE TO SUPPORT THE MAXIMUM WEIGHTS/FORCES SHOWN IN ADDITION TO ANY OTHER LOADS TO THE STRUCTURE. PROVIDE STRENGTHENING OF STRUCTURE AS REQUIRED.
 - b.) CHECK THAT THE FLOOR OR DECK ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. ALL MINIMUM EDGE DISTANCE AND SPACING SHALL BE IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC REPORT FOR THE SPECIFIC ANCHORS USED ON THIS OPM. (SEE TABLE ON NEXT PAGE FOR ANCHOR MINIMUM SPACING & EDGE DISTANCE REQUIREMENTS)
 - c.) CHECK THAT THE INSTALLATION, SUPPORT AND ATTACHMENTS OF THE UNIT COMPLIES WITH THE 2013 CALIFORNIA BUILDING CODE AND WITH THE DETAILS SHOWN IN THIS PREAPPROVAL.
 - d.) VERIFY THAT THE ACTUAL EQUIPMENT'S MODEL NUMBER, OPERATING WEIGHT, CENTER OF GRAVITY (C.G.) LOCATION, ANCHOR LAYOUT, MATERIAL & ASTM GRADE OF THE EQUIPMENT IS THE SAME SHOWN ON THIS OPM PREAPPROVAL.
 - e.) VERIFY THAT THE PROJECT SPECIFIC Sds AND z/h VALUES RESULT IN SEISMIC FORCES DOES NOT EXCEED THE VALUES SHOWN ON THIS OPM.
- 6.) EXPANSION ANCHORS SHALL BE HILTI KB-TZ ANCHORS (PER ICC-ESR- 1917 DATED AUGUST 2014) AND INSTALLED IN NORMAL WEIGHT CONCRETE & SAND LIGHT WEIGHT CONCRETE (AT UNDERSIDE OF DECK) WHERE OCCURS. CARBON STEEL FOR INDOOR APPLICATIONS. MINIMUM EMBEDMENT OF ALL BOLTS AND TEST LOADS (UNLESS NOTED OTHERWISE ON DETAIL) IS SHOWN ON THE NEXT PAGE:



SECHRIST - HYPERBARIC CHAMBERS			
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)		JOSEPH L. LA BRIE Structural Engineer No. SE 3566	
Title: GENERAL NOTES	DATE 02/16/15	PAGE GN-1	 <small>hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947</small>
02/19/2015	OPM-00152-13: Reviewed for Code Compliance by William Staehlin		Page 5 of 15

TEST LOAD FOR HILTI KB-TZ EXPANSION ANCHORS

(f'c = 3000 PSI MIN.) – INSTALLED TO N.W. CONCRETE SLAB & SAND L.W. CONCRETE AS OCCURRING

BOLT SIZE	APPLICATION	MINIMUM EFFECTIVE EMBEDMENT (INS.)	MINIMUM ANCHOR SPACING (INS.)	MINIMUM EDGE DISTANCE (INS.)	DIRECT TENSION (LBS.)	INSTALLATION TORQUE (LBS.-IN)
5/8"	TOP OF SLAB/DECK	4"	3"	14"	2,837	720
1/2"	UNDERSIDE OF DECK	2"	6"	6"	1,190	480



- a.) WHEN INSTALLING DRILLED-IN ANCHORS IN EXISTING NON-PRESTRESSED REINFORCED CONCRETE, USE CARE AND CAUTION TO AVOID CUTTING OR DAMAGING THE THE EXISTING REINFORCING BARS. MAINTAIN A MINIMUM CLEARANCE OF ONE INCH BETWEEN THE REINFORCEMENT AND THE DRILLED-IN ANCHOR.
- b.) AFTER A MIN. OF (24) HRS. HAVE ELAPSED, ALL POST INSTALLED ANCHORS (LOADED IN EITHER PULL OUT OR SHEAR) SHALL BE TORQUED OR TENSION TESTED. WHEN POST-INSTALLED ANCHORS ARE USED FOR NON STRUCTURAL APPLICATIONS SUCH AS EQUIPMENT ANCHORAGE, 50% OR ALTERNATE BOLTS IN A GROUP, INCLUDING AT LEAST ONE-HALF THE ANCHORS IN EACH GROUP SHALL BE TESTED. IF THERE ARE ANY FAILURES, THE IMMEDIATELY ADJACENT BOLTS MUST THEN ALSO BE TESTED.
- c.) THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:

HYDRAULIC RAM METHOD: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE APPLICABLE TEST LOAD. FOR WEDGE AND SLEEVE TYPE ANCHORS, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER UNDER THE NUT BECOMES LOOSE. FOR ADHESIVE ANCHORS, WHERE OTHER THAN BOND IS BEING TESTED, THE DEVICE SHALL NOT RESTRICT THE CONCRETE SHEAR CONE TYPE FAILURE MECHANISM FROM OCCURRING


TORQUE WRENCH METHOD: THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS:
 WEDGE OR SLEEVE TYPE: ONE-HALF (1/2) TURN OF THE NUT ONE-QUARTER (1/4) TURN OF THE NUT FOR 3/8" SLEEVE ANCHOR ONLY.

- d.) OWNER'S REPRESENTATIVE IS RESPONSIBLE FOR ALL ANCHOR TESTING.
- e.) ALL TESTS SHALL BE PERFORMED IN THE PRESENCE OF THE INSPECTOR OF RECORD AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO THE ENFORCEMENT AGENCY.


7.) BOLTS THROUGH CONCRETE ON METAL DECK:

- A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
- B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE
- C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST INSTALLED ANCHORS.

SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91706 PHN: (626) 445-0366 FAX: (626) 445-2947
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: GENERAL NOTES	DATE 02/16/15	PAGE GN-2		
02/19/2015		OPM-00152-13: Reviewed for Code Compliance by William Staehlin		Page 6 of 15



EQUIPMENT SIDE ELEVATION



EQUIPMENT FRONT ELEVATION



OPM-0152-13

BY: William Staehlin


MODEL NO.	EQUIP. WEIGHT (LBS.)	MAX. PATIENT WEIGHT (LBS.)	A (INCHES.)	B (INCHES.)	C (INCHES.)	D (INCHES.)	E (INCHES.)	F (INCHES.)	G (INCHES.)
3200	1,982	500	35"	21.5"	43"	38"	47"	105.5"	68.75"
3300H(S)	2,030	500	36.6"	22.25"	44.5"	39.5"	48"	106"	66.25"
3600H(S)	2,220	700	38.8"	22.25"	44.5"	39.5"	48"	106"	66.25"
4100H(S)	3,500	700	43.25"	23.38"	46.75"	41.75"	46.5"	106"	64"



SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947	
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)					
Title: HYPERBARIC CHAMBER DIMENSIONS		DATE	02/16/15	PAGE	1.0
02/19/2015		OPM-00152-13: Reviewed for Code Compliance by William Staehlin			Page 7 of 15




OPM-0152-13
SEISMIC TIEDOWN BRACKET
 $\angle 6 \times 2 \frac{1}{2} \times \frac{7}{8}$
 (A36; $F_y = 36$ ksi & $F_u = 58$ ksi)
 DATE: 02/19/2015




SPACER
 $1 \frac{1}{2} \times 2 \times \frac{1}{8}$ (ASTM A108)


SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: SEISMIC TIEDOWN BRACKET & SPACER DETAIL	DATE	02/16/15	PAGE	2.0
	<small>hospital building design professionals 555 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947</small>			



TOP VIEW



ISOMETRIC VIEW



FRONT VIEW

SIDE VIEW

SEISMIC TIEDOWN BASE

∠6 x 6 x 3/4"

(A36; Fy = 36 ksi & Fu = 58 ksi)





SECHRIST - HYPERBARIC CHAMBERS			
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)		JOSEPH L. LABRIE Structural Engineer No. SE 3566	
Title: SEISMIC TIEDOWN BASE DETAIL		DATE 02/16/15	PAGE 2.1
02/19/2015		OPM-00152-13: Reviewed for Code Compliance by William Staehlin	
		Page 9 of 15	



SLAB ON GRADE CONDITION

[MODEL Nos. 3200, 3300H(S), 3600H(S), & 4100H(S) : (z/h = 0); Sds = 2.2 MAX.])




SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566	
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)		Scale: N.T.S.	 hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947
Title: HYPERBARIC CHAMBER ELEVATION-SLAB ON GRADE	DATE 02/16/15	PAGE 3.0	

SLAB ON GRADE CONDITION


MODEL NOS. 3200, 3300H(S), 3600H(S), & 4100H(S) : (z/h = 0); Sds = 2.2 MAX

LEGEND:


- (A) SEISMIC TIEDOWN BRACKET L6x 2 1/2 x 7/8" (A36, Fy = 36 ksi, Fu = 58 ksi)
- (B) SEISMIC TIEDOWN BASE L6x 6 x 3/4" (A36, Fy = 36 ksi, Fu = 58 ksi)
- (C) 3/8" FLANGED HEX BOLT GRADE 8 ZINC PLATED MEDIUM STRENGTH STEEL (Fu = 150,000 PSI)
- (D) (2) 5/8" HILTI KB-TZ EXPANSION ANCHORS W/ 4" EFFECTIVE EMBEDMENT (PER ICC-ES ESR-1917)
- (E) 1 1/2" x 2" x 1/8" ASTM A108 ZINC PLATED STEEL SPACER (W/ 7/16" HOLE SIZE), AS NEEDED FOR A LEVEL BRACKET INSTALLATION (SEE FRONT VIEW)
- (F) ASTM A653 11 GA. GALVANIZED SHEET PANEL (Fy = 38 KSI, Fu = 50 KSI) W/ L3 1/2 x 1 1/2 x 1/4" THICK DOUBLER PLATE STIFFENER (A36, Fy = 36 ksi, Fu = 58 ksi) SEE FRONT VIEW




**MODEL 3300H(S) & 4100H(S)
BRACKET SIDE VIEW**



**MODEL 3600H(S)
BRACKET SIDE VIEW**



**MODEL 3200
BRACKET SIDE VIEW**



**TYPICAL
BRACKET FRONT VIEW**

NOTE:
1. WELDING ELECTRODES : E70XX



SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566	
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)			
Title: BRACKET ANCHORAGE DETAIL - SLAB ON GRADE	DATE 02/16/15	PAGE 3.1	 hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947


UPPER LEVEL CONDITION

[MODEL nos. 3200, 3300H(S), & 3600H(S) : (z/h <=1); Sds = 2.5 MAX]



[MODEL no. 4100H(S) : (z/h <=1); Sds = 1.8 MAX]

SEISMIC TIEDOWN BRACKET
AND BASE ANCHORED TO
UPPER LEVEL DECK OR SLAB
(@ 4 PLACES)
SEE SHEET 4.1 FOR DETAILS

NW CONCRETE FLAT SLAB (f'c = 3,000 PSI)
OR 3 1/4" MIN. NWC OR SAND LWC CONC.
(f'c=3,000 PSI) OVER 20 GA METAL DECK
(SEE SHEET 5.0 & 5.1 FOR MIN. DECK REQ'D)



EQUIPMENT SIDE ELEVATION



EQUIPMENT FRONT ELEVATION

SEISMIC TIEDOWN BRACKET
AND BASE ANCHORED TO
UPPER LEVEL DECK OR SLAB
(@ 4 PLACES)
SEE SHEET 4.1 FOR DETAILS


L2 1/2 X 2 1/2 X 1/2 ANGLE
SUPPORT (SEE SHEET 4.1)



SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: HYPERBARIC CHAMBER ELEVATION - UPPER LEVEL		DATE 02/16/15	PAGE 4.0	hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947
02/19/2015	OPM-00152-13: Reviewed for Code Compliance by William Staehlin		Page 12 of 15	

UPPER LEVEL CONDITION
MODEL 3200, 3300H(S), & 3600H(S) : (z/h <=1); Sds = 2.5 MAX.
MODEL 4100H(S) : (z/h <=1); Sds = 1.8 MAX.

- LEGEND:**
- (A) SEISMIC TIEDOWN BRACKET L6x 2 1/2 x 7/8" (A36, Fy = 36 ksi, Fu = 58 ksi)
 - (B) SEISMIC TIEDOWN BASE L6x 6 x 3/4" (A36, Fy = 36 ksi, Fu = 58 ksi)
 - (C) 3/8" FLANGED HEX BOLT GRADE 8 ZINC PLATED MEDIUM STRENGTH STEEL (Fu = 150,000 PSI)
 - (D) (2) 5/8" A36 ALL-THREAD ROD (A36, Fy= 36 ksi, Fu = 58 ksi) SEE SHEET 5.0 & 5.1 FOR BALANCE OF INFORMATION
 - (E) ANGLE SUPPORT L2 1/2 x 2 1/2 x 1/4" (A36, Fy = 36 ksi, Fu = 58 ksi) TO SPAN 2 LOWER FLUTES MIN.
 - (F) 1 1/2" x 2" x 1/8" ASTM A108 ZINC PLATED STEEL SPACER (W/ 7/16" HOLE SIZE), AS NEEDED FOR A LEVEL BRACKET INSTALLATION (SEE FRONT VIEW)
 - (G) ASTM A653 11 GA. GALVANIZED SHEET PANEL (Fy= 38 KSI, Fu= 50 KSI) W/ L3 1/2 x 1 1/2 x 1/4" THICK DOUBLER PLATE STIFFENER (A36, Fy=36 ksi, Fu=58 ksi) SEE FRONT VIEW



MODEL 3300H(S) & 4100H(S)
BRACKET SIDE VIEW

MODEL 3600H(S)
BRACKET SIDE VIEW

MODEL 3200
BRACKET SIDE VIEW

(2) 1/2" HILTI KB-TZ W/ 2" EFFECTIVE EMBED TO UNDERSIDE OF DECK (TYP.)
 SEE SHEET 5.0 & 5.1 FOR BALANCE OF INFORMATION

NOTE:
 1. WELDING ELECTRODES : E70XX

3/8" BOLT TO UNIT FRAME (@ 4 PLACES)

ASTM A653 11 GA. GALVANIZED SHEET FRAME (Fy = 38 KSI & Fu = 50 KSI) & L3 1/2 x 1 1/2 x 1/4" DOUBLER PLATE STIFFENER W/ 7/16" HOLE

1 1/2" x 2" x 1/8" STEEL SPACER AS NEEDED FOR A LEVEL INSTALLATION


HYPERBARIC CHAMBER

NW CONCRETE FLAT SLAB (f'c = 3,000 PSI) OR 3 1/4" MIN. NWC OR SAND LWC CONC. (f'c=3,000 PSI) OVER 20 GA METAL DECK (SEE SHEET 5.0 & 5.1 FOR MIN. DECK REQ'D)

TYPICAL
BRACKET FRONT VIEW




SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: BRACKET ANCHORAGE DETAIL - UPPER LEVEL	DATE 02/16/15	PAGE 4.1	hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91706 PHN: (626) 445-0366 FAX: (626) 445-2947	




HILTI KB-TZ (UNDERSIDE OF TYP. "W" DECK ANCHORAGE)

NOTE:

1. HEX NUT TOP & BOTTOM OF ANGLE (TYP.); WHERE NUT CANNOT BE PROVIDED AT TOP SIDE OF ANGLE, PROVIDE TAPPED HOLE THROUGH ANGLE OR DOUBLE NUTS BELOW ANGLE W/ TOP NUT WELDED TO UNDERSIDE OF ANGLE.


SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: CONCRETE OVER METAL DECK REQUIREMENTS		DATE 02/16/15	PAGE 5.0	hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947



HILTI KB-TZ (UNDERSIDE OF "B" DECK ANCHORAGE)

NOTES:

1. HEX NUT TOP & BOTTOM OF ANGLE (TYP.); WHERE NUT CANNOT BE PROVIDED AT TOP SIDE OF ANGLE, PROVIDE TAPPED HOLE THROUGH ANGLE OR DOUBLE NUTS BELOW ANGLE W/ TOP NUT WELDED TO UNDERSIDE OF ANGLE.

SECHRIST - HYPERBARIC CHAMBERS		JOSEPH L. LA BRIE Structural Engineer No. SE 3566		
Model nos.: 3200, 3300H(S), 3600H(S), 4100H(S)				
Title: CONCRETE OVER METAL DECK REQUIREMENTS		DATE 02/16/15	PAGE 5.1	hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947