

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT **FACILITIES DEVELOPMENT DIVISION**

OFFICE USE ONLY PPLICATION FOR OSHPD PREAPPROVAL OF **MANUFACTURER'S CERTIFICATION (OPM) APPLICATION #:** OPM-0256-13 **OSHPD Preapproval of Manufacturer's Certification (OPM)** Type: ☐ New Renewal Update to Pre-CBC 2013 OPA Number: **Manufacturer Information** Manufacturer: Skytron LLC. Manufacturer's Technical Representative: Bob Vreeland Mailing Address: Skytron US, 5085 Corporate Exchange Blvd SE, Grand Rapids, MI 49512 Telephone: (616) 656-1189 Email: byreeland@skytron.us **Product Information** Product Name: Skytron Ergon 3 Hospital Ceiling Mounted Booms 10256-13 Product Type: 100 Series, 110 Series, 120 Series, 200 Series, Product Model Number: 300 Series, 310 Series, 320 Series, 330 Series, Stellar XL Surgical Lights General Description: Skytron's Ergon 3 ceiling mounted booms are mounted from a single point of attachment. Multiple booms at different elevations are used for equipment, monitors and lights **Applicant Information** Applicant Company Name: ISAT Seismic Bracing

Contact Person: William V Joerger

Mailing Address: 1020 Crews Road, Suite Q, Matthews NC 28105

Telephone: 510-714-0216 Email: wvjoerger@isatsb.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

WVJaces Signature of Applicant: Date: 01/22/2016

Title: Principal Structural Engineer Company Name: ISAT Seismic Bracing





"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs

STATE OF CALIFORNIA - HEALTH AND HUMAN SERVICES AGENCY OSH-FD-700 (REV 5/30/13)



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations								
Company ISAT Seismic Bracing Name:								
Name: William V Joerger California License Number: S4545								
Mailing Address: 1020 Crews Rd, Matthews NC 28105								
Telephone: 510-714-0216 Email: wvjoerger@isatsb.com								
OSHPD Special Seismic Certification Preapproval (OSP)								
 □ Special Seismic Certification is preapproved under OSP-(Separate application for OSP is required) □ Special Seismic Certification is no preapproved 								
Certification Method(s)								
☐ Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-10 ☐ Other* (Please Specify):								
Equipment is considered to be rugged. OPM is for anchorage to concrete slabs.								
*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing. BY: Jeffrey Enzler								
✓ Analysis								
Experience Data DATE: 12/19/2016								
☐ Combination of Testing, Analysis, and/or Experience Data (Please Specify):								
List of Attachments Supporting the Manufacturer's Certification								
 ☐ Test Report ☐ Other(s) (Please Specify): Calculations ☐ Manufacturer's Catalog								
OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY								
Signature: Date: 12-19-2016								
Print Name: Jeffrey Enzler								
Title: DSE								
Condition of Approval (if applicable):								

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

Page 2 of 2



A Division of Tomarco Contractor Specialties

Submittal Documents

OSHPD OPM-0256-13

SUPPORT AND ATTACHMENT OPM CONSTRUCTION DRAWINGS SKYTRON ERGON 3 BOOMS WITH STELLAR XL SURGICAL LIGHTS

SKYTRON, LLC

BUILDING

ISAT 1020 Crews Road Suite Q Matthews, NC 28105 704-841-4080

"Empowered by Experience"

REV 6

OSHPD OPM-0256-13 DWG - i

WUJ 15 NOVIG

FILE NO.: CLT-0515-070



OSHPD OPM-0256-13

DRAWING INDEX

INDEX

Cover Page	рi
ndex Page	p ii
General Notes Attachment Notes	p 1
	p 2
Skytron Ergon 3 at Metal Deck Slabs OPM-0256-13 Skytron Ergon 3 at Solid Concrete Slabs	р3
Skytron Ergon 3 at Solid Concrete Slabs	p 4
Skytron Ergon 3 Miscellaneous Steel Page 1 effrey Enzler	p 5
Skytron Ergon 3 Miscellaneous Steel Page 2: 12/19/2016	p 6
Skytron Ergon 3 Miscellaneous Steel Page 3	p 7
Skytron Ergon 3 Top Plate Details	p 8
Skytron Ergon 3 Mounting Plate Details	p 9
Skytron Ergon 3 Forces and Moments	p 10
Skytron Ergon 3 Stellar XL Forces and Moments	p 11
Skytron Equipment Mounting and Anchorage Forces	p 12

OSHPD OPM-0256-13

MANUFACTURE: SKYTRON
EQUIPMENT TYPE: CEILING MOUNTED BOOM
AND SURGICAL LIGHTS (NON-MOTORIZED)

GENERAL NOTES:

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.
- 2. SEISMIC CRITERIA USED: $S_{DS} = 2.5$ $I_P = 1.5$ ap = 1.0 Rp = 1.5 $z/h \le 1.0$ FpH = 3.00 AND FpV = 0.50.
- 3. SUPPORT AND ATTACHMENT FORCES ARE DETERMINED USING ASCE 7-10 CHAPTER 13 "SEISMIC DESIGN REQUIREMENTS FOR NONSTRUCTURAL COMPONENTS". AN OVERSTRENGTH FACTOR Ω_0 = 1.5 IS USED FOR CONCRETE MATERIALS PER ASCE 7-10 SUPPLEMENT 1 TABLE 13.6-1. LOADS SHOWN ARE STRENGTH DESIGN LOADS PER CBC 2013 SECTION 1909A.
- 4. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 5. STEEL MATERIALS: ANGLES AND PLATE ASTM A36, ROUND HSS ASTM A500 GRADE B (FY = 42 KSI), RECTANGULAR HSS ASTM A500 GRADE B (FY = 46 KSI), WELDING ELECTROIDS E70XX. ALL THREAD ROD ASTM A36, NUTS ASTM A563, WASHERS ASTM F436. BOLTS SUPPLIED BY SKYTRON ARE ASTM A193 B7.
- 6. CONCRETE SLABS:
 - a. FOR ELEVATED SOLID CONCRETE SLABS: 6" THICKNESS OF NORMAL WEIGHT CONCRETE WITH 4000 PSI MINIMUM STRENGTH.
 - b. METAL DECK: 3" DEEP COMPOSITE STEEL DECK, 20 GAGE MINIMUM, 4 1/2 INCH MINIMUM BOTTOM FLUTE WIDTH AND FLUTE SPACING IS 12", WITH 3 1/4 INCH SAND LIGHT WEIGHT CONCRETE COVER AT 4000 PSI MINIMUM STRENGTH.
- 7. POST-INSTALLED CONCRETE ANCHORS: HILTI KWIK BOLT TZ (ESR-1917) 5/8" DIAMETER x 4 3/4" MIN. HOLE DEPTH (4" EFFECTIVE EMBEDMENT) AND 60 FT-LBS INSTALLATION TORQUE; 12" EDGE DISTANCE. ANCHOR SPACING IS SHOWN ON THE DRAWINGS INCLUDED IN THIS OPM.
- 8. STELLAR XL SURGICAL LIGHTS THAT ARE PART OF THIS OPM ARE ATTACHED TO THE ERGON 3 BOOMS AND HUB.

OPM-0256-13

SKYTRON ERGON 3 GENERAL NOTES



OSHPD OPM-0256-13 DWG - 1

SCALE

N.T.S.

PAGE

GEN NOTES

OSHPD OPM-0256-13

MANUFACTURE: SKYTRON

EQUIPMENT TYPE: CEILING MOUNTED BOOM AND SURGICAL LIGHTS (NON-MOTORIZED)

ATTACHMENT NOTES:

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.
- 2. BRACE ARM INCLINATION MAY VARY FROM 30° TO 45° FROM HORIZONTAL.
- 3. PERIODIC SPECIAL INSPECTION PER CBC 2013 SECTION 1705A AND TABLE 1705A.3 INCLUDING VERIFICATION OF ANCHOR TYPE. ANCHOR DIMENSIONS, CONCRETE TYPE, CONCRETE COMPRESSIVE STRENGTH, ANCHOR SPACING, EDGE DISTANCES, CONCRETE MEMBER THICKNESS, TIGHTENING TORQUE, HOLE DIMENSIONS, ANCHOR EMBEDMENT AND ADHERENCE TO THE MANUFACTURER'S PRINTED INSTALLATION INSTRUCTIONS. IN ADDITION, FOLLOW THE PROVISIONS OF THE 2013 CALIFONIA BUILDING CODE SECTION 1913A.7.2 BY CONFIRMING THE INSTALLATION TORQUE SPECIFIED BY THE MANUFACTURER. TESTING IS NOT TO OCCUR UNTIL A MINIMUM OF 24 HOURS HAS ELAPSED AFTER THE INSTALLATION OF THE SUBJECT ANCHORS. TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR. TEST 50% OF THE ANCHORS FOR EACH PIECE OF EQUIPMENT. USING A CALIBRATED TORQUE WRENCH VERIFY THE INSTALLATION TORQUE IS OBTAINED WITHIN 1/2 TURN OF THE NUT. REPORT OF TEST RESULTS ARE TO BE SUBMITTED TO THE ENGORCEMENT AGENGY. THE SEOR SHALL POVIDE REMEDIAL ANCHORAGE DETAILS IN THE EVENT THAT AN ANCHOR FAILS TO MEET THE TEST REQUIREMENTS.
- 4. WELDS ARE TO BE VISUALLY INSPECTED BY A QUALIFIED WELDING INSPECTOR.
- 5. STRENGTH DESIGN WAS USED FOR ANCHOR FORCE CALCULATIONS INCLUDING Ω_0 PER ACI 318-11 WHERE REQUIRED FOR ATTACHMENT TO CONCRETE.
- 6. EXCERCISE DUE CARE WHEN DRILLING POST-INSTALLED ANCHORS TO AVOID DAMAGING CONCRETE REINFORCEMENT OR
- 7. PROVIDE FULL THREAD ENGAGEMENT OF NUT AND WASHER.

RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD

- 1. CONFIRM THE MATERIAL PROPERTIES AND THICKNESS OF THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ATTACHED MEETS THE REQUIREMENTS OF THIS OPM.
- 2. PROVIDE A PLAN FOR INSPECTION OF SUPPORTS AND ATTACHMENTS AND VERIFY ITS IMPLEMENTATION.

OSI

- 3. CONFIRM THE SPECIFIED MINIMUM CONCRETE EDGE DISTANCES ARE MAINTAINED BASED ON THE ACTUAL EQUIPMENT LOCATION. VERIFY THAT EXISTING OR NEW ANCHORS ARE AN ADEQUATE DISTANCE FROM THIS UNIT'S ATTACHMENT.
- 4. VERIFY THAT THE EXISITING STRUCTURE IS ADEQUATE FOR THE IMPOSED DEAD, LATERAL AND TENSION FORCES SHOWN IN ADDTION TO ALL OTHER LOADS.
- 5. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH CBC 2013 AND WITH THE OPM-0256-13 DETAILS INCLUDING MATERIALS AND DIMENSIONS OF THE SUPPORT WHERE THE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN.
- 6. VERIFY THAT THE PROJECT SPECIFIC SDS AND Z/h VALUES RESULT IN SEISMIC FORCES (Eh AND EV) DO NOT EXCEED THE VALUES SHOWN IN THESE DETAILS.

SKYTRON ERGON 3 ATTACHMENT NOTES OPM-0256-13



A Division of Tomarco Contractor Specialties 1020 Crews Road, Suite Q, Matthews, NC 28105 704-841-4080 www.isatsb.com




DRAWN BY: WV.J **DATE**: 08/10/15


REVISED BY: WV.J **DATE:** 03/01/16


REV NO: 2


SCALE


PAGE N.T.S. ATTACH NOTES




















		Seismic Design Forces at Ultimate Loads							
Model Number	Weight at 1g (lbs.)	Vertical Force (lbs.)	Lateral force (lbs.)	Mx (ftlbs.)	My (ftlbs.)				
100 SERIES	216	367 C O D	648	1,481	2,401				
110 SERIES	494	840	C1,482	4,477	6,166				
120 SERIES	639	1,086	1,917	5,454	7,221				
200 SERIES	597	615	1,791	3,118	7,148				
300 SERIES	625 47	1,068	1,884	7,012	10,422				
310 SERIES	839	DPM- 4.426 6-13	2,517	6,137	12,460				
320 SERIES	1006	1,710	3,018	9,557	14,543				
330 SERIES	1168	1,986	3,504	12,059	19,546				

- 1. WEIGHTS AND MOMENTS ARE FACTORED LOADS USING STRENGTH DESIGN AND INCLUDE THE FOLLOWING FACTORS: DL = 1.2, FpV = 0.5 AND FpH = 3.0. MAXIMUM LOADING AFFECT IS USED FOR DESIGN.
- 2. Mx ARE MOMENTS IN THE HORIZONTAL AXIS DUE TO COMPONENT WEIGHT TIMES THE MAXIMUM HORIZONTAL ECCENTRICITY.
- 3. My are moments in the horizontal axis due to component weight plus vertical seismic acceleration times the maximum vertical eccentricity.

OPM-0256-13 SKYTRON ERGON 3 FORCES AND MOMENTS



A Division of Tomarco Contractor Specialties
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com




DRAWN BY: WVJ
DATE: 08/10/15
REVISED BY: WVJ

DATE: 11/15/16

REV NO: 3

SCALE PAGE
N.T.S. FORCES



Equipment Data						Seismic Design Forces at Ultimate Loads				
Model Series	Model Number	Hub	Bolts at	Bolts at	Torque	Weight at 1g	Vertical Force	Lateral Force	Mx	Му
	Used for Design	Conf. No.	at Hub	Mount. Pl.	ftlbs.	lbs.	lbs.	lbs.	ftlbs.	ftlbs.
STANDARD	STX555	2	(6) 1/2"	(4) 7/8"	32	353	600	1,059	2,653	3,960
EXTENDED ARM SINGLE	STX5-EA-2000	1	(4) 7/8"	(4) 7/8"	166	159	270	477	1,599	1,432
EXTENDED ARM DUAL	STX55-EA-2000/1500	1	(4) 7/8"	(4) 7/8"	166	315	536	945	1,601	1,389
EXTENDED HUB & ARM	STX55-EH-2000/1500	3	(6) 7/8"	(6) 7/8"	166	421	716	1,263	3,469	6,503
EXTENDED HUB	STX555-EH	3	(6) 7/8"	(6)7/8"	166	466	792	1,398	3,075	8,352
LFS	LFS STX555	1	(4) 7/8"	(4) 7/8"	166	562	955	1,686	3,656	6,136
LFSLFS	LFSLFS STX55	1	(4) 7/8"	(4) 7/8"	166	579	984	1,737	3,493	5,641

- 1. WEIGHTS AND MOMENTS ARE FACTORED LOADS USING STRENGTH DESIGN AND INCLUDE THE FOLLOWING FACTORS: DL = 1.2, FpV = 0.5 AND FpH = 3.0.
- 2. Mx ARE MOMENTS IN THE HORIZONTAL AXIS DUE TO COMPONENT WEIGHT TIMES THE MAXIMUM HORIZONTAL ECCENTRICITY.
- 3. My ARE MOMENTS IN THE HORIZONTAL AXIS DUE TO COMPONENT LATERAL FORCE TIMES THE MAXIMUM VERTICAL ECCENTRICITY.
- 4. MODEL NUMBER TABULATED YIELDS THE WORST CASE LOADING FOR THE MODEL SERIES.
- 5. EXTENDED HUB MODELS DO NOT HAVE A UNIVERSAL MOUNTING PLATE.
- 6. SEE PAGE "ATTACH NOTES" FOR A LISTING OF THE MODEL NUMBERS INCLUDED WITHIN A GIVEN MODEL SERIES.

ERIES MODELS INCLUDED IN OPM-0256-13

STANDARD STX5-895, STX5-1090, STX5-1285, STX5-895/700, STX55-1090/895, STX55-1285/1090, STX555, STX5-HY-895, STX5-HY-1090

LFS LFSSTXS, LFSSTXSS, LFSSTXSS 1 E: 12/19/2016

EXTENDED ARM - SINGLE STX5-EA-1500, STX5-EA-2000

EXTENDED ARM - DUAL STX5-HY-1285, STX-HY-1500, STX55-EA-1500/1300, STX55-EA-2000/1500

EXTENDED HUB/EXTENED ARM STX55-EH-2500/1300, STX55-EH-2000/1500

EXTENDED HUB STX5-EH-895, STX5-EH-1090, STX5-EH-1285, STXE-EH-1500, STX5-EH-2000, STX55-895/700, STX55-EH-1090/895, STX55-EH-1285/1090, STX55-EH

OPM-0256-13

SKYTRON ERGON 3 WITH STELLAR XL SURGICAL LIGHTS FORCES AND MOMENTS



International Seismic Application Technology
A Division of Tomarco Contractor Specialties
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com



DRAWN BY: WVJ DATE: 01/22/16

REVISED BY: WVJ DATE: 11/15/16

REV NO: 2

SCALE N.T.S. S

PAGE STELLAR XL

EQUIPMENT DATA AND ATTACHMENT AND ANCHORAGE FORCES ³									
Model Series	Bolts at Hub	Torque	Diagonal Brace	Top Plate Max. 2		Brace Anchoragge ²			
	and Mount. Pl.	ftlbs.	Anchorage 1	Tu Lbs	Vu Lbs	Tu Lbs	Vu Lbs		
ERGON 3									
100 SERIES	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	300	2260	1210	1290		
110 SERIES	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	300	2260	1210	1290		
120 SERIES	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	300	2260	1210	1290		
200 SERIES	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	300	2260	1210	1290		
300 SERIES	(6) 1 1/4"	500	(6) Hilti KB TZ 5/8"	500	2950	1100	1460		
310 SERIES	(6) 1 1/4"	500	(6) Hilti KB TZ 5/8"	500	2950	1100	1460		
320 SERIES	(6) 1 1/4"	500	(6) Hilti KB TZ 5/8"	500	2950	1100	1460		
330 SERIES	(6) 1 1/4"	500	(6) Hilti KB TZ 5/8"	500	2950	1100	1460		
	STELLAR XL SURGICAL LIGHTS								
STANDARD	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	230	2160	870	1310		
EXTENDED ARM SINGLE	(6) 1 1/4"	500	(4) Hilti-KB TZ 5/8"	230	2160	870	1310		
EXTENDED ARM DUAL	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	230	2160	870	1310		
EXTENDED HUB & ARM	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	230	2160	870	1310		
EXTENDED HUB	(6) 4 1/4"	500	(4) Hilti KB TZ 5/8"	230	2160	870	1310		
LFS	(6) 1 1/4"	500	(4) Hilti KB TZ 5/8"	230	2160	870	1310		
LFSLFS	(6) 1 1/4"OPI	I- 500 5	(4) Hilti KB TZ 5/8"	230	2160	870	1310		

- 1. SEE PAGE "GEN NOTE" GENERAL NOTE 7 HILTEKWIK BOLT TZ (ESR-1917) WITH 4" EFFECTIVE EMBEDMENT AND 60 FT-LBS INSTALLATION TORQUE.
- 2. INCLUDES AN OVERSTRENGTH FACTOR OF $\Omega_0 = 1.50_{\circ} / 2.016$
- 3. FORCES SHOWN ARE FOR EACH ANCHOR.

OPM-0256-13 SKYTRON ERGON 3 WITH STELLAR XL SURGICAL LIGHTS EQUIPMENT MOUNTING AND ANCHOR FORCES



International Seismic Application Technology
A Division of Tomarco Contractor Specialties
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com



DRAWN BY: WVJ
DATE: 02/26/16
REVISED BY: WVJ

DATE: 11/15/16 **REV NO:** 3

SCALE

N.T.S. ANCH FORCES

PAGE