



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL OF
MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY

APPLICATION #: OPM-0259-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: [X] New [] Renewal [] Update to Pre-CBC 2013 OPA Number:

Manufacturer Information

Manufacturer: Grifols Diagnostic Solutions, Inc.

Manufacturer's Technical Representative: Adam Ostergren

Mailing Address: 4560 Horton Street, Emeryville, CA 94608

Telephone: 510-923-5412 Email: Adam.ostergren@grifols.com

Product Information

Product Name: Erytra Analyzer

Product Type: Blood Analyzer for immunohematology testing

Product Model Number: None

General Description: Instrument for performing pre-transfusion compatibility tests

Applicant Information

Applicant Company Name: Grifols Diagnostic Solutions, Inc.

Contact Person: Adam Ostergren

Mailing Address: 4560 Horton Street, Emeryville, CA 94608

Telephone: 510-923-5412 Email: Adam.ostergren@grifols.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

Signature of Applicant: [Handwritten Signature] Date: 24 AUG 15

Title: Assoc. Dir. North America Company Name: GRIFOLS

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

Registered Design Professional Preparing Engineering Recommendations

Company Name: CYS Structural Engineers, Inc.

Name: Dieter T. Siebald California License Number: S4346

Mailing Address: 2495 Natomas Park Drive, Suite #650, Sacramento, CA 95833

Telephone: 916-920-2020 Email: dieters@cyseng.com

OSHPD Special Seismic Certification Preapproval (OSP)

Special Seismic Certification is preapproved under OSP-
(Separate application for OSP is required)

Special Seismic Certification is not preapproved

Certification Method(s)

Testing in accordance with: ICC-ES AC156 FM 1950-10

Other* (Please Specify): _____

*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing.

Analysis

Experience Data

Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

List of Attachments Supporting the Manufacturer's Certification

Test Report Drawings Calculations Manufacturer's Catalog

Other(s) (Please Specify): _____

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY

Signature:  Date: 12-14-2015

Print Name: Jeffrey Kikumoto

Title: SSE

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"



TABLE OF CONTENTS
OPM-0259-13

	PAGE
GENERAL NOTES	2
ABBREVIATIONS & DESIGN CRITERIA	5
EQUIPMENT DIMENSIONS	6
ANCHORAGE BRACKET LOCATIONS	7
FRONT FLOOR BRACKET DETAIL	8
REAR FLOOR BRACKET DETAIL	9
REAR EQUIPMENT BRACKET DETAIL	10
ATTACHMENT DETAILS TO EQUIPMENT FRAME	11
ANCHORAGE DETAIL TO CONCRETE FILL OVER METAL DECK (CASE 1)	12
ANCHORAGE DETAIL TO CONCRETE SLAB (CASE 2)	15

- NOTES:**
1. THESE DRAWINGS ARE PREPARED FOR GRIFOLS DIAGNOSTIC SOLUTIONS, INC., EMERYVILLE, CALIFORNIA.
 2. THE CONTRACTOR AND INSPECTOR OF RECORD SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OSHPD PRE-APPROVAL PROGRAMS WEBSITE.
 3. THIS PRE-APPROVAL COVERS THE SUPPORTS & ATTACHMENTS OF THE UNIT TO THE SUPPORTING STRUCTURE. THE UNIT, ANCHORAGE BRACKETS & CONNECTION HARDWARE TO UNIT IS SUPPLIED & INSTALLED BY GRIFOLS. THRU-BOLTS, STRUT HARDWARE & ATTACHMENTS AT SOFFIT UNDER METAL DECK & EXPANSION BOLTS SHOWN ON PAGES 12 TO 15 SHALL BE SUPPLIED & INSTALLED BY THE CONTRACTOR.



SHEET TITLE: TABLE OF CONTENTS



CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833


TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	1 of 15

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:51pm Login:meyerhoferm Dimscale:1 L1Scale:6

GENERAL NOTES:

1. THIS OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2013. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2013.
2. IT IS THE RESPONSIBILITY OF THE STRUCTURAL ENGINEER OF RECORD FOR A SITE SPECIFIC PROJECT TO VERIFY:
 - A. THE ADEQUACY OF THE NEW OR EXISTING STRUCTURE TO RESIST THE FORCES & WEIGHT SPECIFIED FOR EACH EQUIPMENT IN ADDITION TO ALL OTHER LOADS. PROVIDE & DESIGN SUPPLEMENTARY MEMBERS AS REQUIRED.
 - B. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS.
 - C. THAT THE FLOOR ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. THE SPACING SHOWN IN THE TEST LOADS TABLE ON PAGE 3 IS THE REQUIRED MINIMUM SPACING OF THE GIVEN DIAMETER ANCHORS. THE REQUIRED SPACING FROM ANCHORS OF OTHER DIAMETERS & EMBEDMENTS MAY VARY & SHALL BE EVALUATED BY THE SEOR.
 - D. THAT THE INSTALLATION IS IN CONFORMANCE WITH THE CBC 2013 & WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL.
 - E. THAT THE ACTUAL EQUIPMENT'S WEIGHT, CENTER OF GRAVITY (CG) LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS, & THE MATERIAL & GAGE OF THE EQUIPMENT WHERE ATTACHMENTS ARE MADE, AGREE WITH THE INFORMATION SHOWN ON THE PRE-APPROVAL DOCUMENTS.
3. EXPANSION ANCHORS INSTALLED IN NORMAL WEIGHT OR SAND-LIGHTWEIGHT CONCRETE SHALL BE CARBON STEEL HILTI KB-TZ EXPANSION ANCHORS COMPLYING WITH ESR-1917 REISSUED MAY 2015.
 - A. INSTALLATION: INSTALL THE EXPANSION ANCHORS IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC EVALUATION REPORT FOR THE SPECIFIC ANCHOR & THE PARAMETERS GIVEN IN THE TABLE ON PAGE 3. PROVIDE FULL THREAD ENGAGEMENT FOR NUT & WASHER.
 - B. JOB TESTING: FOR VERIFYING SATISFACTORY INSTALLATION WORKMANSHIP, PERFORM JOB SITE TESTING IN ACCORDANCE WITH THE TEST LOAD TABLE PROVIDED IN THIS DOCUMENT. TEST 50% OF THE INSTALLED ANCHORS. THE TEST LOAD MAY BE APPLIED BY ANY METHOD THAT WILL EFFECTIVELY MEASURE THE TENSION OR TORQUE IN THE ANCHOR SUCH AS DIRECT PULL WITH A HYDRAULIC JACK OR CALIBRATED SPRING LOADING DEVICES OR CALIBRATED TORQUE WRENCH METHOD. ALL TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE INSPECTOR OF RECORD (IOR). IF ANY ANCHOR FAILS THE TEST, TEST ALL ANCHORS. THE TEST SHALL BE PERFORMED 24 HOURS OR MORE AFTER INSTALLATION. TESTING MAY BE DONE PRIOR TO EQUIPMENT INSTALLATION. ALSO REFER TO CBC 1913A.7 "FIELD TESTS FOR POST-INSTALLED ANCHORS IN CONCRETE".
 - C. FAILURE/ACCEPTANCE CRITERIA: THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:
 - HYDRAULIC RAM METHOD: APPLY & HOLD TEST LOAD FOR A MINIMUM OF 15 SECONDS. THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE APPLICABLE TEST LOAD WHERE WASHERS ARE USED. FOR WEDGE TYPE ANCHORS, SUCH AS HILTI KB-TZ, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER UNDER THE NUT BECOMES LOOSE.
 - TORQUE WRENCH METHOD: THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: ONE-HALF (1/2) TURN OF THE NUT.




SHEET TITLE: GENERAL NOTES

 <p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	<p>TEL (916) 920-2020 www.cyseng.com</p>	<p>Job No: 15077 Date: 12/10/2015 Page: 2 of 15</p>
--	--	---

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:52pm Login:meyerhoferm Dimscale:1 L1Scale:6

GENERAL NOTES CONTINUED:

3D. TEST VALUES: APPLY TEST LOADS TO ANCHORS WITHOUT REMOVING THE NUT IF POSSIBLE, SEE TABLE BELOW.



ANCHOR DIA (INCH) d_a	INSTALLATION EMBED (INCH) h_{nom}	EFFECTIVE EMBED (INCH) h_{ef}	HOLE DEPTH (INCH) h_o	MIN CONC THICKNESS (INCH) h_{min}	MIN CONC EDGE DISTANCE (INCH)	MIN AB SPACING (INCH)	TEST LOAD		CONDITION OF ANCHORAGE
							TENSION LOAD (LBS)	TORQUE (FT-LBS)	
1/2	2 3/8	2	2 5/8	4	12	4	1660	40	CASE 2

4. BOLTS THROUGH CONC ON MTL DECK:

- A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUT AFTER SNUG TIGHT CONDITION IS ACHIEVED, UNO. THE SNUG TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT.
- B. THRU-BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16")
- C. THRU-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION & TESTING IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS. THRU-BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TESTING.

5. SCREW ANCHORS TO BOTT OF CONC FILL OVER MTL DECK: 5

- A. HILTI KH-EZ (ICC ESR-3027) TENSION TEST LOAD FOR CASE 1.

ANCHOR DIA (INCH) d_a	INSTALLATION EMBED (INCH) h_{nom}	EFFECTIVE EMBED (INCH) h_{ef}	HOLE DEPTH (INCH) h_o	MIN CONC THICKNESS (INCH) h_{min}	MIN CONC EDGE DISTANCE (INCH)	MIN AB SPACING (INCH)	TENSION TEST LOAD (LBS)
1/4	1 5/8	1.18	2	3/4	1 1/4*	10*	400

* SEE PG 13 IN THIS OPM & FOOTNOTE 2, TABLE 2 IN ESR-3027




SHEET TITLE: GENERAL NOTES (CONTINUED)

 <p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020	Job No: 15077
	www.cyseng.com	Date: 12/10/2015
		Page: 3 of 15

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:53pm Login:meyerhoferm Dimscale:1 LITScale:6

GENERAL NOTES CONTINUED:


6. TWO (2) CASES OF ATTACHMENT ARE SPECIFIED & PRESENTED IN THIS PRE-APPROVAL:



CASE 1: ATTACHMENT DETAILS LOCATED AT UPPER FLRS ABV THE BASE OF A BLDG ($z/h \leq 1$). THE FLRS ARE ASSUMED TO BE BUILT OF A MIN 3/4" SLWC TOPPING OVER 20 GA MTL DECK ($f'c = 3000$ PSI, MIN).

CASE 2: ATTACHMENT DETAILS LOCATED AT OR BLW THE BASE OF A BLDG ($z/h=0$). THE FLRS ARE ASSUMED TO BE BUILT OF A MIN 4" NWC SLAB ($f'c = 3000$ PSI, MIN).

7. THIS PRE-APPROVAL MAY BE USED AT ANY GEOGRAPHICAL LOCATION IN THE STATE OF CALIFORNIA WHERE S_{DS} IS LESS THAN OR EQUAL TO 2.50.



SHEET TITLE: GENERAL NOTES (CONTINUED)

 <p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	<p>TEL (916) 920-2020 www.cyseng.com</p>	<p>Job No: 15077 Date: 12/10/2015 Page: 4 of 15</p>
--	--	---

L:\Jobs\15\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:54pm Login:mmyerhoferm DimScale:1 LITScale:6

ABBREVIATIONS:

@	AT	LBS	POUNDS
ABV	ABOVE	LLH	LONG LEG HORIZONTAL
ASTM	AMERICAN SOCIETY FOR TESTING & MATERIALS	LLV	LONG LEG VERTICAL
BLDG	BUILDING	LRFD	LOAD AND RESISTANCE FACTOR DESIGN
BLW	BELOW	MAX	MAXIMUM
CBC	CALIFORNIA BUILDING CODE	MFR	MANUFACTURER
CG	CENTER OF GRAVITY	MIN	MINIMUM
☉	CENTERLINE	MTL	METAL
CONC	CONCRETE	NO. (#)	NUMBER OR POUNDS
COORD	COORDINATE	NWC	NORMAL WEIGHT CONCRETE
DBL	DOUBLE	OSHPD	OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT
DIA (φ)	DIAMETER	PSI	POUNDS PER SQUARE INCH
(E)	EXISTING CONDITION	SLWC	SAND-LIGHTWEIGHT CONCRETE
ELEV	ELEVATION	SS	STAINLESS STEEL
EQUIP	EQUIPMENT	STD	STANDARD
f'c	MINIMUM ULTIMATE COMPRESSIVE STRENGTH OF CONCRETE	STL	STEEL
FLR	FLOOR	Tu	ANCHORAGE TENSION REACTION DUE TO SEISMIC FORCE AT LRFD
FT (')	FOOT/FEET	THRD	THREAD OR THREADED
Fy	SPECIFIED MINIMUM YIELD STRESS OF STEEL	TYP	TYPICAL
GA	GAUGE	Vu	ANCHORAGE SHEAR REACTION DUE TO SEISMIC FORCE AT LRFD
IN (")	INCH	W/	WITH
KSI	KIPS PER SQUARE INCH	Wp	OPERATING WEIGHT
		WT	WEIGHT

DESIGN CRITERIA

SUPPORT & ATTACHMENT DESIGN IS PER 2013 CBC AT LRFD LEVEL FORCES.

OTHER MECHANICAL OR ELECTRICAL COMPONENTS PER TABLE 13.6-1 OF ASCE 7-10 SUPPLEMENT #1:

$$\alpha_p = 1.0 \quad R_p = 1.5 \quad I_p = 1.5 \quad \Omega_0 = 1.5 \text{ (CONC ANCHORS)}$$

Wp AS NOTED ON EQUIP DRAWING SHOWN ON PG 6

FOR CASE 1 – UPPER FLRS ABV THE BASE, z/h = 1

$$S_{Ds} = 2.50 \quad F_p = 3.00 W_p$$


FOR CASE 2 – SLAB AT OR BLW BASE, z/h = 0

$$S_{Ds} = 2.50 \quad F_p = 1.13 W_p$$

LOAD COMBINATIONS

$$(0.9 - 0.2 S_{Ds}) D - \Omega_0 F_p \text{ (FOR MAXIMUM TENSION)}$$


$$(1.2 + 0.2 S_{Ds}) D + \Omega_0 F_p \text{ (FOR MAXIMUM COMPRESSION)}$$



SHEET TITLE: ABBREVIATIONS & DESIGN CRITERIA


 <p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020	Job No: 15077
	www.cyseng.com	Date: 12/10/2015
		Page: 5 of 15

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:54pm Login:meyerhoferm DimScale:1 L1Scale:6



FRONT ELEV

SIDE ELEV



PLAN VIEW

OPM-0259-13

BY: Jeffrey Y. Kikumoto

DATE: 12/14/2015



SHEET TITLE: EQUIPMENT DIMENSIONS




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	6 of 15



NOTE:

SEE PG 11 FOR ATTACHMENT DETAILS TO EQUIP FRAME.
SEE PGS 12 & 15 FOR ANCHORAGE DETAILS TO SUPPORTING FLR.



SHEET TITLE: ANCHORAGE BRACKET LOCATIONS




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com


Job No:	15077
Date:	12/10/2015
Page:	7 of 15

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:55pm Login:meyerhoferm Dimscale:1 L1Scale:6



ELEV VIEW

SECTION



PLAN VIEW



SHEET TITLE: FRONT FLOOR BRACKET DETAIL




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833


TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	8 of 15


L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:56pm Login:meyerhoferm Dimstyle:1 LTScale:6




ELEV VIEW



SECTION



PLAN VIEW



SHEET TITLE: REAR FLOOR BRACKET DETAIL




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833


TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	9 of 15


L:\Jobs\15\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:56pm Login:myerhoferm Dimscale:1 LTSscale:6




ELEV VIEW



SECTION



PLAN VIEW



SHEET TITLE: REAR EQUIPMENT BRACKET DETAIL




CYS STRUCTURAL ENGINEERS, INC.


2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	10 of 15



DETAIL A
ATTACHMENT AT FRONT OF EQUIP



DETAIL B
ATTACHMENT AT REAR OF EQUIP



SHEET TITLE: ATTACHMENT DETAILS
TO EQUIPMENT FRAME



CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com


Job No:	15077
Date:	12/10/2015
Page:	11 of 15

L:\Jobs\15\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:57pm Login:meyerhoferm Dimscale:1 LTScale:6

MAX ANCHOR FORCES
AT LRFD AT EA AB

	Tu	Ω_o Vu
CASE 1 $z/h \leq 1$	1545#	594#

OVERSTRENGTH FACTOR (Ω_o)
INCLUDED WHERE NOTED.



SHEET TITLE: ANCHORAGE DETAIL
TO CONCRETE FILL OVER METAL DECK (CASE 1)




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	12 of 15

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:57pm Login:myerhoferm DimScale:1 LITScale:6



SHEET TITLE: STRUT PLATE DETAIL




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	15077
Date:	12/10/2015
Page:	13 of 15

EXTEND STRUT LENGTH TO NEXT ADJ LOW FLUTE IF AB'S ARE LESS THAN 2" FROM STRUT ϕ SUPPORT SCREWS
LENGTH SHALL ENGAGE 2 LOW FLUTES MIN



STRUT ϕ SUPPORT SCREWS:
HILTI KH-EZ $\frac{1}{4}$ " ϕ x $2\frac{5}{8}$ " CONC SCREWS (ICC ESR-3027) 1 EA END OF ϕ . USE 18 FT-LBS MAX INSTALLATION TORQUE USING A CALIBRATED TORQUE WRENCH, TYP

AB PER CASE 1

$\frac{5}{16}$ " STRUT ϕ , TYP. FIELD VERIFY DECK ORIENTATION FOR COORD W/ AB LAYOUT


1.00" MAX EITHER SIDE OF FLUTE ϕ , TYP

PLAN VIEW

ANCHORS PERPENDICULAR TO FLUTES

EXTEND STRUT LENGTH TO NEXT ADJ LOW FLUTE IF AB'S ARE LESS THAN 2" FROM STRUT ϕ SUPPORT SCREWS

LENGTH SHALL ENGAGE 2 LOW FLUTES MIN



STRUT ϕ SUPPORT SCREWS:
HILTI KH-EZ $\frac{1}{4}$ " ϕ x $2\frac{5}{8}$ " CONC SCREWS (ICC ESR-3027) 1 EA END OF ϕ


AB PER CASE 1

$\frac{5}{16}$ " STRUT ϕ . FIELD VERIFY DECK ORIENTATION FOR COORD W/ AB LAYOUT

BOTT OF MTL DECK (LOW FLUTES SHADED FOR CLARITY)

PLAN VIEW

ANCHORS PARALLEL TO FLUTES



SHEET TITLE: STRUT PLATE DETAIL



CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com


Job No:	15077
Date:	12/10/2015
Page:	14 of 15

L:\Jobs\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:58pm Login:myerhoferm DimScale:1 L1Scale:6

MAX ANCHOR FORCES
AT LRFD AT EA AB

	$\Omega_o Tu$	$\Omega_o Vu$
<u>CASE 2</u> $z/h = 0$	1331#	359#

OVERSTRENGTH FACTOR (Ω_o)
INCLUDED WHERE NOTED.



SHEET TITLE: ANCHORAGE DETAIL
TO CONCRETE SLAB (CASE 2)

 <p>CYS STRUCTURAL ENGINEERS, INC. 2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 15077 Date: 12/10/2015 Page: 15 of 15
--	--------------------------------------	---

L:\Jobs\15\15077 Grifols Diagnostics- Erytra Analyzer OPM\STRU\S1.dwg Time:Dec10,2015-04:58pm Login:meyerhoferm Dimscale:1 LITScale:6