

APPLICATION FOR OSHPD PREAPPROVAL

## OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

OF MANUFACTURER'S CERTIFICATION (OPM)  APPLICATION #: OPM-0290-13										
OSHPD Preapproval of Manufacturer's Certification (OPM)										
Type: ☐ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number:										
Manufacturer Information										
Manufacturer: _Draeger Medical, Inc.										
Manufacturer's Technical Representative: Merouane Djerbal										
Mailing Address: 3135 Quarry Road, Telford, PA. 18969										
Telephone: (215) 412-5868 Email: Merouane.Djerbal@draeger.com										
Product Information										
Product Name: Agila-Movita System										
Product Type: Overhead OPM-0290-13										
Product Model Number: N/A										
General Description: Medical gas, electrical and misc ceiling supported supply arm(s)										
DATE: 10/11/2019										
Applicant Information										
Applicant Company Name: EASE Co.										
Contact Person: Jonathan Roberson, S.E.										
Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709										
Telephone: (909) 606-7622 Email: J.Roberson@EASECo.com  I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2016.										
Signature of Applicant: Date: 12/8/15										
Title: Principal Engineer Company Name: EASE Co.										

14/1/WW



STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY OSH-FD-700 (REV 1/24/13)

Page 1 of 2

OFFICE USE ONLY

<sup>&</sup>quot;Access to Safe Quality Healthcare Environments that Meet California's



### OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations									
Company Name: EASE Co.									
Name: Jonathan Roberson, S.E. California License Number: S4197									
Mailing Address: _5877 Pine Ave. Suite 210, Chino Hills, CA. 91709									
Telephone: 909-606-7667 Email: J.Roberson@EASECo.com									
OSHPD Special Seismic Certification Preapproval (OSP)									
<ul> <li>□ Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required)</li> <li>□ Special Seismic Certification is not preapproved</li> </ul>									
Certification Method(s)									
☐ Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-16 ☐ Other* (Please Specify):									
*Use of test criteria other than those adopted by the California Building Standards Code, 2016 (CBSC 2016) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2016 may be used when approved by OSHPD prior to testing.  BY: Jeffrey Kikumoto									
Experience Data  DATE: 10/11/2019									
Combination of Testing, Analysis, and/or Experience Data (Please Specify):									
List of Attachments Supporting the Manufacturer's Certification									
☐ Test Report ☐ Drawings ☐ Calculations ☐ Manufacturer's Catalog ☐ Other(s) (Please Specify):									
OFFICE USE ONLY - OSHPD APPROVAL VALID FOR CBC 2016 & ALL PRE-2016 CODE BASED PROJECTS									
Signature: Date:									
Print Name:									
Title: Structural Engineer									
Condition of Approval (if applicable):									

"Access to Safe Quality Healthcare Environments that Meet California's







5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development
PREAPPROVAL OF MANUFACTURER'S CERTIFICATION
OPM-0290-13

THIS PREAPPROVAL CONFORMS TO THE 2016 CALIFORNIA BUILDING CODE

MANUFACTURER: DRAI

DRAEGER

**EQUIPMENT NAME:** 

**AGILA-MOVITA SYSTEM** 

Sheet: 1 of 10 Date: 10/10/19

#### **GENERAL NOTES**

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2016 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2016 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2016 CALIFORNIA BUILDING CODE.
- 4. FORCES PER ASCE 7-10 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE SDS = 1.10,  $\mathbf{a}_P = 2.5$ ,  $\mathbf{l}_P = 1.5$ ,  $\mathbf{R}_P = 2.5$ ,  $\mathbf{z}/h \le 1$ .
- 5. THE DETAILS IN THIS PREAPPROVAL MAY BE USED AT ANY LOCATION IN THE STATE OF CALIFORNIA, WHERE SDS IS NOT GREATER THAN 1.10.
- 6. ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 8. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING
  - A. PROVIDE SUPPORTING STRUCTURE REQUIRED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN ADDITION TO ALL OTHER LOADS.
  - B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2016 CBC AND WITH THE DETAILS SHOWN IN THIS PREAPPROVAL. VERIFY THAT THE ACTUAL EQUIPMENT'S WEIGHT, CG LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS AND THE MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
  - C. VERIFY THAT THE COMBINATION OF SDS & z/h RESULT IN SEISMIC FORCES (Eh , Ev) THAT ARE NOT GREATER THAN THE VALUES ON THE DETAILS.



www.EquipmentAnchorage.com

#### DRAEGER

DES. J. ROBERSON

JOB NO.

DATE

11-1456

10/10/19

SHEETS

SHEET


#### AGILA-MOVITA SYSTEM

#### 10. EXPANSION ANCHORS:


A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension Test
1/2"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	3.25"	8"	24"	6"	40 FT-LB	3026 lb
3/4"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	4.75"	12"	24"	8"	110 FT-LB	5410 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 24" AWAY MINIMUM (i.e. - CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING OF EXPANSION ANCHORS PER 2016 CBC, 1910A.5: TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO OSHPD
  - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
  - (ii) ACCEPTANCE CRITERIA:
    - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
    - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
  - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.



TYPICAL CONCRETE EDGE DETAIL



DES. J. ROBERSON

www.EquipmentAnchorage.com

DRAEGER


11-1456 JOB NO.

SHEET

AGILA-MOVITA SYSTEM

10/10/19 DATE

SHEETS



1. FORCES ARE DETERMINED PER 2016 CALIFORNIA BUILDING CODE AND ASCE 7-10

STRENGTH DESIGN IS USED.(Sps = 1.10,  $a_p$  = 2.5,  $I_p$  = 1.5,  $R_p$  = 2.5,  $z/h \le 1$ )

HORIZONTAL FORCE (En) = 1.98 Wp VERTICAL FORCE (Ev) = 0.22 Wp

- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. SEE GENERAL NOTES: SHEET 1



## EASE

#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

**DRAEGER** 

DES. J. ROBERSON

11-1456

4

SHEET

AGILA-MOVITA SYSTEM


DATE 10/10/19

JOB NO.

<sub>F</sub> 10 <sub>SHEETS</sub>

#### SEISMIC SUPPORTS & ATTACHMENTS

#### COLUMN FLANGE TO INTERFACE PLATE





## EASE

#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

**DRAEGER** 

DES. J. ROBERSON

11-1456

SHEET 5

AGILA-MOVITA SYSTEM


DATE 10/10/19


JOB NO.

F 10 SHEETS

#### SEISMIC SUPPORTS & ATTACHMENTS

#### PATIENT LIFT ADAPTER PLATE TO INTERFACE PLATE





## **DRAEGER**

#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

DES. J. ROBERSON

11-1456 JOB NO.

DATE


10/10/19


SHEETS

SHEET

AGILA-MOVITA SYSTEM

SEISMIC SUPPORTS & ATTACHMENTS





# EASE

#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

DES. J. ROBERSON

www.EquipmentAnchorage.com

**DRAEGER** 

**ЈОВ NO.** 11-1456

SHEET 7


AGILA-MOVITA SYSTEM


DATE 10/10/19

OF 10 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

BRACING DETAILS





#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING** www.EquipmentAnchorage.com SHEET

### **DRAEGER**

DES. J. ROBERSON

11-1456 JOB NO.

SHEETS

### AGILA-MOVITA SYSTEM

10/10/19 DATE SEISMIC SUPPORTS & ATTACHMENTS CEILING MOUNTED USE (4)- HSS3x3x3/I6" (A500 GR B, Fy=46 KSI MIN) VERTICAL MEMBERS E70XX TYP / 3/16 USE (6)- MI8 (CL 10.9) BOLTS (TORQUE TO 249 FT-LB) STRUCTURAL (SECTION A-A) USE (8)- HSS3x3x3/16" (A500 GR B, Fu=46 KSI MIN) BRACING MEMBERS (6)- 20mm DIA HOLE 0 BRACE TO COLUMN STRUCTURAL PLATE 27" X I" X 2'-3" (A529 GR 50, Fy=50 KSI MIN) E70XX No. 4197 BRACE TO COLUMN 3/16 EXP. 6-30-2020 STRUCTURAL PLATE DETAIL

www.EquipmentAnchorage.com

DES. J. ROBERSON

11-1456 JOB NO.


DATE

SHEET

10/10/19 SHEETS

#### AGILA-MOVITA SYSTEM

**DRAEGER** 



www.EquipmentAnchorage.com

**DRAEGER** 

DES. J. ROBERSON

10/10/19

11-1456 JOB NO.

DATE

SHEET


SHEETS

CEILING MOUNTED

AGILA-MOVITA SYSTEM

SEISMIC SUPPORTS & ATTACHMENTS DIRECTION OF DECK FLUTES USE (4)- 1/2"Φ HILTI KB-TZ EXPANSION ANCHORS E70XX (MIN. EMBED. (her) = 3.25") (4 PER BRACE, 32 TOTAL) 3/8" PL (A36, Fy=50 KSI MIN) 13.75" BRACING MEMBER OPM-0290-13 ANCHOR AT BRACE PLATE  $T_{U1} = 778 LB/BOLT (MAX)$ Vu<sub>1</sub> = 654 LB/BOLT (MAX) BUILDING

BRACE PLATE DETAIL



(VALUES INCLUDE  $\Omega_{\circ}$ )