

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT **FACILITIES DEVELOPMENT DIVISION**

APPLICATION FOR OSHPD PREAPPROVAL OF	OFFICE USE ONLY								
MANUFACTURER'S CERTIFICATION (OPM)	APPLICATION #: OPM-0332								
OSHPD Preapproval of Manufacturer's Certification (OPM)									
Type: New X Renewal/Update									
Manufacturer Information									
Manufacturer: Shimadzu Medical Systems									
Manufacturer's Technical Representative: Jim Mekker									
Mailing Address: 25101 Chagrin Blvd., Suite 240, Beechwood, OH 44122									
Felephone: (216) 288-0709 Email: mekker@shimadzu-usa.com									
EOR CODE CO	1 .								
Product Information OSHPD	· ·								
Product Name: BK-200 X-RAY TABLE									
Product Type: Other Electrical & Mechanical Components	CH								
Product Model Number: BK-200									
General Description: Patient Table									
DATE: 09/29/2021	2013								
Applicant Information	â.								
Applicant Company Name: EASE LLC.	10 N								

Contact Person: Tiffany Tonn

Mailing Address: 1515 FAIRVIEW AVE, STE 205, MISSOULA, MT 59801

Telephone: (406) 541-3273 Email: tiffany@easeco.com

Title: Office Manager

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA - HEALTH AND HUMAN SERVICES AGENCY







OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations								
Company Name: EASE								
Name: Jonathan Roberson California License Number: S4197								
Mailing Address: 5877 Pine Ave., Suite 210, Chino Hills, CA 91709								
Telephone: (951) 295-1892 Email: jon@EASECo.com								
OSHPD Special Seismic Certification Preapproval (OSP)								
Special Seismic Certification is preapproved under OSP OSP Number:								
OR CODE								
Certification Method								
Testing in accordance with:								
Other(s) (Please Specify):								
*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHPD prior to testing.								
X Analysis BY: Hayne Kim								
Experience Data DATE: 09/29/2021								
Combination of Testing, Analysis, and/or Experience Data (Please Specify):								
TO THE CONTRACT OF THE CONTRAC								
OSHPD Approval BUILDING								
Date: 9/29/2021								
Name: Hayne Kim Title: Senior Structural Engineer								
Condition of Approval (if applicable):								

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"









EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development

PREAPPROVAL OF MANUFACTURER'S CERTIFICATION OPM-0332

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER: SHIMADZU MEDICAL SYSTEMS

EQUIPMENT NAME: **BK-200 X-RAY TABLE**

Sheet: 1 of 8 Date: 9/27/21

GENERAL NOTES

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 1.90, 2.00 & 2.30. SEE DETAIL FOR APPLICABILITY
- 4. FORCES PER ASCE 7-16 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3,
 - WHERE SDS = 1.90, a_p = 1.0, I_p = 1.5, R_p = 1.5, z/h = 0 AT CONCRETE SLAB. SEE FOLLOWING SHEETS FOR Ω_0
 - WHERE SDS = 2.30, a_p = 1.0, I_p = 1.5, R_p = 1.5, z/h = 0 AT CONCRETE SLAB. SEE FOLLOWING SHEETS FOR Ω_0
 - WHERE SDS = 2.00, a_p = 1.0, I_p = 1.5, R_p = 1.5, z/h < 1 AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_0
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING, (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT OR BELOW GRADE. (i.e. z/h = 0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR REPORT AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SHIMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

9/27/21

JOB NO. 11-2101

DATE

of 8 sheets


SHEET

10. EXPANSION ANCHORS:

A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension Test
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	6.75"	12"	See Detail "A"	30 FT-LB	N/A
1/2"	Normal Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	8"	14"	4"	50 FT-LB	1983 lb
1/2	Normal Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	3.25"	10"	24"	6"	50 FT-LB	3325 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, (SEE SCHEDULE) AWAY MINIMUM (i.e. CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING AND SPECIAL INSPECTION OF EXPANSION ANCHORS SHALL
 BE PERFORMED BY AN APPROVED INDEPENDENT AGENCY
 EMPLOYED BY THE FACILITY OWNER PER CBC 1704A & 1910A.5_ 0 3 3 2
 AND CAC 7-149. ALL REPORTS SHALL BE SENT TO THE INSPECTOR
 OF RECORD, OWNER AND THE ARCHITECT OR ENGINEER IN
 RESPONSIBLE CHARGE.
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION,
 DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF 2 0 2 1
 THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:
 - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
 - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
 - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
 - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
 - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL (SLAB ON GRADE ONLY)



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SHIMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

9/27/21

JOB NO. 11-2101


DATE

of 8 SHEETS

SHEET

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB



NOTES:

- 1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16 STRENGTH DESIGN IS USED. ($\alpha_p = 1.0$, $p_p = 1.5$, $p_p = 1.$
- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. BK-200 X-RAY TABLE (UUT-13) HAS OBTAINED SPECIAL SEISMIC CERTIFICATION REFER TO OSP-0319-10.
- 5. SEE GENERAL NOTES: SHEET 1 AND 2.

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SHIMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

9/27/21

JOB NO. 11-2101

DATE


4


<u>of 8 sheets</u>

SEISMIC SUPPORTS & ATTACHMENTS

MAX Sps ≤ 1.90

CONCRETE SLAB





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SHIMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

JOB NO. 11-2101

DATE 9/27/21

SHEET


5


OF **8** SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

1.90 < MAX Sps < 2.30

CONCRETE SLAB





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SHIMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

11-2101 JOB NO.


9/27/21 DATE

SHEET

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK



NOTES:

FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 2.00, $\Delta p = 1.0$, lp = 1.5, Rp = 1.5, $\Omega_0 = 1.5$, z/h < 1)

HORIZONTAL FORCE (Eh) = 2,40 Wp

HORIZONTAL FORCE (Emh) = 3.60 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (Ev) = 0.40 Wp

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. BK-200 X-RAY TABLE (UUT-13) HAS OBTAINED SPECIAL SEISMIC CERTIFICATION REFER TO OSP-0319-10.
- 5. SEE GENERAL NOTES: SHEET 1 AND 2.



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SHÍMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

11-2101 JOB NO.


9/27/21 DATE


SHEET

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

OF

SHIMADZU MEDICAL SYSTEMS

BK-200 X-RAY TABLE

DES. J. ROBERSON

JOB NO. 11-2101

DATE 9/27/21


SHEET

8 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

THREADED ROD OCCURS AT FLUTE)

CONCRETE DETAILS



MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL

ORWIA BUILDING

