

## DEPARTMENT OF HEALTH CARE ACCESS AND INFORMATION FACILITIES DEVELOPMENT DIVISION

all Mil No.	
<b>APPLICATION FOR HCAI PREAPPROVAL C</b>	OFFICE USE ONLY
MANUFACTURER'S CERTIFICATION (OPM)	
HCAI Preapproval of Manufacturer's Certification (OPM	
Type: New X Renewal/Update	
Manufacturer Information	
Manufacturer: Advanced Sterilization Products	
Manufacturer's Technical Representative: Curtis Felt	
Mailing Address: 33 Technology Drive, Irvine, CA 92618	
Telephone: (949) 507-6186 Email: curtis	.felt@asp.com
LED TO THE PARTY OF THE PARTY O	ONID
Product Information	
Product Name: EVOTECH OPM-0	0340
Product Type: Other Mechanical and Electrical Components	
Product Model Number: EVOTECH BY: Mohamm	nad Aliaari
General Description: Endoscope cleaner and reprocessor	4/0000
DATE: U6/1	4/2022
Applicant Information	3 4.
Applicant Company Name: EASE LLC.	CO
Contact Person: Tiffany Tonn	ING
Mailing Address: 1515 FAIRVIEW AVE, STE 205, MISSOULA, M	MT 59801

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





Telephone: (406) 541-3273

Title: Office Manager

Email: tiffany@easeco.com



# DEPARTMENT OF HEALTH CARE ACCESS AND INFORMATION FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations							
Company Name: EASE							
Name: Jonathan Roberson California License Number: S4197							
Mailing Address: 5877 Pine Ave., Suite 210, Chino Hills, CA 91709							
Telephone:         (951) 295-1892         Email:         jon@EASECo.com							
HCAI Special Seismic Certification Preapproval (OSP)							
Special Seismic Certification is preapproved under OSP OSP Number:							
2000							
Certification Method							
Testing in accordance with:							
Other(s) (Please Specify):							
*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by HCAI prior to testing.							
X Analysis							
Experience Data  C DATE: 06/14/2022							
Combination of Testing, Analysis, and/or Experience Data (Please Specify):							
COSE.							
HCAI Approval							
Date: 6/14/2022							
Name: Mohammad Aliaari Title: Senior Structural Engineer							
Condition of Approval (if applicable):							

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"







#### EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

The Department of Health Care Access and Information
PREAPPROVAL OF MANUFACTURER'S CERTIFICATION
OPM-0340

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER: ADVANCED STERILIZATION PRODUCTS

**EQUIPMENT NAME: EVOTECH** 

Sheet: 1 of 8

Date: 6/7/22

#### **GENERAL NOTES**

- 1. THIS HCAI PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 2.30.
- 4. FORCES PER ASCE 7-16 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE SDS = 2.30,  $\mathbf{a}_p$  = 1.0,  $\mathbf{I}_p$  = 1.5,  $\mathbf{r}_p$  = 1.5,
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 6. ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT OR BELOW GRADE. (i.e. z/h = 0)

#### 9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev ) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR REPORT. AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

### ADVANCED STERILIZATION PRODUCTS

DES. J. ROBERSON

11-2206

2 2

DATE 6/7/22

JOB NO.

OF 8 SHEETS


### **EVOTECH**

#### 10. EXPANSION ANCHORS:


A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension Test
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	6.75"	12"	See Detail "A"	30 FT-LB	N/A
3/8"	Normal Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	3"	12"	4"	30 FT-LB	1488 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 12" AWAY MINIMUM (i.e. - CORNER).
  SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE
  EDGE DISTANCES.
- C. TESTING AND SPECIAL INSPECTION OF EXPANSION ANCHORS SHALL
  BE PERFORMED BY AN APPROVED INDEPENDENT AGENCY
  EMPLOYED BY THE FACILITY OWNER PER CBC 1704A & 1910A.5
  AND CAC 7-149. ALL REPORTS SHALL BE SENT TO THE INSPECTOR
  OF RECORD, OWNER AND THE ARCHITECT OR ENGINEER IN
  RESPONSIBLE CHARGE.
  - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
  - (ii) ACCEPTANCE CRITERIA:
    - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
    - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED
       WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
  - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
  - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
  - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
  - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL (SLAB ON GRADE ONLY)



#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

### ADVANCED STERILIZATION PRODUCTS

**EVOTECH** 

DES. J. ROBERSON

11-2206 JOB NO.


6/7/22 DATE

SHEET

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB



(VALUES INCLUDE  $\Omega$ )

#### NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 2.30,  $\Delta_p = 1.0$ ,  $I_p = 1.5$ ,  $R_p = 1.5$ ,  $\Omega_0 = 2.0$ , z/h = 0)

HORIZONTAL FORCE (Eh) = 1.035 Wp

HORIZONTAL FORCE (Emh) = 2.07 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (Ev) = 0.46 Wp

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. EVOTECH IS USED WITH ASP'S WATER FILTRATION SYSTEM (OPM-0341-13)
- 5. SEE GENERAL NOTES: SHEET 1 AND 2



# EASE

#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

### ADVANCED STERILIZATION PRODUCTS

**EVOTECH** 

DES. J. ROBERSON

JOB NO. 11-2206


4

SHEET

DATE 6/7/22

ог 8 знеетз

SEISMIC SUPPORTS & ATTACHMENTS





#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

# ADVANCED STERILIZATION PRODUCTS DES. J. ROBERSON

**EVOTECH** 

11-2206

SHEET


6/7/22 DATE

JOB NO.

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK




NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 2.30,  $\Delta p = 1.0$ , p = 1.5, Rp = 1.5,  $\Omega_0 = 2.0$ , z/h < 1)

HORIZONTAL FORCE (En) = 2.76 Wp HORIZONTAL FORCE (Emh) = 5.52 Wp (FOR CONCRETE ANCHORAGE) VERTICAL FORCE (E<sub>V</sub>) = 0.46 W<sub>p</sub>

- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. EVOTECH IS USED WITH ASP'S WATER FILTRATION SYSTEM (OPM-0341-13)
- 5. SEE GENERAL NOTES; SHEET 1 AND 2



# EASE

#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

CONCRETE SLAB ON METAL DECK

## ADVANCED STERILIZATION PRODUCTS

**EVOTECH** 

DES. J. ROBERSON

JOB NO. 11-2206

DATE 6/7/22

SHEET

6

OF **8** SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

- CG WT= 578 LB MAX

UNIT FRAME MATERIAL

- AT CONNECTION

(13 GA, ASTM A1008 20 KSI MIN)

SEISMIC BRACKET (BY ASP)
(II GA, ASTM A240 TYPE 304, 30 KSI MIN)

W/ (4) - 5/16" \( \text{(A2-70) BOLTS} \)

TO UNIT FRAME (2 PLACES, 8 TOTAL) (TORQUE BOLTS TO 7 FT-LB)

(SEE BRACKET DETAIL "B")

USE 4- 3/8"φ (A36) THREADED RODS THRU FLOOR (2 PLACES, & TOTAL)

N.W. OR SAND L.W. CONC. (3000 PSI MIN.)

►VU

REFER TO MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL "A"

FOR ADDITIONAL DETAILS



Tu

# EASE

### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING**

www.EquipmentAnchorage.com

OF

### ADVANCED STERILIZATION PRODUCTS

**EVOTECH** 

DES. J. ROBERSON

**JOB NO.** 11-2206


DATE 6/7/22

NEET

8 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE DETAIL



L5 X 3 X I/4" X I'-2" (LLH) MIN
(A36) AT EACH ANCHOR
(EXTEND ANGLE TO
ADJACENT FLUTE WHEN
THREADED ROD OCCURS AT FLUTE)


FLANGE (TYP) AT CONDITIONS WHERE

NUT CANNOT BE PROVIDED AT TOP
SIDE OF STRUT, PROVIDE TAPPED
HOLE THROUGH STRUT FLANGE.
(TORQUE NUTS TO 13 FT-LB)

DATE: 06/14/2022

PVIA BUILDING COS

MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL



#### **EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING** www.EquipmentAnchorage.com SHEET DES. J. ROBERSON ADVANCED STERILIZATION PRODUCTS 11-2206 JOB NO. **EVOTECH** 6/7/22 DATE SHEETS BRACKET DETAILS SEISMIC SUPPORTS & ATTACHMENTS 24" **BRACKET** (II GA, ASTM A240 STAINLESS STEEL TYPE 304, Fy=30 KSI MIN) 0 0 (4) O.5Ι"Φ HOLES PLAN (4) 0.375" x 0.94" 19. 2.76" SLOTTED HOLES TYP TYP 2.6 Mohammad Aliaari 06/14/2022 FRONT SIDE