



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY
APPLICATION #: OPM-0441-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: [X] New [] Renewal [] Update to Pre-CBC 2013 OPA Number:

Manufacturer Information

Manufacturer: Handicare

Manufacturer's Technical Representative: Chris Huisman

Mailing Address: 10888 Metro Court, St Louis MO 63043

Telephone: 610-266-5260 x147 Email: Chris.huisman@handicare.com

Product Information

Product Name: Prism Medical Lifts

Product Type: Patient Lifts

Product Model Number: C-300, C-450, C-625, C-800, & C-1000

General Description: Patient lifts - Overhead supported & Wall post supported

Applicant Information

Applicant Company Name: Prism Medical

Contact Person: Chris Huisman

Mailing Address: 10888 Metro Court, St Louis MO 63043

Telephone: 610-266-5260 x147 Email: Chris.huisman@handicare.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2016.

Signature of Applicant: [Signature] Date: 8/18/2017

Title: Project manager Company Name: Handicare

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

Registered Design Professional Preparing Engineering Recommendations

Company Name: Degenkolb Engineers

Name: Robert Graff California License Number: SE 5113

Mailing Address: 235 Montgomery St, Suite 500, San Francisco CA 94104

Telephone: 415-392-6952 Email: rgraff@degenkolb.com

OSHPD Special Seismic Certification Preapproval (OSP)

- Special Seismic Certification is preapproved under OSP-
(Separate application for OSP is required)
- Special Seismic Certification is not preapproved

Certification Method(s)

- Testing in accordance with: ICC-ES AC156 FM 1950-16
- Other* (Please Specify): _____

OPM-0441-13

*Use of criteria other than those adopted by the California Building Standards Code, 2016 (CBSC 2016) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2016 may be used when approved by OSHPD prior to testing.

- Analysis
- Experience Data
- Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

DATE: 02/01/2018

List of Attachments Supporting the Manufacturer's Certification

- Test Report Drawings Calculations Manufacturer's Catalog
- Other(s) (Please Specify): Manufacturer product info provided within the calculation package.

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2016 & ALL PRE-2016 CODE BASED PROJECTS

Signature: *William Staehlin* Date: 02-01-2018

Print Name: William Staehlin

Title: SSE

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"



GENERAL NOTES

I. GENERAL

- THIS OSHPD PRE-APPROVAL OF MANUFACTURE'S CERTIFICATION (OPM) IS BASED ON THE CBC 2016. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2016.
- THIS PRE-APPROVAL IS VALID FOR THE EQUIPMENT DESCRIBED IN THESE DRAWINGS THROUGHOUT THE STATE OF CALIFORNIA, AND IS VALID FOR EQUIPMENT INSTALLED AT ANY HEIGHT WITHIN THE BUILDING.
- PROVIDE LABELING ON LIFTS WITH THE DESIGN LIFT CAPACITY IDENTIFIED.

II. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD

- VERIFY MATERIALS AND WORKMANSHIP TO CONFORM WITH THE 2016 EDITION OF THE CALIFORNIA BUILDING CODE AND THE REQUIREMENTS OF THIS PRE-APPROVAL DOCUMENT.
- VERIFY THE ADEQUACY OF THE EXISTING FRAMING TO SUPPORT THE LOADS INDICATED ON THIS SHEET, IN ADDITION TO ALL OTHER LOADS.
- VERIFY ANCHORS ARE ADEQUATE DISTANCES FROM OPENINGS AND EDGES OF SLABS.
- VERIFY ANCHORS ARE ADEQUATE DISTANCES FROM NEW OR EXISTING ANCHORS.
- DESIGN ANY SUPPLEMENTARY MEMBER AND THEIR ATTACHMENTS OTHER THAN THOSE DETAILED WITHIN THIS PRE-APPROVAL.
- VERIFY THE EQUIPMENTS WEIGHT, LOCATION. ANCHOR LOCATIONS AND ANCHOR DETAILS AGREE WITH THE INFORMATION SHOWN IN THIS PRE-APPROVAL.

III. STRUT FRAMING

- CHANNEL FRAMING COMPONENTS AND CONNECTORS MANUFACTURED BY MASON WEST CORPORATION.
- CHANNEL FRAMING TO CONFORM TO ASTM A1011 SS, GRADE 33.
- INSTALL BRACING WITH NO MORE THAN 5 DEGREE +/- PLAN DEVIATION.
- STRUT TYPE: SOLID SECTIONS ONLY.
- ALL STRUT NUTS AND BOLTS ARE 1/2" AND ARE TO BE TORQUED TO 50 FT-LBS UON.

IV. MECHANICAL ANCHORS

- WEDGE ANCHORS INTO CONCRETE: USE ZINC PLATED CARBON STEEL HILTI KB-TZ (ICC ESR-1917 ISSUED MAY 2015. INSTALL ANCHORS IN ACCORDANCE WITH ICC REPORT IF REINFORCEMENT IS ENCOUNTERED DURING DRILLING, ABANDON AND SHIFT THE HOLE LOCATION TO AVOID THE REINFORCEMENT. PROVIDE A MINIMUM OF 2 ANCHOR DIAMETERS OR 1 INCH, WHICHEVER IS LARGER, OF SOUND CONCRETE BETWEEN THE DOWEL AND THE ABANDONED HOLE. FILL THE ABANDONED HOLE WITH NON-SHRINK GROUT APPROVED BY THE ENGINEER OF RECORD. NOTIFY THE ENGINEER OF RECORD IF ANY REINFORCING IS DAMAGED.
- ANCHORS WILL BE PROOF-TESTED BY OWNER'S TESTING AND INSPECTION AGENCY. WITH A REPORT OF THE TEST RESULTS SUBMITTED TO OSHPD.
- IF ANY ANCHOR FAILS TESTING, REPLACE ANCHOR AND TEST ADDITIONAL ANCHORS OF THE SAME CATEGORY NOT PREVIOUSLY TESTED UNTIL TWENTY (20) CONSECUTIVE PASS, THEN RESUME INITIAL TESTING FREQUENCY.
- TEST ANCHORS NO SOONER THAN 24 HOURS AFTER INSTALLATION.
- TEST 50% WEDGE ANCHORS PER THE FOLLOWING METHOD:
 - TORQUE WRENCH METHOD: TEST ANCHORS TO THE TORQUE LOAD INDICATED IN THE TABLE BELOW WITHIN THE FOLLOWING LIMITS:
 - ONE-HALF TURN OF THE NUT.

WEDGE	
ANCHOR DIA. (IN)	TORQUE LOAD (FT-LBS)
3/8	25
1/2	40
5/8	60

V. ROUGH CARPENTRY

- FRAMING LUMBER: DOUGLAS FIR (COAST REGION) GRADED AND MARKED IN ACCORDANCE WITH THE STANDARD GRADING RULES NO. 17 OF THE WEST COAST LUMBER INSPECTION BUREAU (W.C.L.I.B.) OR WESTERN LUMBER GRADING RULES, OF THE WESTERN WOOD PRODUCTS ASSOCIATION (W.W.P.A.). USE LUMBER WITH A MINIMUM GRADE OF D.F. #2, U.O.N..
- ROUGH HARDWARE:
 - NAILS: COMMON WIRE NAILS, ASTM F1667, STANDARD LENGTHS U.O.N. USE HOT-DIPPED ZINC-COATED GALVANIZED NAILS WHEN PENETRATING PRESSURE TREATED OR FIRE-RETARDANT LUMBER.
 - BOLTS AND THREADED RODS: ASTM A307, SQUARE OR HEXAGONAL HEAD MACHINE BOLTS WITH ASTM A563 NUTS. USE MALLEABLE IRON WASHERS UNDER HEAD AND NUT WHEN IN CONTACT WITH WOOD.
 - CREWS: ASTM A307, ANSI/ASME STANDARD B18.6.1. USE CADMIUM-PLATED PAN OR ROUND HEADED SCREWS AT STEEL TO WOOD AND WOOD TO WOOD CONNECTIONS.
 - MISCELLANEOUS STEEL: ASTM A36.
- BOLT INSTALLATION:
 - DRILL BOLT HOLES A MAXIMUM OF 1/16 INCH LARGER IN DIAMETER THAN THE BOLT NOMINAL DIAMETER.

VI. STRUCTURAL STEEL

- STRUCTURAL STEEL TO CONFORM TO THE FOLLOWING UNLESS OTHERWISE NOTED:

SECTIONS	TYPE
COLD FORMED HOLLOW STRUCTURAL SECTION (HSS)	ASTM A500 GRADE B
MACHINE BOLTS	ASTM A307
THREADED AND HANGER ROD	ASTM A36
NUTS FOR BOLTS AND MACHINE BOLTS	ASTM A563
PLAIN WASHERS	ANSI B18.22.1

- HOT DIPPED GALVANIZED IN ACCORDANCE WITH ASTM A123 AND ASTM A153 STRUCTURAL STEEL AND FASTENERS THAT ARE PERMANENTLY EXPOSED TO WEATHER.

VII. STRUCTURAL TESTS, INSPECTIONS, AND OBSERVATIONS

- AN INDEPENDENT TESTING AGENCY AND SPECIAL INSPECTORS WILL BE RETAINED BY THE OWNER TO PERFORM THE FOLLOWING TESTS AND INSPECTION. PROVIDE ACCESS AND FURNISH SAMPLES TO THE AGENCY AS REQUIRED.
- THE FOLLOWING ITEMS REQUIRE TESTS AND INSPECTIONS IN ACCORDANCE WITH THE REQUIREMENTS OF THE CHAPTER "STRUCTURAL TESTS AND INSPECTIONS" OF THE CODE.
- MECHANICAL ANCHORS:
 - VERIFY TYPE OF ANCHOR, ANCHOR DIMENSIONS, CONCRETE TYPE AND COMPRESSIVE STRENGTH, PREDRILLED HOLE DIMENSIONS, ANCHOR SPACING, EDGE DISTANCE, SLAB THICKNESS AND ANCHOR EMBEDMENT.
 - PROOF-TEST AS INDICATED IN THE MECHANICAL ANCHORS SECTION OF THESE GENERAL NOTES.

VIII. DESIGN CRITERIA

- APPLICABLE CODE: 2016 CALIFORNIA BUILDING CODE.
- SEISMIC DESIGN:

SEISMIC FORCE	$F = 3.00 W_p$	$E_v = 0.50 W_p$	$R_p = 4.5$
WHERE:			$a_p = 2.5$
Sds = 250% G	WORST CASE ACCEL.		$\Omega = 2.0$
Ip = 1.5	FOR NON-ESSENTIAL EQUIP.		
Z/h = 1.0	FOR ANY FLOOR		
- CRANE LOADING PER AISC

TRANSVERSE LOADING	= 0.2 (DL+LL)
LONGITUDINAL LOADING	= 0.1 (DL+LL)

IX. HOW TO USE THIS PRE-APPROVAL

- REVIEW AND UNDERSTAND ALL GENERAL NOTES AND FIGURES BEFORE PROCEEDING.
- FOR THE SELECTED LIFT AND TRACK DETERMINE THE MAX HANGER SPACING AND BRACE SPACING FROM THE TABLES ON S2.
- BASED ON THE LIFT AND STRUCTURE TYPES SELECT A HANGER CONNECTION FROM THE TABLE ON S3.
- BASE ON THE LIFT AND STRUCTURE TYPES SELECT A BRACE CONNECTION FROM THE TABLE ON S5.
- DETERMINE THE MAXIMUM DEMANDS ON THE EXISTING STRUCTURE FROM THE NEW UNIT FROM THE TABLE ON THIS SHEET, AND VERIFY THE ADEQUACY OF THE EXISTING STRUCTURE WITH THE ENGINEER OF RECORD FOR THE BUILDING.

SHEET LIST

- S1 GENERAL NOTES
- S2 LAYOUT & TRACK SECTIONS
- S3 SECTION & TRAPEZE
- S4 HANGER CONNECTION DETAILS
- S5 HANGER CONNECTION DETAILS
- S6 BRACE CONNECTION DETAILS
- S7 BRACE CONNECTION DETAILS
- S8 WALL POST INSTALLATION

LOADS IMPOSED ON STRUCTURE

MAX LRFD LOADS W/OMEGA		
	T Hanger	P Brace
Lift	lbs	lbs
C300	1216	707
C450	1611	1471
C625	1977	1471
C800	2360	1506
C1000	2779	1506



DEGENKOLB ENGINEERS
 375 Beale Street, Suite 500
 San Francisco, CA 94105
 415.392.6952 Phone
 415.981.3157 Fax
 www.degenkolb.com



HANDICARE PATIENT LIFTS


Title: GENERAL NOTES

Sheet Number


Drawn	QL	Job number:	B5535011.00
Design:	RMG	Rev:	
Check:	RMG	Scale:	N.T.S.
Date	11/17/17		

S1


OF Sheets




SINGLE



MULTIPLE SECTION




CURVED TRACK



STANDARD TRACK

6063 T6 ALUMINIUM


LIFT	MAXIMUM TRACK SPANS (INCHES)		
	STD	SUPER	PLUS
C300	51	113	130
C450	40	96	108
C625	34	84	94
C800	30	74	83
C1000	26	66	74



SUPER TRACK

6005 T5 ALUMINIUM

LIFT	MAXIMUM TRACK CANTILEVER (INCHES)		
	STD	SUPER	PLUS
C300	12	28	32
C450	10	24	26
C625	8	21	23
C800	7	18	20
C1000	6	16	18



TRACK PLUS

6063 T6 ALUMINIUM

1 TRACK LAYOUT PLANS
N.T.S.

2 TRACK SECTIONS
N.T.S.

Degenkolb
DEGENKOLB ENGINEERS
375 Beale Street, Suite 500
San Francisco, CA 94105
415.392.6952 Phone
415.981.3157 Fax
www.degenkolb.com



HANDICARE PATIENT LIFTS


Title: LAYOUT & TRACK SECTIONS

Sheet Number

Drawn	Author	Job number: B5535011.00
Design:	Designer	Rev:
Check:	Checker	Scale: As indicated
Date	11/17/17	


S2

OF Sheets




MAXIMUM OFFSET	
LIFT	INCHES
C300	2
C450	1.625
C625	1.5
C800	1.375
C1000	1.125

1 SECTION @ BRACE
N.T.S.



LIFT	STRUT SPAN (INCHES)				
	STD	DOUBLE	DEEP	DOUBLE DEEP	HSS 4X4X3/16
C300	28	76	60	126	216
C450	20	58	49	114	183
C625	16	48	45	102	162
C800	13	39	37	93	146
C1000	11	32	31	87	134



2 TRAPEZE
N.T.S.



DEGENKOLB ENGINEERS
375 Beale Street, Suite 500
San Francisco, CA 94105
415.392.6952 Phone
415.981.3157 Fax
www.degenkolb.com



HANDICARE PATIENT LIFTS


Title: SECTION & TRAPEZE

Sheet Number


Drawn: QL Job number: B5535011.00
Design: RMG Rev:
Check: RMG Scale: As indicated
Date: 11/17/17

S3

OF Sheets




H-MDS# CONC. ON DECK SINGLE ANCHOR




H-MDD# CONC. ON DECK DOUBLE ANCHOR


HANGER CONNECTION TABLE																				
LIFT	CONCRETE FLAT SLAB						COMP METAL DECK				CONC JOIST			STEEL BEAM		STEEL JOIST	SOLID LUMBER		WOOD I-JOIST	
	H-CFS1	H-CFS2	H-CFD1	H-CFD1N	H-CFD2	H-CFD3	H-MDS1	H-MDS2	H-MDD1	H-MDD2	H-CPJ1	H-CPJ2	H-CPJ3	H-SB1	H-SB2	H-OWJ1	H-DL1	H-DL2	H-IJ1	H-IJ2
C300	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
C450	OK	OK	OK	OK	OK	OK	NO	OK	NO	OK	NO	OK	OK	OK	OK	OK	OK	OK	OK	OK
C625	NO	OK	NO	NO	OK	OK	NO	OK	NO	OK	NO	OK	OK	NO	OK	OK	OK	OK	OK	OK
C800	NO	NO	NO	NO	OK	OK	NO	NO	NO	OK	NO	OK	OK	NO	OK	OK	NO	OK	NO	OK
C1000	NO	NO	NO	NO	NO	OK	NO	NO	NO	OK	NO	NO	OK	NO	OK	OK	NO	OK	NO	OK




H-CFS# CONCRETE SLAB SINGLE ANCHOR




H-CFD# CONCRETE SLAB DOUBLE ANCHOR




HANDICARE PATIENT LIFTS
 Title: HANGER CONNECTIONS
 Sheet Number: S4
 Drawn: MSA Job number: B5535011.00
 Design: RMG Rev:
 Check: RMG Scale: N.T.S.
 Date: 11/17/17
 OF Sheets




H-DL# WOOD




SB2 ABOVE STEEL BEAM FLANGE




H-CPJ# CONCRETE JOIST




H-IJ# WOOD I - JOIST



H-OWJ1 OPEN WEB STEEL JOIST



SB1 BELOW STEEL BEAM



DEGENKOLB ENGINEERS
 375 Beale Street, Suite 500
 San Francisco, CA 94105
 415.392.6952 Phone
 415.981.3157 Fax
 www.degenkolb.com



HANDICARE PATIENT LIFTS

Title: HANGER CONNECTION DETAILS


Sheet Number

Drawn	QL	Job number:	B5535011.00
Design:	RMG	Rev:	
Check:	RMG	Scale:	As indicated
Date	11/17/17		


S5

OF Sheets


BRACE CONNECTION TABLE											
	CONCRETE FLAT SLAB		COMP METAL DECK		CONC JOIST	STEEL BEAM			STEEL JOIST	STEEL LUMBER	WOOD I-JOIST
LIFT	B-CFS1	B-CFD1	B-MDS1	B-MDD1	B-CPJ1	B-SB1	B-SB2	B-SB3	B-OWJ1	B-DL1	B-IJ1
C300	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
C450	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
C625	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
C800	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
C1000	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK




B-MDS-1 CONCRETE ON DECK SINGLE ANCHOR




B-CFS-1 CONCRETE SLAB SINGLE ANCHOR



H-CPJ 1 CONCRETE JOIST




B-MDD-1 CONCRETE ON DECK DOUBLE ANCHOR



B-CFD-1 CONCRETE SLAB DOUBLE ANCHOR

Degenkolb
DEGENKOLB ENGINEERS
 375 Beale Street, Suite 500
 San Francisco, CA 94105
 415.392.6952 Phone
 415.981.3157 Fax
 www.degenkolb.com




HANDICARE PATIENT LIFTS


Title: **BRACE CONNECTIONS** Sheet Number

Drawn	QL	Job number:	B5535011.00
Design:	RMG	Rev:	
Check:	RMG	Scale:	N.T.S.
Date	11/17/17		


S6
OF Sheets




B-DL1 WOOD




B-IJ1 WOOD I - JOIST




B-SB 3 BEAM CLAMP



B-SB 1 BOLT BEAM



B-OWJ 1 OPEN WEB JOIST



B-SB 2 WELD BEAM



DEGENKOLB ENGINEERS
 375 Beale Street, Suite 500
 San Francisco, CA 94105
 415.392.6952 Phone
 415.981.3157 Fax
 www.degenkolb.com



HANDICARE PATIENT LIFTS


Title:
 BRACE CONNECTIONS

Sheet Number


Drawn	Author	Job number:	B5535011.00
Design:	Designer	Rev:	
Check:	Checker	Scale:	As indicated
Date	11/17/17		

S7


OF Sheets




1 LAYOUT PLAN
N.T.S.



3 WALL POST AT METAL STUD WALL
N.T.S.



2 WALL POST AT CONC WALL
N.T.S.



Degenkolb
DEGENKOLB ENGINEERS
375 Beale Street, Suite 500
San Francisco, CA 94105
415.392.6952 Phone
415.981.3157 Fax
www.degenkolb.com



HANDICARE PATIENT LIFTS		Sheet Number
Title: WALL POST INSTALLATION		
Drawn: JQS	Job number: B5535011.00	S8
Design: RMG	Rev:	
Check: RMG	Scale: As indicated	
Date: 11/17/17		
		OF Sheets