

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL	OFFICE USE ONLY									
OF MANUFACTURER'S CERTIFICATION (OPM)	APPLICATION #:	OPM-0540-19								
	<u>-</u>									
OSHPD Preapproval of Manufacturer's Certification (OPM)										
Type: ⊠ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number:										
Manufacturer Information										
Manufacturer: Beckman Coulter, Inc.										
Manufacturer's Technical Representative: Matthew Sasse										
Mailing Address: 322 Lake Hazeltine Dr., Chaska, MN. 55318										
Telephone: On File Email: On File										
Product Information										
Product Name: UniCel Dxl 800/DXA Access Immunoassay System										
Product Type: Other Electrical and Mechanical Components 40-19	Z									
Product Model Number: Dxl 800/DXA	(((())))									
General Description: Blood Analyzer Blood Analyzer										
DATE: 04/20/2019	767									
F	7									
Applicant Information	Ŷ									
Applicant Company Name: EASE Co.	301									
Contact Person: Jonathan Roberson, S.E.										
Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709										
Telephone: (909) 606-7622 Email: J.Robe I hereby agree to reimburse the Office of Statewide Health F accordance with the California Administrative Code, 2019.	rson@EASECo.com Planning and Dev									

Signature of Applicant:

Principal Engineer

Title:

EASE Co.

Company Name:

Date: 6/12/19



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations									
Company Name: EASE Co.									
Name: Jonathan Roberson, S.E. California License Number: S4197									
Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709									
Telephone: 909-606-7622 Email: J.Roberson@EASECo.com									
OSHPD Special Seismic Certification Preapproval (OSP)									
 □ Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required) □ Special Seismic Certification is not preapproved 									
Certification Method(s)									
 ☐ Testing in accordance with: ☐ Other* (Please Specify): ☐ ICC-ES AC156 ☐ FM 1950-16 ☐ FM 1950-16 									
OPM-0540-19									
*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHPD prior to testing.									
Analysis DATE: 04/20/2019									
 Experience Data Combination of Testing, Analysis, and/or Experience Data (Please Specify): 									
List of Attachments Supporting the Manufacturer's Certification									
 ☐ Test Report ☐ Other(s) (Please Specify): 									
OFFICE USE ONLY - OSHPD APPROVAL VALID FOR CBC 2019 & ALL PRE-2019 CODE BASED PROJECTS									
Signature: Name And Amount Date: 4/20/2020									
Print Name: Haeseong Lim Title: Senior Structural Engineer									
Condition of Approval (if applicable):									



5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development

PREAPPROVAL OF MANUFACTURER'S CERTIFICATION OPM-0540-19

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER: BECKMAN COULTER
EQUIPMENT NAME: UniCel Dxl 800/DXA


Sheet: 1 of 9 Date: 12/5/19

GENERAL NOTES

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 1.80 & 2.20. SEE DETAIL FOR APPLICABILITY
- 4. FORCES PER ASCE 7-16 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE SDS = 1.80, a_p = 1.0, l_p = 1.5, R_p = 1.5, z/h = 0 AT CONCRETE SLAB. SEE FOLLOWING SHEETS FOR Ω_o WHERE SDS = 2.2, a_p = 1.0, l_p = 1.5, R_p = 1.0, z/h = 0 AT CONCRETE SLAB & $z/h \le 1$ AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_o
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 6. ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION AT OR BELOW GRADE. (i.e. z/h = 0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR. AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



www.EquipmentAnchorage.com

BECKMAN COULTER

UniCel Dxl 800/DXA

DES. J. ROBERSON

JOB NO. 11-1908

DATE 12/5/19

SHEET 2

9 SHEETS

10. EXPANSION ANCHORS:

A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.


Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension Test
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	2"	4"	3"	See Detail "A"	25 FT-LB	1825 lb
5/8"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	3.125"	12"	24"	5"	60 FT-LB	3135 lb
5/8"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	4"	12"	21"	6"	60 FT-LB	4540 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 21", 24" AWAY MINIMUM (i.e. CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING OF EXPANSION ANCHORS PER 2019 CBC, 1910A.5:

 TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR AND A REPORT OF THE TEST RESULTS SHALL BE 4 0 1 9 SUBMITTED TO OSHPD
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:

DATE: 04/20/2019

- DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO
 OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY
 TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER
 BECOMES LOOSE.
- TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
- (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
 - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
 - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL



www.EquipmentAnchorage.com

BECKMAN COULTER

UniCel Dxl 800/DXA

DES. J. ROBERSON

12/5/19

JOB NO. 11-1908


DATE

3

9 SHEETS


SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB



NOTES:

- 1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16 STRENGTH DESIGN IS USED. ($a_p = 1.0$, $|_p = 1.5$, $R_p = 1.5$, $\Omega_0 = 1.5$, z/h = 0)
- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.



www.EquipmentAnchorage.com

BECKMAN COULTER

UniCel DxI 800/DXA

DES. J. ROBERSON

ЈОВ NO. 11-1908

DATE 12/5/19

SHEET


4


9 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

MAX Sps ≤ I.80

CONCRETE SLAB





www.EquipmentAnchorage.com

BECKMAN COULTER

UniCel Dxl 800/DXA

DES. J. ROBERSON

JOB NO. 11-1908

DATE 12/5/19


5

9 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

1.80 < MAX Sps ≤ 2.20

CONCRETE SLAB





www.EquipmentAnchorage.com

BECKMAN COULTER

UniCel Dxl 800/DXA


DES. J. ROBERSON

JOB NO. 11-1908

DATE 12/5/19

8HEET

9 SHEETS





www.EquipmentAnchorage.com

BECKMAN COULTER

UniCel Dxl 800/DXA

DES. J. ROBERSON

JOB NO. 11-1908


DATE 12/5/19


7

9 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK





BECKMAN COULTER

UniCel DxI 800/DXA

www.EquipmentAnchorage.com

DES. J. ROBERSON

11-1908 JOB NO.

FLUTE DETAIL

12/5/19 DATE

SHEETS OF

SHEET

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE DETAILS 40" MIN **EQUIPMENT** EDGE DISTANCE N.W. OR SAND L.W. 3" CONC. (3000 PSI MIN.) Θ 12" (MIN) TYP I" MAX **OFFSET** Vu STRUT <u>w</u> ₹ MIN. 20 GA (TYP) STEEL W-DECK 4.5" FLUTE I" MIN USE 3/8" PHILTI KB-TZ TYP MIN Ł EXPANSION ANCHORS

HEX NUT TOP & BOT OF FLANGE

(TYP) AT CONDITIONS WHERE NUT CANNOT BE PROVIDED AT TOP SIDE

OF STRUT, PROVIDE TAPPED HOLE

L3 X 3 X I/4" X I'-2" MIN (A36) AT EACH ANCHOR (EXTEND ANGLE TO ADJACENT FLUTE WHEN THREADED ROD OCCURS AT FLUTE)

(MIN. EMBED. (hef) = 2")

(2 ANCHORS MIN PER STRUT)

DECK REQUIREMENTS AND STRUT DETAIL

ORWIA BUILDING

BY: HJHROUGH STRUT FLANGE.





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING www.EquipmentAnchorage.com SHEET DES. J. ROBERSON BECKMAN COULTER 11-1908 JOB NO. UniCel DxI 800/DXA 12/5/19 DATE SHEETS SEISMIC SUPPORTS & ATTACHMENTS BRACKET DETAILS 5" 0.50" 3" (2) HOLES 0.44"Φ 5.5" PART: C44090 COLD ROLLED ASTM A36 CS 3.5" SIDE **BRACKET** (I) HOLE 0.687" Φ FRONT BRACKET CHASSIS COMPONENT: SEISMIC MOUNT ·(0.313" THK MIN, A1008 SS, GRADE 30, Fy=29 KSI MIN) 3/8" \$\Phi BOLT (AISI 304, Fu=75 KSI MIN) PLAN **BRACKET** BRACKET DETAIL SIDE **ASSEMBLY** No. 4197

EXP. 6-30-2020