

# OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT **FACILITIES DEVELOPMENT DIVISION**

| APPLICATION FOR OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) |
|---|
|   |
| OSHPD Preapproval of Manufacturer's Certification (OPM)                 |

OFFICE USE ONLY

| MANUFACTURER 5 CERTIFICATION (OPM)   | APPLICATION #: OPM-0549  |
|--|--|
| OSHPD Preapproval of Manufacturer's Certification (OPM)  |  |
| Type: New X Renewal/Update   |  |
|  |  |
| Manufacturer Information   |  |
| Manufacturer: Abbott Laboratories  |  |
| Manufacturer's Technical Representative: Claudia Moreno  |  |
| Mailing Address: 1921 Hurd Drive, MS 2-33, Irving, TX 75038  |  |
| Telephone: (972) 518-7691 Email: Claudia.Moreno@a  | abbott.com   |
| FOR CODE CON   |  |
| Product Information OSHPD  | The state of the s |
| Product Name: Accelerator® a3600 Automation Track 3rd – Party Interfaces   | 7  |
| Product Type: Interfaces for 3rd -party unified automated diagnostic processing  | ng laboratory instruments  |
| Product Model Number: a3600  |  |
| General Description: modular system designed to automate pre/post-analytica allow for 3rd -party clinical analyzers (XN-9000 and Stag Accelerator® a3600 modular system.   |  |
| THE STATE OF THE S |  |
| Applicant Information  | D <sub>A</sub>   |
| Applicant Company Name:  |  |
| Contact Person: Dieter Siebald   |  |
| Mailing Address: 2495 Natomas Park Drive, Suite 650, Sacramento, CA 9583   | 3  |
| Telephone: (916) 920-2020 Email: dieters@cyseng.co   | m  |

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA - HEALTH AND HUMAN SERVICES AGENCY





Title: Structural Engineer/Project Manager



# OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

| Registered Design Professonal Preparing Engineering Recommendations  |  |  |  |  |  |  |
|--|--|--|--|--|--|--|
| Company Name: CYS STRUCTURAL ENGINEERS, INC.   |  |  |  |  |  |  |
| Name: Dieter Siebald California License Number: S4346  |  |  |  |  |  |  |
| Mailing Address: 2495 Natomas Park Drive, Suite 650, Sacramento, CA 95833  |  |  |  |  |  |  |
| Telephone: (916) 920-2020 Email: dieters@cyseng.com  |  |  |  |  |  |  |
|  |  |  |  |  |  |  |
| OSHPD Special Seismic Certification Preapproval (OSP)  |  |  |  |  |  |  |
| Special Seismic Certification is preapproved under OSP OSP Number:   |  |  |  |  |  |  |
|  |  |  |  |  |  |  |
| OR CODE  |  |  |  |  |  |  |
| Certification Method   |  |  |  |  |  |  |
| Testing in accordance with:   ICC-ES AC156   FM 1950-16  |  |  |  |  |  |  |
| Other(s) (Please Specify):   |  |  |  |  |  |  |
| *Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHPD prior to testing. |  |  |  |  |  |  |
| X Analysis BY: Jeffrey Kikumoto  |  |  |  |  |  |  |
| Experience Data  DATE: 11/20/2020  |  |  |  |  |  |  |
| Combination of Testing, Analysis, and/or Experience Data (Please Specify):   |  |  |  |  |  |  |
| CODE   |  |  |  |  |  |  |
| OSHPD Approval  BUILDING   |  |  |  |  |  |  |
| Date: 11/20/2020   |  |  |  |  |  |  |
| Name: Jeffrey Kikumoto Title: Senior Structural Engineer   |  |  |  |  |  |  |
| Condition of Approval (if applicable):   |  |  |  |  |  |  |

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"







GENERAL NOTES

THESE DRAWINGS ARE PREPARED FOR ABBOTT LABORATORIES, AN ILLINOIS CORPORATION, ABBOTT PARK, ILLINOIS.

1. THE CONTRACTOR AND THE INSPECTOR SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OSHPD WEBSITE.

2. THIS PRE-APPROVAL COVERS THE SUPPORTS AND ATTACHMENTS OF THE LABORATORY EQUIPMENT TO THE STRUCTURE.

SHEET TITLE: SHEET INDEX



ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



# CYS STRUCTURAL ENGINEERS, INC.

GENERAL INFORMATION

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

TEL (9 WW۱

|                |  | Date: | 6/30/2020 |
|----------------|--|-------|-----------|
|                |  | By:   | CYS       |
| (916) 920-2020 |  |       |           |
| ww.cyseng.com  |  | Page: | 1 of 25   |

Date

Rev Description

19090

Job No:


Doto.

PAGE


2. IT IS THE RESPONSIBILITY OF THE STRUCTURAL ENGINEER OF RECORD FOR A SITE SPECIFIC PROJECT TO VERIFY:

- A. THE ADEQUACY OF THE NEW OR EXISTING STRUCTURE TO RESIST THE FORCES AND WEIGHT SPECIFIED FOR EACH COMPONENT IN ADDITION TO ALL OTHER LOADS. PROVIDE AND DESIGN SUPPLEMENTARY MEMBERS AS REQUIRED.
- B. THAT THE ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS.
- C. THAT THE ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. THE SPACING SHOWN IN THE TEST VALUES TABLE ON THIS PAGE IS THE REQUIRED SPACING FROM ANCHORS OF OTHER DIAMETERS AND EMBEDMENTS WILL VARY.
- D. THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL.
- E. THAT THE ACTUAL EQUIPMENT'S WEIGHT, CENTER OF GRAVITY (CG) LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS, AND THE MATERIAL AND GAGE OF THE EQUIPMENT WHERE ATTACHMENTS ARE MADE, AGREE WITH THE INFORMATION SHOWN ON THE PRE—APPROVAL DOCUMENTS.
- 3A. EXPANSION ANCHORS INSTALLED IN NORMAL WEIGHT OR SAND-LIGHTWEIGHT CONCRETE SHALL BE STAINLESS STEEL HILTI KB-TZ EXPANSION ANCHORS COMPLYING WITH ICC-ES ESR-1917 REISSUED MAY 2019. ADHESIVE ANCHORS INSTALLED IN NORMAL WEIGHT CONCRETE SHALL BE ASTM F593 CW1 (316) INSTALLED USING HILTI HIT-RE 500 V3 ADHESIVE COMPLYING WITH ICC-ES ESR-3814 REVISED APRIL 2019.
- B. INSTALLATION: INSTALL THE POST-INSTALLED DRILLED-IN CONCRETE ANCHORS IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC EVALUATION REPORT FOR THE SPECIFIC ANCHOR AND THE PARAMETERS GIVEN IN THE TABLES ON THIS PAGE.
- C. TESTING:
  - JOB TESTING: FOR VERIFYING SATISFACTORY INSTALLATION WORKMANSHIP, PERFORM JOB SITE TESTING IN ACCORDANCE WITH THE TEST LOAD TABLE PROVIDED IN THIS DOCUMENT. TEST 50% OF THE INSTALLED ANCHORS. FOR TENSION TESTING, THE TEST LOAD MAY BE APPLIED BY ANY METHOD THAT WILL EFFECTIVELY MEASURE THE TENSION IN THE ANCHOR SUCH AS DIRECT PULL WITH A HYDRAULIC JACK OR CALIBRATED SPRING LOADING DEVICES. FOR TORQUE TESTING, THE TEST LOAD SHALL BE APPLIED WITH A CALIBRATED TORQUE WRENCH. ALL TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE INSPECTOR OF RECORD. IF ANY ANCHOR FAILS THE TEST, TEST ALL ANCHORS. THE TEST SHALL BE PERFORMED 24 HOURS OR MORE AFTER INSTALLATION. TESTING MAY BE DONE PRIOR TO EQUIPMENT INSTALLATION. ALSO REFER TO CBC 1910A.5.5 TESTS FOR POST—INSTALLED ANCHORS IN CONCRETE".
  - FAILURE/ACCEPTANCE CRITERIA: THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:
    - •• HYDRAULIC RAM METHOD: APPLY AND HOLD TEST LOAD FOR A MINIMUM OF 15 SECONDS.

      THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE APPLICABLE TEST LOAD WHERE WASHERS ARE USED. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER UNDER THE NUT BECOMES LOOSE OR BY A CONTINUOUS LOSS OF JACKING PRESSURE.
  - •• TORQUE WRENCH METHOD (EXPANSION ANCHORS ONLY): THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: ONE-HALF (½) TURN OF THE NUT.
- D. TEST VALUES: APPLY TEST LOADS TO ANCHORS WITHOUT REMOVING THE NUT.



ALL-THRD ROD



MECHANICAL ANCHOR

|   | POST-INSTALLED <u>ADHESIVE</u> ANCHOR SCHED |   |                          |      |                                      |   |                                 |                            |  |
|---|---|---|--------------------------|------|--------------------------------------|---|---------------------------------|----------------------------|--|
| CONDITION ANCHOR TYPE OF & DIA (INCH) ANCHORAGE |   |   | HOLE DIA<br>(INCH)<br>do |      | MIN CONC<br>THICKNESS<br>(INCH)<br>h | MIN CONC<br>EDGE<br>DISTANCE<br>(INCH)<br>c | MIN AB<br>SPCG<br>UNO<br>(INCH) | TENSION TEST<br>LOAD (LBS) |  |
|   | CASE 2                                      | ½"ø HILTI HAS-R<br>(ASTM F593 CW1 316 SS) | 0.5625                   | 2.75 | 4                                    | 12  | 6.75                            | 2590                       |  |

| 1               | POST-INSTALLED <u>MECHANICAL</u> ANCHOR SCHED |                       |                    |               |                       |                    |  |                       |                    |  |
|-----------------|---|-----------------------|--------------------|---------------|-----------------------|--------------------|--|-----------------------|--------------------|--|
| CONDITION       | ANCHOR<br>TYPE                                | INSTALLATION<br>EMBED | EFFECTIVE<br>EMBED | HOLE<br>DEPTH | MIN CONC<br>THICKNESS |                    | MIN AB<br>SPCG UNO                     | TEST I                | LOAD               |  |
| OF<br>ANCHORAGE | & DIA<br>(INCH)                               | (INCH)<br>hnom        | (INCH)<br>heff     | (INCH)<br>ho  | (INCH)<br>h           | DISTANCE<br>(INCH) | (INCH)                                 | TENSION<br>LOAD (LBS) | TORQUE<br>(FT-LBS) |  |
|                 | KB-TZ 304<br>SS 0.375"ø                       |                       | 2.00               | 2.625         | SEE<br>DTLS           | 12                 | 6.75<br>PARALLEL TO MTL<br>DECK FLUTES | 1350                  | 25                 |  |

No. S4346

No. S4346

No. S4346

Job No:

SHEET TITLE: GENERAL NOTES



ABBOTT LABORATORIES ACCELERATOR  $^{\circledR}$  a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



| CYS STRUCTURAL | ENGINEERS, | INC. |
|----------------|------------|------|

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-2020 www.cyseng.com


Date: 6/30/2020
By: CYS
Page: 2 of 25

Date

Rev Description

### **GENERAL NOTES CONTINUED:**

4. TWO (2) CONDITIONS OF SUPPORTS & ATTACHMENTS ARE SPECIFIED AND PRESENTED IN THIS PRE—APPROVAL:



<u>CASE 1:</u> SUPPORTS & ATTACHMENTS DTLS LOCATED AT UPPER FLRS ABV THE BASE OF A BLDG (z/h <= 1.0), IT IS ASSUMED THAT THE FLRS ARE BUILT OF A MIN  $3\frac{1}{4}$ " NWC OR SLWC TOPPING OVER MTL DECK (f'c = 3000 PSI, MIN).

<u>CASE 2:</u> SUPPORTS & ATTACHMENTS DTLS LOCATED AT OR BLW THE BASE OF A BLDG (z/h=0). THE FLRS ARE ASSUMED TO BE BUILT OF A MIN 4" NWC SLAB (f'c = 3000 PSI, MIN).

- 5. THIS PRE-APPROVAL MAY BE USED AT ANY GEOGRAPHICAL LOCATION IN THE STATE OF CALIFORNIA WHERE  $S_{DS}$  IS LESS THAN OR EQ TO 2.49.
- 6. COORDINATE THE ANCHOR BOLT LAYOUT WITH THE COMPONENT IN THE FIELD PRIOR TO SETTING ANCHOR BOLTS.
- 7. ANCHOR BRACKETS SHALL BE PAINTED WITH A RUST INHIBITIVE PRIMER FOLLOWED BY A COLOR COAT SELECTED BY THE HOSPITAL FACILITY OR MATCH THE COLOR OF THE BASE OF THE EQUIP IF A COLOR IS NOT SPECIFIED BY THE HOSPITAL.
- 8. FASTENERS AND ASSOCIATED HARDWARE SHALL BE FIELD PAINTED TO MATCH BRACKETS AFTER INSTALLATION IS COMPLETE.
- 10. UNIT COMPONENTS, TRACK SUPPORTS AND INTER-CONNECTION OF THE COMPONENTS TO BE PROVIDED BY ABBOTT. CONTRACTOR SHALL FURNISH AND INSTALL THE SEISMIC SUPPORTS AND ATTACHMENTS (INCLUDING SEISMIC BRACKETS, EXPANSION ANCHORS, THRU-BOLTS, STRUT PLATES BELOW SLABS, HIGH STRENGTH BOLTS, ETC.) IN CONJUNCTION WITH COMPONENT SETTING INSTRUCTIONS FROM ABBOTT FIELD INSTALLATION PERSONNEL.

  BY: Jeffrey
- 11. DRAWING SCALES ARE NOT PROVIDED. <u>DO NOT SCALE OFF OF THESE DRAWINGS</u>. THE INTENT OF THESE DRAWINGS IS TO SHOW HOW TO FABRICATE THE SEISMIC BRACKET TO ANCHOR THE EQUIPMENT SPECIFIED. THE REPRESENTATIONS OF THE EQUIPMENT ARE ONLY INTENDED TO SHOW THE COORDINATION WITH THE SEISMIC BRACKETS.
- 12. BOLTS THROUGH CONCRETE ON METAL DECK:
  - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER SNUG TIGHT (THE SNUG TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNO.
  - B. THRU BOLT HOLES SHALL BE  $\frac{1}{16}$ " LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE +  $\frac{1}{16}$ ").
  - C. THRU BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION & TESTING (THRU BOLTS WITH STEEL-TO-STEEL CONNECTION IN TENSION DO NOT REQUIRE TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.
- 13. TAKE CARE TO AVOID DAMAGING REBAR OR POST—TENSIONING TENDONS WHEN INSTALLING ANCHORS TO CONCRETE.

# **GENERAL NOTES CONTINUED:**

14. DRY BOLT AND NUT INSTALLATION TORQUES SHALL BE AS FOLLOWS:

| BOLT OR NUT DIA | TORQUE (FT-LBS) |
|-----------------|-----------------|
| M5              | 5               |
| M8              | 8               |
| M10             | 15              |
| M16             | 67              |
| M20             | 136             |
| 0.25"           | 10              |
| 0.50"           | 40              |

THESE VALUES DO NOT APPLY TO POST-INSTALLED CONCRETE ANCHORS.

15. FUTURE ALTERATIONS TO TRACK SYSTEMS, INCLUSIVE OF BUT NOT LIMITED TO TRACK AND TRACK COMPONENTS, TRACK MODULES AND OTHER INSTRUMENTS ATTACHED TO THE TRACK SYSTEM MUST BE REVIEWED BY OSHPD.

# **WELDING NOTES:**

- 1. WELDING OF <u>SEISMIC BRACKETS</u> SHALL BE PERFORMED BY CERTIFIED WELDERS USING E70XX ELECTRODES (UNO). THE USE OF E70—T4 WELDING WIRE IS NOT ALLOWED FOR ANY APPLICATION. WELDS SHALL BE IN CONFORMITY WITH THE STRUCTURAL WELDING CODE—STEEL OF THE AMERICAN WELDING SOCIETY (AWS D1.1—15). SUBMIT WELDING PROCEDURES AND SPECIFICATIONS TO OWNER'S TESTING LABORATORY FOR REVIEW AND APPROVAL PRIOR TO BEGINNING SEISMIC BRACKET FABRICATION.
- 2. WELD LENGTHS CALLED FOR ON PLANS ARE THE NET EFFECTIVE LENGTH REQUIRED. WHERE FILLET WELD SYMBOL IS GIVEN WITHOUT INDICATION OF SIZE, USE MINIMUM SIZE WELDS AS SPECIFIED IN AISC 360-16, SECTION J2.2b.

# **DESIGN CRITERIA**

DESIGN OF SUPPORTS & ATTACHMENTS FOR ALL EQUIP COMPONENTS IS PER 2019 CBC

ASCE 7-16 TABLE 13.6-1

OTHER MECHANICAL OR ELECTRICAL COMPONENTS (EXCEPT XN-9000 DOUBLE ROBOT CASE 2)

Y: Jeffrey Kikumo Sos = 2.5

 $I_p = 1.5$ 

 $a_p = 1.0$ 

 $R_p = 1.5$ 

 $\Omega_{\rm o} = 1.5$ 

WP AS NOTED ON COMPONENT BASE PLAN & ELEV. SEE PG 7.

SEISMIC LOADS FOR CASE 1 - UPPER FLRS ABV THE BASE, z/h <= 1.0 (LRFD)

 $F_p \ = \ 3.00 \ W_p$ 

 $F_v = 0.50 W_p$ 


SEISMIC LOADS FOR CASE 2 - SLAB AT OR BLW BASE, z/h = 0 (LRFD)

 $F_p = 1.12 W_p$ 

XN-9000 DOUBLE ROBOT:  $S_{DS} = 2.3$ 

 $F_v = 0.50 W_p$ 

 $F_p = 1.04 W_p F_v = 0.46 W_p$ 



## SHEET TITLE: GENERAL NOTES & DESIGN CRITERIA



ABBOTT LABORATORIES
ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



# CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO. CA 95833 TEL (916) 920-2020 www.cyseng.com

|          |  | Date: |
|----------|--|-------|
|          |  | By:   |
| 920-2020 |  |       |
| seng.com |  | Page: |
|          |  |       |

Rev Description

6/30/2020

CYS

3 of 25

|                             | NO  | NO. COMPONENT                       |           | LEVELING FOOT DIMS |         | CC     | LOCATIO | OP WT             | PG  |
|-----------------------------|-----|-------------------------------------|-----------|--------------------|---------|--------|---------|-------------------|-----|
|                             | NO. | COMPONENT                           | WIDTH     | DEPTH              | ×       | ÿ      | z       | (W <sub>p</sub> ) | -6  |
| A-R<br>CE                   |     | TYPICAL TRACK MODULE <sup>(2)</sup> | 90" MAX   | 11.50"             | MIDSPAN | 5.75"  | 26.9"   | 165#              | 7,8 |
| STAGO STA—I<br>INTERFACE    | 1.  | ROBOT MODULE                        | N/A       | N/A                | 7.32"   | 3.27"  | 2.56"   | 33#               | 11  |
| STA(                        | 2.  | ELECTRICAL CONTROL PANEL            | 5.3"      | 14.7"              | 7.35"   | 2.65"  | 10.5"   | 11#               | 12  |
| 00                          | 3.  | OUTPUT LOADER MODULE                | 24.80"    | 18.11"             | 13.34"  | 5.89"  | 24.06"  | 62#               | 13  |
| XN-90<br>FACE               | 4.  | RETURN LANE <sup>(4)</sup>          | 46.5" MAX | 5.7 <b>"</b>       | MIDSPAN | 2.56"  | 27.36"  | 9 PLF             | 14  |
| SYSMEX XN-9000<br>INTERFACE | 5.  | INPUT LOADER MODULE                 | 22.83"    | 18.11"             | 7.72"   | 11.33" | 30.03"  | 84#               | 15  |
| SY                          | 6.  | DOUBLE ROBOT MODULE                 | 20.87"    | 41.77"             | 7.13"   | 20.79" | 29.33"  | 374#              | 17  |


### **NOTES:**

SEE COMPONENT DIMENSION DTL FOR MORE INFO.

COMPONENT OP WT (Wp) INCLUDES WT OF TWO TOP TRACKS, BOTT TRACK, & TWO SUPPORTS FRAMES AS SHOWN ON PG 7.

THE DIMS FURNISHED ABY ARE NOT INTENDED TO BE USED TO LAYOUT THE COMPONENTS. THIS INFORMATION IS BEING FURNISHED FOR USE BY THE SEOR.

FRAMES FOR RETURN LANE WEIGH APPROXIMATELY 7 LBS EA. WT IN TABLE IS FOR RETURN LANE ONLY.



**ABBREVIATIONS:** ΑT AB ANCHOR BOLT **ABV** ABOVE ADJ **ADJACENT AISC** AMERICAN INSTITUTE FOR STEEL CONSTRUCTION **ALUM ALUMINUM ASCE** AMERICAN SOCIETY OF CIVIL ENGINEERS ASD ALLOWABLE STRENGTH DESIGN **ASTM** AMERICAN SOCIETY FOR TESTING & MATERIALS AWS AMERICAN WELDING SOCIETY **BLDG BUILDING** BLW BELOW **BOTT** воттом CBC CALIFORNIA BUILDING CODE CG CENTER OF GRAVITY CJP COMPLETE JOINT PENETRATION **CLR** CLEAR OR CLEARANCE **CLSE** CALIFORNIA LICENSED STRUCTURAL ENGINEER CENTERLINE CONC CONCRETE CONN CONNECTION

COORD COORDINATE CTR CENTER

DBL DOUBLE DIA (ø) DIAMETER DIM DIMENSION DTL DETAIL DWG DRAWING (E) **EXISTING CONDITION** 

OPM-

EA **EACH** EACH END EE ELEV **ELEVATION EQUAL** EQ **EQUIP EQUIPMENT** 

MINIMUM ULTIMATE COMPRESSIVE STRENGTH DATE: 11/20/2020 OF CONCRETE

> FINISHED FLOOR FLG FLANGE FLR **FLOOR FRMG** FRAMING FT (') FOOT/FEET  $F_{v}$

SPECIFIED YIELD STRENGTH OF REINFORCING, PSI OR SPECIFIED MINIMUM YIELD STRESS OF STEEL, KSI

GA GAUGE **GALVANIZED** GALV GENERAL CONTRACTOR GC ICC INTERNATIONAL CODE COUNCIL

IN (") INCH INFO INFORMATION

KSI KIPS PER SQUARE INCH LENGTH

LBS **POUNDS** 

**LRFD** LOAD & RESISTANCE FACTOR DESIGN **LFRS** LATERAL FORCE RESISTING SYSTEM

MAX MAXIMUM MFR

MIN

MTL

**MANUFACTURER** MINIMUM **METAL** 

NO. (#) NUMBER OR POUNDS NTS NOT TO SCALE

NEAR SIDE & FAR SIDE NS&FS NORMAL WEIGHT CONCRETE NWC **OPG OPENING** 

OP WT OPERATING WEIGHT **OSHPD** OFFICE OF STATEWIDE HEALTH PLANNING

& DEVELOPMENT

PG(S) PAGE(S) PLATE

PLF POUNDS PER LINEAL FOOT PSI POUNDS PER SQUARE INCH

**RADIUS** R REQ **REQUIRED** SC SLIP-CRITICAL **SCHED SCHEDULE** 

**SEOR** STRUCTURAL ENGINEER OF RECORD

SIM SIMILAR

SLWC SAND LIGHT WEIGHT CONCRETE

SOG SLAB ON GRADE **SPCG SPACING** 

SQ SQUARE SS STAINLESS STEEL

STIFF **STIFFENER** STL STEEL TOP & BOTTOM T&B

**TEMP TEMPORARY** 

THRD THREAD OR THREADED

MAXIMUM TENSION DUE TO SEISMIC FORCE Tmax TOC TOP OF CONCRETE

ANCHORAGE TENSION REACTION Tu

DUE TO SEISMIC FORCE TYP **TYPICAL** 

UNO UNLESS NOTED OTHERWISE ANCHORAGE SHEAR REAC

**VERT VERTICAL** 

MAXIMUM SHEAR DUE TO SEISMIC FORCE Vmax ANCHORAGE SHEAR REACTION ۷u

DUE TO SEISMIC FORCE

W/ WITH

COMPONENT SELF-WEIGHT

WEIGHT

Rev Description

Job No:

SHEET TITLE: COMPONENT DIMENSIONS SCHEDULE & ABBREVIATIONS



ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES




| CYS STRUCTURAL | . ENGINEERS, | INC. |
|----------------|--------------|------|

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

TEL (916) 920www.cysend

|        |  |  | Date: | 6/30/2020 |
|--------|--|--|-------|-----------|
|        |  |  | Ву:   | CYS       |
| )-2020 |  |  |       |           |
| g.com  |  |  | Page: | 4 of 25   |
|        |  |  | •     |           |

Date



ANCHOR OPTIONS

19090 Date Job No: 6/30/2020 Date: CYS By:

Page:

Rev Description


SHEET TITLE: TYPICAL STRUT DETAILS ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES

SACRAMENTO, CA 95833

|   | CYS STRUCTURAL ENGINEER            | SINC |                |  |
|---|------------------------------------|------|----------------|--|
| > | 2495 NATOMAS PARK DRIVE. SUITE 650 | •    | (916) 920-2020 |  |
|   | SACRAMENTO, CA 95833               |      | ww.cyseng.com  |  |

Abbott

5 of 25



Date Job No: 19090

Date: 6/30/2020

By: CYS

Page:

Rev Description

SHEET TITLE: TYPICAL STRUT DETAILS


ABBOTT LABORATORIES

ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



|   | CYS STRUCTURAL ENGINEERS           | S. INC. |                |  |
|---|------------------------------------|---------|----------------|--|
| > | 2495 NATOMAS PARK DRIVE, SUITE 650 | •       | (916) 920-2020 |  |
|   | SACRAMENTO, CA 95833               | W       | ww.cyseng.com  |  |

6 of 25



SHEET TITLE: MODULE & INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT) TYPICAL TRACK MODULE

Abbott

ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



# CYS STRUCTURAL ENGINEERS, INC.


2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

TEL (916) 920-2020 www.cyseng.com

Job No: 6/30/2020 Date: CYS By: lPage: 7 of 25

Date

Rev Description



SHEET TITLE: MODULE & INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT) TYPICAL TRACK MODULE



ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



| CYS STRUCTURAL ENGINEERS, 2495 NATOMAS PARK DRIVE, SUITE 650 | NC. |      |
|--|-----|------|
| 2495 NATOMAS PARK DRIVE, SUITE 650                           | TEL | (916 |


2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

|   | (916) 920-2020 |  |
|---|----------------|--|
|   | www.cyseng.com |  |
| - |                |  |

Rev Description

Date

Job No:




- 1. COMPONENTS ARE SHOWN ONLY TO DESCRIBE THE INTERFACE SYSTEMS. NUMBER & LOCATION OF COMPONENTS MAY VARY BASED ON SITE SPECIFIC NEEDS.
- 2. COMPONENT SUB-ASSEMBLY CONNS SHALL BE SUPPLIED & PERFORMED BY ABBOTT, NOT BY THE GENERAL CONTRACTOR.
- 3. SEE PG 4 FOR COMPONENT UNIT WTS & DIMS.
- 4. CLADDING PANELS NOT SHOWN FOR CLARITY.

|   | Rev | Description | Date | Job No: | 19090     |
|---|-----|-------------|------|---------|-----------|
|   |     |             |      | Date:   | 6/30/2020 |
| CYS STRUCTURAL ENGINEERS, INC.                        |     |             |      | Ву:     | CYS       |
| 2495 NATOMAS PARK DRIVE, SUITE 650 TEL (916) 920-2020 |     |             |      |         |           |
| SACRAMENTO, CA 95833 www.cyseng.com                   | 1   |             |      | Page:   | 9 of 25   |



ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES Abbott


XN-9000 INTERFACE W/ INPUT LOADER LOAD BAY MODULE



NOTES:

- 1. COMPONENTS ARE SHOWN ONLY TO DESCRIBE THE INTERFACE SYSTEMS. NUMBER & LOCATION OF COMPONENTS MAY VARY BASED ON SITE SPECIFIC NEEDS.
- 2. COMPONENT SUB-ASSEMBLY CONNS SHALL BE SUPPLIED & PERFORMED BY ABBOTT, NOT BY THE GENERAL CONTRACTOR.
- 3. SEE PG 4 FOR COMPONENT UNIT WTS & DIMS.
- 4. CLADDING PANELS NOT SHOWN FOR CLARITY.

|                     |                                    |     |                | Rev | Description | Date | Job No: | 19090     |
|---------------------|------------------------------------|-----|----------------|-----|-------------|------|---------|-----------|
|                     |                                    |     |                |     |             |      | Date:   | 6/30/2020 |
|                     | CYS STRUCTURAL ENGINEERS, I        | NC. |                |     |             |      | By:     | CYS       |
| <b>&gt;&gt;&gt;</b> | 2495 NATOMAS PARK DRIVE, SUITE 650 | TEL | (916) 920-2020 |     |             |      |         |           |
|                     | SACRAMENTO, CA 95833               | V   | vww.cyseng.com |     |             |      | Page:   | 10 of 25  |
|                     |                                    |     |                |     |             |      | -       |           |







ABBOTT LABORATORIES

ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES


45x45 OPEN PROFILE




45x45 CLOSED PROFILE




25x40x40 FRMG BRACKET DETAIL



40x40x40 FRMG BRACKET DETAIL



TYP ALUM FRMG BRACKET DETAIL



# EXTRUDED ALUM PROFILE SECTION PROPERTIES:

| PROFILE            | A<br>(MM <sup>2</sup> ) | Sx (MM <sup>3</sup> ) | Sy<br>(MM <sup>3</sup> ) | lx<br>(MM <sup>4</sup> ) | y<br>(MM <sup>4</sup> ) |
|--------------------|-------------------------|-----------------------|--------------------------|--------------------------|-------------------------|
| 45x45 OPEN FRAME   | 721                     | 5.88                  | 5.88                     | 13.25                    | 13.25                   |
| 32x45 OPEN FRAME   | 533                     | 3.55                  | 4.28                     | 5.80                     | 8.70                    |
| 45x45 CLOSED FRAME | 585                     | 5.20                  | 5.20 <sub>I</sub>        | 11.73                    | 11.73/                  |


# TYP LEVELING LEG CONN DETAIL



20mmø AISI 304 SS PIN 20 (ASTM A240 Fy MIN= 30 KSI) W/ TAPPED HOLE FOR LEVELING LEG

> ALUM CLADDING PANEL, WHERE OCCURS

> > 45x45 OPEN PROFILE EXTRUDED ALUM FRMG, SEE PG 11



10mmø MIN LEVELING LEG W/FOOT. SEISMIC BRACKET PER PG 22 (NOT SHOWN FOR CLARITY)

# PROFESS /ON ALL TO SHOW THE OF CALIFORNIA

Joh No:

# 32x45 OPEN PROFILE

# EXTRUDED ALUM PROFILE CROSS SECTIONS

### NOTES:

- 1. ALL DIMS SHOWN W/ OUT UNITS ARE IN MM.
- 2. ATTACHMENTS TO FRMG USING M8 UNI 5931-8.8 BOLTS W/ SPRING NUTS, TYP. Fy= 96 KSI MIN. INSTALL TORQUE 8 FT-LBS. ATTACHMENT BY ABBOTT.
- 3. MATERIAL: ALUM EN AW-6060-T5 Fy=  $14.5 \text{ KSI } (100 \text{ Nmm}^2) \text{ MIN}$

SHEET TITLE: INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT)
XN-9000 SUBCOMPONENTS TYPICAL CONNECTIONS



ABBOTT LABORATORIES ACCELERATOR  $^{\circledR}$  a3600 AUTOMATION TRACK 3RD PARTY INTERFACES




# CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-202 www.cyseng.com

|        | I/G A | Description | Date | 300 NO. | 19090     |
|--------|-------|-------------|------|---------|-----------|
|        |       |             |      | Date:   | 6/30/2020 |
|        |       |             |      | By:     | CYS       |
| )-2020 |       |             |      |         |           |
| g.com  |       |             |      | Page:   | 11 of 25  |


Date


Pay Description



NOTES:

- 1. COMPONENTS ARE SHOWN ONLY TO DESCRIBE THE INTERFACE SYSTEMS. NUMBER & LOCATION OF COMPONENTS MAY VARY BASED ON SITE SPECIFIC NEEDS.
- 2. COMPONENT SUB-ASSEMBLY CONNS SHALL BE SUPPLIED & PERFORMED BY ABBOTT, NOT BY THE GENERAL CONTRACTOR.
- 3. SEE PG 4 FOR COMPONENT UNIT WTS & DIMS.
- 4. CLADDING COVER PANELS NOT SHOWN FOR CLARITY.








ABBOTT LABORATORIES ACCELERATOR  $^{\circledR}$  a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



|    | CYS STRUCTURAL ENGINEERS,          | INC |
|----|------------------------------------|-----|
| >> | CAOS NATONAC DADIC DDIVE CULTE CEC | -   |

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO. CA 95833 TEL (916) 920-2020 www.cyseng.com

Rev Description






ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



|                                    |      |                | Rev | Description | Date | Job No: | 19090     |
|------------------------------------|------|----------------|-----|-------------|------|---------|-----------|
|                                    |      |                |     |             |      | Date:   | 6/30/2020 |
| CYS STRUCTURAL ENGINEERS,          | INC. |                |     |             |      | By:     | CYS       |
| 2495 NATOMAS PARK DRIVE. SUITE 650 | TEL  | (916) 920-2020 |     |             |      |         |           |
| SACRAMENTO, CA 95833               | V    | vww.cyseng.com |     |             |      | Page:   | 13 of 25  |

Abbott



SHEET TITLE: MODULE & INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT)

COMPONENT 2: OUTPUT LOADER MODULE

ABBOTT LABORATORIES

ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES

CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650

TEL (916) 920-2020


11/20/2020

Abbott

Page:

www.cyseng.com

14 of 25



SHEET TITLE: MODULE & INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT)
COMPONENT 3: RETURN LANE

ABBOTT LABORATORIES
ACCELERATOR <sup>®</sup> a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



| CYS STRUCTURAL ENGINEERS,         | INC |
|-----------------------------------|-----|
| 2495 NATOMAS PARK DRIVE SHITE 650 | TF  |

, TEL

TEL (916) 920-2020 www.cyseng.com

SACRAMENTO, CA 95833

ROBOT SUPPORT ARM (ALUM UNI 3051 Fy= 80Mpa MIN) ATTACHED TO INPUT LOADER MODULÉ FRAME W/ M8 CLASS 8.8 BOLTS & SPRING NUTS, TYP OF 4

TYP TRACK MODULE TOP RAIL (NOT SHOWN), SEÈ <u>PGS 7 & 8</u> TYP ES OF PANEL. SEE PG 11 ELECTRICAL CONTROL PANEL. ATTACH W/ 2- M8 SCREWS MIN W/ SPRING NUTS. SIM TO BY: Jeffrey Kikumot PANEL SHOWN ON PG 13 TYP TRACK MODULE -(ATTACHMENT BY ABBOTT) BOTT RAIL (NOT SHOWN), DATE: 11 SEE PGS 7 & 8 ALUM SUPPORT FRAME, SEE PLAN VIEW ALUM FRMG BRACKET, SEE PG 11. ATTACH TO TOP & BOTT RAIL OF TYP TRACK (NOT SHOWN)

POLYCARBONATE COVERS

Job No:

SEE SPLICE DTL

ON PG 15 FOR

RETURN LANE PER <u>PG 15</u>

CONN

**ROBOT ARM ATTACHMENT** 

SHEET TITLE: MODULE & INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT) COMPONENT 4: INPUT LOADER MODULE & ELECTRICAL CONTROL PANEL



ABBOTT LABORATORIES ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES

CYS STRUCTURAL ENGINEERS, INC.


2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833


TEL (916) 920-2020 www.cyseng.com

Rev Description

6/30/2020 Date: CYS By: lPage: 16 of 25

Date







SHEET TITLE: MODULE & INTERFACE SUB-ASSEMBLY DETAILS (CONNECTION BY ABBOTT)


COMPONENT 5: DOUBLE ROBOT



ABBOTT LABORATORIES
ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



| ) |                                       |    |                | Rev | Description | Date | Job No: | 19090     |
|---|---------------------------------------|----|----------------|-----|-------------|------|---------|-----------|
|   |                                       |    |                |     |             |      | Date:   | 6/30/2020 |
|   | CYS STRUCTURAL ENGINEERS, INC         | ). |                |     |             |      | Ву:     | CYS       |
|   | 2495 NATOMAS PARK DRIVE, SUITE 650 TE |    | (916) 920-2020 |     |             |      |         |           |
|   | SACRAMENTO, CA 95833                  | w  | ww.cyseng.com  |     |             |      | Page:   | 17 of 25  |



Abbott

**BASE PLAN & ELEVATIONS** 

ABBOTT LABORATORIES

ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES

2495 NATOMAS PARK DRIVE, SUITE 650

CYS STRUCTURAL ENGINEERS, INC.

TEL (916) 920-2020

www.cyseng.com

Date:


Page:

By:

6/30/2020

CYS

18 of 25



- 1. ECCENTRICITY & PRYING ACTION CONSIDERED BASED ON THE SEISMIC BRACKET CONFIGURATION SHOWN.
- 2. INCLUDES OVERSTRENGTH FACTOR  $(\Omega_0)$ .
- 3. OVERSTRENGTH FACTOR  $(\Omega)$  MUST BE APPLIED FOR ANCHORAGE TO CONC.

No. S4346

No. S4346

No. S4346

Job No:

SHEET TITLE: XN - 9000 SPUR BASE PLAN W/ INPUT LOADER MODULE

Abbott

ABBOTT LABORATORIES ACCELERATOR  $^{(8)}$  a3600 AUTOMATION TRACK 3RD PARTY INTERFACES




CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-2020 www.cyseng.com

|      |  | Date: | 6/30/2020 |
|------|--|-------|-----------|
|      |  | By:   | CYS       |
| 2020 |  |       |           |
| .com |  | Page: | 19 of 25  |

Date

Rev Description



- 1. ECCENTRICITY & PRYING ACTION CONSIDERED BASED ON THE SEISMIC BRACKET CONFIGURATION SHOWN.
- 2. INCLUDES OVERSTRENGTH FACTOR  $(\Omega_0)$ .
- 3. OVERSTRENGTH FACTOR  $(\Omega)$  MUST BE APPLIED FOR ANCHORAGE TO CONC.

PROFESS /ON TO TO SHARE TO SHA

Job No:

SHEET TITLE: XN - 9000 SPUR BASE PLAN W/ DOUBLE ROBOT

Abbott

ABBOTT LABORATORIES
ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



| CYS STRUCTURAL ENGINEERS, INC. |
|--------------------------------|
|--------------------------------|

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833


TEL (916) 920-2020 www.cyseng.com

|        |  | Date: | 6/30/2020 |
|--------|--|-------|-----------|
|        |  | By:   | CYS       |
| 0-2020 |  |       |           |
| ng.com |  | Page: | 20 of 25  |


Date

Rev Description

# XN-9000 DOUBLE ROBOT SEISMIC BRACKET DETAIL:



"H" VARIES TO ACCOMMODATE VERT ADJUSTMENT OF THE COMPONENT FOR LEVELING PURPOSES AS MEASURED FROM THE FLR TO THE BOTT OF THE COMPONENT PER THE CASE 1 & CASE 2 FLR TO COMPONENT CLEARANCES AS FOLLOWS:



|         |          | 14             |      |  |
|---------|----------|----------------|------|--|
| BRACKET | "H" (IN) | CLEARANCE (IN) |      |  |
| DRACKET |          | MIN            | MAX  |  |
| Α       | 4.15     | 4.15           | 4.55 |  |
| В       | 3.75     | 3.75           | 4.15 |  |
| С       | 3.35     | 3.35           | 3.75 |  |
| D       | 2.95     | 2.95           | 3.35 |  |
| E       | 2.55     | 2.55           | 2.95 |  |
| F       | 2.15     | 2.15           | 2.55 |  |
| G       | 1.75     | 0,0            | 2.15 |  |
|         |          |                |      |  |

<u>NOTES:</u>

- . FOR INFO NOT SHOWN OR NOTED, SEE TYP SEISMIC BRACKET DETAIL.
- 2. FOR CASE 1 & CASE 2 ANCHORAGE TO FLR, SEE PG 25

SHEET TITLE: TYPICAL TRACK MODULE & XN-9000 DOUBLE ROBOT SEISMIC BRACKET FABRICATION DETAIL

ABBOTT LABORATORIES ACCELERATOR  $^{(8)}$  a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



# CYS STRUCTURAL ENGINEERS, INC.

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833

PM - 0549

rey Kikumoto

TEL (916) 920-2020 www.cyseng.com


|     |  | Date: | 6/30/2020 |
|-----|--|-------|-----------|
|     |  | By:   | CYS       |
| 020 |  |       |           |
| om  |  | Page: | 21 of 25  |
|     |  |       |           |


Date

Rev Description


Abbott

Job No:





OVERSTRENGTH FACTOR (126) MUST BE APPLIED TO Vu FOR ANCHORAGE TO CONC.



INCLUDES OVERSTRENGTH FACTOR  $(\Omega_0)$ 

USE HARDENED ASTM F436

IS TOO SMALL FOR A NUT

MAX

4.75"

Date

±Vu

CLR


WASHERS TO SHIM GAP IF GAP

EMBED (hef)

(E) 4" MIN SOG

(f'c = 3000 PSI)


Ζ



STRUT PL PER PGS 5&6

CASE 1 - SUSPENDED FLR W/ THRU BOLTS

ACCELERATOR ® a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



SEE <u>PG 2</u> FOR TYPE & INSTALLATION REQUIREMENTS

CASE 2 - SLAB ON GRADE

(SLAB AT OR BLW GRADE)

SHEET TITLE: TYPICAL TRACK MODULE
SUPPORT & ATTACHMENT DETAILS

ABBOTT LABORATORIES

CYS STRUCTURAL ENGINEERS, INC.

 CYS STRUCTURAL ENGINEERS, INC.
 By:
 CYS

 2495 NATOMAS PARK DRIVE, SUITE 650
 TEL (916) 920-2020 www.cyseng.com
 Page:
 23 of 25


Abbott

**OCCURS** 

EDGE DISTANCE

/2020

Job No:



Abbott


ABBOTT LABORATORIES ACCELERATOR  $^{\circledR}$  a3600 AUTOMATION TRACK 3RD PARTY INTERFACES



| CYS STRUCTURAL | ENGINEERS, | INC. |
|----------------|------------|------|

2495 NATOMAS PARK DRIVE, SUITE 650 SACRAMENTO, CA 95833 TEL (916) 920-2020 www.cyseng.com

Date: 6/30/202
By: CYS
Page: 24 of 25



SACRAMENTO, CA 95833

25 of 25

lPage:

www.cyseng.com