

# OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

altiful I hals								
APPLICATION FOR OSHPD F	OFFICE USE ONLY APPLICATION #: OPM-0593							
MANUFACTURER'S CERTIFI								
OSHPD Preapproval of Manufacturer's Certification (OPM)								
Type: X New Renewal/Update								
Manufacturer Information								
Manufacturer: Carestream Health								
Manufacturer's Technical Representative:	Christopher Kralles							
Mailing Address: 1049 Ridge Road West, Rochester, NY 14615								
Telephone: (800) 328-2910	Email: christopher.kralles@carestreamhealth.com							
	FOR CODE CO.							
Product Information	OSHPD							
Product Name: DRX-Excel Plus Generator	Cabinet	Y						
Product Type: Generator	OPM-0593	CH						
Product Model Number: N/A	BY: Haeseong Lim							
General Description: Provides Power to M	ledical Equipment							
A	DATE: 10/06/2021	502						
Applicant Information	TO A							
Applicant Company Name: EASE LLC.	(0)	<b>Y</b>						
Contact Person: Tiffany Tonn	BUILDING							
Mailing Address: 1515 FAIRVIEW AVE, S	TE 205, MISSOULA, MT 59801							

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

Telephone: (406) 541-3273 Email: tiffany@easeco.com





Title:



# OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations								
Company Name: EASE								
Name: Jonathan Roberson	California License Number: S4197							
Mailing Address: 5877 Pine Ave., Suite 210, Chino Hills, CA 91709								
Telephone: (951) 295-1892	Email: jon@EASECo.com							
OSHPD Special Seismic Certification Preap	proval (OSP)							
X Special Seismic Certification is preapproved ur	nder OSP OSP Number: OSP-0675							
	OR CODE							
Certification Method	EON CON							
Testing in accordance with: ICC-ES AC156	6 FM 1950-16							
Other(s) (Please Specify):	THE STATE OF THE S							
and attachments are not permitted. For distribution	lifornia Building Standards Code, 2019 (CBSC 2019) for component supports system, interior partition wall, and suspended ceiling seismic bracings, test may be used when approved by OSHPD prior to testing.							
X Analysis	Y: Haeseong Lim							
C C	ATE: 10/06/2021							
Combination of Testing, Analysis, and/or Experience Data (Please Specify):								
RAZ	CODE							
OSHPD Approval	BUILDING							
Date: 10/6/2021								
Name: Haeseong Lim	Title: Senior Structural Engineer							
Condition of Approval (if applicable):								

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"







5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development PREAPPROVAL OF MANUFACTURER'S CERTIFICATION OPM-0593

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER: CARESTREAM HEALTH. INC

Sheet: 1 of 9

EQUIPMENT NAME: DRX EXCEL PLUS DIGITAL R/F SYSTEM GENERATOR CABINET Date: 8/25/21

#### **GENERAL NOTES**

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE Sos IS NOT GREATER THAN 2.30.
- 4. FORCES PER ASCE 7-16 SECTION 13.3.1, EQUATIONS 13.3-1 13.3-2 & 13.3-3, WHERE SDS = 2.30,  $\mathbf{a}_p$  = 2.5,  $\mathbf{l}_p$  = 1.5,  $\mathbf{R}_p$  = 6.0,  $\mathbf{z}/h$  = 0 AT CONCRETE SLAB &  $\mathbf{z}/h$  < 1 AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR  $\Omega_o$
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 6. ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION AT OR BELOW GRADE. (i.e. z/h = 0)

#### 9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev ) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

# DRX EXCEL PLUS DIGITAL R/F SYSTEM GENERATOR CABINET

DES. J. ROBERSON

**JOB NO.** 11-2002

DATE 8/25/21

SHEET

9 SHEETS

#### 10. EXPANSION ANCHORS:


A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension Test
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	6.75"	12"	See Detail "A"	30 FT-LB	N/A
1/2"	Normal Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	4"	14"	4"	50 FT-LB	1626 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES 14" AWAY MINIMUM (i.e. - CORNER).
  SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE
  EDGE DISTANCES.
- C. TESTING AND SPECIAL INSPECTION OF EXPANSION ANCHORS SHALL
  BE PERFORMED BY AN APPROVED INDEPENDENT AGENCY
  EMPLOYED BY THE FACILITY OWNER PER CBC 1704A & 1910A.5
  AND CAC 7-149. ALL REPORTS SHALL BE SENT TO THE INSPECTOR
  OF RECORD, OWNER AND THE ARCHITECT OR ENGINEER IN M 0 5 9 3
  RESPONSIBLE CHARGE.
  - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, im DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
  - (ii) ACCEPTANCE CRITERIA:
    - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
    - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
  - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK

A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.

- B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE
- C. (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
- D. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL



www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX EXCEL PLUS DIGITAL R/F SYSTEM **GENERATOR CABINET**


DES. J. ROBERSON

11-2002 JOB NO.

8/25/21 DATE

SHEET

SHEETS



#### NOTES:

#### 1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16

STRENGTH DESIGN IS USED. (SDS = 2.30 ap = 2.5, p = 15, p = 6.0, p = 6.0

HORIZONTAL FORCE (En) = 1.035 Wp HORIZONTAL FORCE (Emh) = 2.07 Wp (FOR CONCRETE ANCHORAGE) VERTICAL FORCE (E<sub>V</sub>) = 0.46 W<sub>p</sub>

- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. GENERATOR CABINET (UUT-2001-4) HAS OBTAINED SPECIAL SEISMIC CERTIFICATION REFER TO OSP-0675.
- 5. SEE GENERAL NOTES: SHEETS 1 AND 2.



www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX EXCEL PLUS DIGITAL R/F SYSTEM GENERATOR CABINET

DES. J. ROBERSON

JOB NO. 11-2002

DATE 8/25/21

SHEET

4

OF 9 SHEETS

CONCRETE SLAB

SEISMIC SUPPORTS & ATTACHMENTS FRONT

23.12" 11.84" II.28" 9.3" 8.75" HOLE NOT USED (SIM ON OTHER BRACKET) 0 O <u>4</u> 0 O 0 Ö C.G. WT. = 227 LB. HAREARY Lim (Y = 19.5")PLAN AT BASE

BRACKET MATERIAL
(0.13" THK, ASTM 1008, Fy= 29 KSI MIN)
W/(3)- 1/2"\$\phi\$ HILTI KB-TZ2
EXPANSION ANCHORS
(MIN. EMBED. (het) = 2")
(EA BRACKET, 6 TOTAL)
(SEE BRACKET DETAIL "B")

BRACKET MATERIAL
(O.16" THK, ASTM 1008, Fy= 29 KSI MIN)
W (2)- MI2 (CL 8.8) BOLTS
W STANDARD WASHERS
BRACKET TO UNIT MATERIAL (II GA ALUM)
(EA SIDE, 4 TOTAL)
(SEE BRACKET DETAILS "C")

No. 4197

EXP. 6-30-2022

S. 8/25/21

PUCTURE

OF CALKER

OF CALKE

www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX EXCEL PLUS DIGITAL R/F SYSTEM GENERATOR CABINET

DES. J. ROBERSON


JOB NO. 11-2002

DATE 8/25/21

SHEET

5

f 9 sheets




#### NOTES:

#### 1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (Sps = 2.30, 2p = 2.5, 2p = 1.5, 2p = 6.0, 2p = 2.0, 2p = 2.0, 2p = 1.5, 2p = 6.0, 2p = 2.0, 2p

HORIZONTAL FORCE (En) = 1.73 Wp HORIZONTAL FORCE (Emh) = 3.46 Wp (FOR CONCRETE ANCHORAGE) VERTICAL FORCE (Ev) = 0.46 Wp

- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE ESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. GENERATOR CABINET (UUT-2001-4) HAS OBTAINED SPECIAL SEISMIC CERTIFICATION REFER TO OSP-0675.
- 5. SEE GENERAL NOTES: SHEETS 1 AND 2.



www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX EXCEL PLUS DIGITAL R/F SYSTEM GENERATOR CABINET

DES. J. ROBERSON

JOB NO. 11-2002

DATE 8/25/21

SHEET

6

9 <sub>SHEETS</sub>

SEISMIC SUPPORTS & ATTACHMENTS CONCRETE SLAB ON METAL DECK FRONT 23.12" BRACKET MATERIAL 11.84" 11.28" (0.13" THK, ASTM 1008, Fy= 29 KSI MIN) W/(3)- 1/2" (A36) THREADED RODS THRU FLOOR (EA BRACKET, 6 TOTAL) 9.3" 8.75" (SEE BRACKET DETAIL "B") HOLE NOT USED (SIM ON OTHER BRACKET) O ₽ 0 O 0 BRACKET MATERIAL (0.16" THK, ASTM 1008, Fy= 29 KSI MIN) 0 W (2)- MI2 (CL 8.8) BOLTS W/ STANDARD WASHERS BRACKET TO UNIT MATERIAL (II GA ALUM) (EA SIDE, 4 TOTAL) C.G. WT. = 227 LB (SEE BRACKET DETAILS "C") REAR Lim (Y = 19.5")DAPILEA NO AGE / BRASE



www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX EXCEL PLUS DIGITAL R/F SYSTEM **GENERATOR CABINET**

DES. J. ROBERSON

11-2002 JOB NO.

8/25/21 DATE

MIN 20 GA

STEEL W-DECK

SHEET

SHEETS

(TYP)

FLUTE

Œ

FLUTE DETAIL

SEISMIC SUPPORTS & ATTACHMENTS CONCRETE DETAIL 18" MIN **EQUIPMENT** EDGE DISTANCE N.W. OR SAND L.W.  $\Theta$ CONC. (3000 PSI MIN.) 12" (MIN) TYP I" MAX **OFFSET** Vu STRUT

USE 3/8" PHILTI KB-TZ2 EXPANSION ANCHORS (MIN. EMBED. (het) = 2") (2 ANCHORS MIN PER STRUT)

 $\bar{\sigma}$ 

L3 X 3 X I/4" X I'-2" MIN (A36) AT EACH ANCHOR (EXTEND ANGLE TO ADJACENT FLUTE WHEN THREADED ROD OCCURS AT FLUTE)

HEX NUT TOP & BOT OF FLANGE (TYP) AT CONDITIONS WHERE NUT CANNOT BE PROVIDED AT TOP SIDE OF STRUT, PROVIDE TAPPED HOLE

I" MIN

TYP

4.5"


MIN

BY: HJHROUGH STRUT FLANGE.

MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL

OPWIA BUILDING





www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX-EXCEL PLUS DIGITAL R/F SYSTEM GENERATOR CABINET

DES. J. ROBERSON

**JOB NO.** 11-2002


DATE 8/25/21

8HEET

<u> 9 sнеетв</u>

SEISMIC SUPPORTS & ATTACHMENTS

BRACKET DETAILS





www.EquipmentAnchorage.com

# CARESTREAM HEALTH, INC

## DRX EXCEL PLUS DIGITAL R/F SYSTEM **GENERATOR CABINET**

PLAN

I

DES. J. ROBERSON

11-2002 JOB NO.

8/25/21 DATE

SHEET

SHEETS BRACKET DETAILS

SEISMIC ANCHORAGE

BRACKET MATERIAL 23.23" (0.16" THK, ASTM 1008, Fy= 29 KSI MIN) 3.15" TYP 0.16" 1.10" 1.52" 0 0 (2)- 0.53" HOLES 5.31 (BRACKET TO BRACKET CONN.) 0 C 1.1" ACCESS HOLE 5.50" (ALIGNS WITH I ½"中 ACCES HOLES) SIDE

> (2)- 0.53" HOLES (CONNECTION TO UNIT)

FRONT

