

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT **FACILITIES DEVELOPMENT DIVISION**

APPLICATION FOR OSHPD PREAPPROVAL OF	OFFICE USE ONLY APPLICATION #: OPM-0613							
MANUFACTURER'S CERTIFICATION (OPM)								
OSHPD Preapproval of Manufacturer's Certification (OPM)								
Type: X New Renewal/Update								
Manufacturer Information								
Manufacturer: AGFA Healthcare Corporation								
Manufacturer's Technical Representative: Brendan Conley								
Mailing Address: 10 South Academy Street, Greenville, SC 29601								
elephone: (214) 548-7214 Email: brendan.conley@agfa.com								
EOR CODE COM								
Product Information OSHPD								
Product Name: DR 800 Non-Tilting Wall Stand	Y							
Product Type: Other Mechanical & Electrical Components	CH							
Product Model Number: N/A BY: Keh-Shin Chi								
General Description: Wall Stand for the DR 800 System								
DATE: 10/08/2021	70							
	N							
Applicant Information								
Applicant Company Name: EASE LLC.								
Dr. C.								

Contact Person: Tiffany Tonn

Mailing Address: 1515 FAIRVIEW AVE, STE 205, MISSOULA, MT 59801

Telephone: (406) 541-3273

Email: tiffany@easeco.com

Title: Office Manager

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA - HEALTH AND HUMAN SERVICES AGENCY







OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations							
Company Name: EASE							
Name: Jonathan Roberson	California License Number: S4197						
Mailing Address: 5877 Pine Ave., Su	ite 210, Chino Hills, CA 91709						
Telephone: (951) 295-1892	Email: jon@EASECo.com						
OSHPD Special Seismic Certific	ation Preapproval (OSP)						
X Special Seismic Certification is p	reapproved under OSP OSP Number: 0526						
	OR CODE						
Certification Method	EOIL						
Testing in accordance with:	CC-ES AC156 FM 1950-16						
Other(s) (Please Specify):							
and attachments are not permitted. For	ted by the California Building Standards Code, 2019 (CBSC 2019) for component supports or distribution system, interior partition wall, and suspended ceiling seismic bracings, test e CBSC 2019 may be used when approved by OSHPD prior to testing.						
X Analysis	BY: Keh-Shin Chi						
Experience Data	DATE: 10/08/2021						
Combination of Testing, Analysis	, and/or Experience Data (Please Specify):						
	CODE CODE						
OSHPD Approval	BUILDING						
Date: 10/8/2021							
Name: Keh-Shin Chi	Title: Senior Structural Engineer						
Condition of Approval (if applicable):							

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"









5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development
PREAPPROVAL OF MANUFACTURER'S CERTIFICATION
OPM-0613

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER: AGFA HE

AGFA HEALTHCARE CORP

AFD

Sheet: <u>1 of 7</u>

EQUIPMENT NAME:

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER

Date: 7/26/21

GENERAL NOTES

- THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 2.0.
- 4. FORCES PER ASCE 7-16 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE SDS = 2.00, \mathbf{a}_p = 1.0, \mathbf{l}_p = 1.5, \mathbf{r}_p = 1.5, \mathbf{z}/h = 0 AT CONCRETE SLAB & \mathbf{z}/h < 1 AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_o
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION AT OR BELOW GRADE. (i.e. z/h = 0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS,
 MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE
 PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



www.EquipmentAnchorage.com

AGFA HEALTHCARE CORP

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER

DES. J. ROBERSON

7/26/21

ЈОВ NO. 11-2011

DATE

SHEET 2


7 SHEETS

10. EXPANSION ANCHORS:

A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	1.5"	6"	12"	3.25"	30 FT-LB	773 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 12" AWAY MINIMUM (i.e. CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING AND SPECIAL INSPECTION OF EXPANSION ANCHORS SHALL BE PERFORMED BY AN APPROVED INDEPENDENT AGENCY EMPLOYED BY THE FACILITY OWNER PER CBC 1704A & 1910A.5 AND CAC 7-149. ALL REPORTS SHALL BE SENT TO THE INSPECTOR OF RECORD, OWNER AND THE ARCHITECT OR ENGINEER IN RESPONSIBLE CHARGE.
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST AT LEAST 50% OF THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:
 - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO
 OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY
 TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER
 BECOMES LOOSE.
 - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
 - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.



TYPICAL CONCRETE EDGE DETAIL



www.EquipmentAnchorage.com

ÁGFA HEALTHCARE CORP

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER


DES. J. ROBERSON

11-2011 JOB NO.

7/26/21 DATE

SHEET

SHEETS



NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16

STRENGTH DESIGN IS USED. (SDS = 2.00, Δp = 1.0, |p| = 1.5, Rp = 1.5, Ω_0 = 2.0, z/h < 1)

HORIZONTAL FORCE (En) = 2.40 Wp HORIZONTAL FORCE (Emh) = 4.80 Wp (FOR CONCRETE ANCHORAGE) VERTICAL FORCE (E_V) = 0.40 W_D

- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. DR 800 WALL STAND HAS OBTAINED SPECIAL SEISMIC CERTIFICATION REFER TO OSP-0526.
- SEE GENERAL NOTES: SHEETS 1 AND 2.



www.EquipmentAnchorage.com

OF

AGFA HEALTHCARE CORP

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER

DES. J. ROBERSON


JOB NO. 11-2011

DATE 7/26/21

SHEET

4

SHEETS





(f'c = 3000 PSI MIN.)

www.EquipmentAnchorage.com

ÁGFA HEALTHCARE CORP

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER

DES. J. ROBERSON

JOB NO. 11-2011


DATE 7/26/21

SHEET 5

OF 7 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

BRACKET DETAIL



NOTE: VERIFIED THROUGH SHAKE TABLE TESTING PER 0SP-0526



www.EquipmentAnchorage.com

ÁGFA HEALTHCARE CORP

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER

DES. J. ROBERSON

11-2011 JOB NO.


7/26/21 DATE

SHEET

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS


BRACKET DETAIL



ATE

ADAPTER BUILDING

NOTE: VERIFIED THROUGH SHAKE TABLE TESTING PER OSP-0526



www.EquipmentAnchorage.com

AGFA HEALTHCARE CORP

DR 800 NON-TILTING WALLSTAND WITH AND WITHOUT SPACER

DES. J. ROBERSON

JOB NO. 11-2011


DATE 7/26/21

7

_F 7 _{энеетэ}

SEISMIC SUPPORTS & ATTACHMENTS

BRACKET DETAIL



BUILDING

NOTE: VERIFIED THROUGH SHAKE TABLE TESTING PER OSP-0526

