

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

WWW DAY				
APPLICATION FOR OSHPD PREAPPROVAL OF	OFFICE USE ONLY			
MANUFACTURER'S CERTIFICATION (OPM)	APPLICATION #: OPM-0618			
OSHPD Preapproval of Manufacturer's Certification (OPM)				
Type: X New Renewal/Update				
Manufacturer Information				
Manufacturer: Steris Corporation				
Manufacturer's Technical Representative: Kyle Irish				
Mailing Address: 5900 Heisley Rd., Mentor, OH 44060				
Telephone: (440) 392-7292 Email: Kyle_Irish@steris.com	ı			
FOR CODE COM				
Product Information OSHPD				
Product Name: AMSCO 600-37 Steam Sterilizer	T.			
Product Type: Other Electrical & Mechanical Components	CH C			
Product Model Number: 42 SD, 54 SD, 66 SD, 42 DD, 54 DD, 66 DD				
General Description: Sterilzer used to sanitize medical instruments, gowns etc				
DATE: 09/13/2021	201			
Applicant Information				
Applicant Company Name: EASE LLC.				

Contact Person: Tiffany Tonn

Mailing Address: 1515 FAIRVIEW AVE, STE 205, MISSOULA, MT 59801

Telephone: (406) 541-3273 Email: tiffany@easeco.com

Title: Office Manager

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY







OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations							
Company Name: EASE							
Name: Jonathan Roberson California License Number: S4197							
Mailing Address: 5877 Pine Ave., Suite 210, Chino Hills, CA 91709							
Telephone: (951) 295-1892 Email: jon@EASECo.com							
OSHPD Special Seismic Certification Preapproval (OSP)							
Special Seismic Certification is preapproved under OSP OSP Number:							
OR CODE							
Certification Method							
Testing in accordance with:							
Other(s) (Please Specify):							
*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHPD prior to testing.							
X Analysis BY: William Staehlin							
Experience Data DATE: 09/13/2021							
Combination of Testing, Analysis, and/or Experience Data (Please Specify):							
CODE CODE							
OSHPD Approval BUILDING							
Date: 9/13/2021							
Name: William Staehlin Title: Senior Structural Engineer							
Condition of Approval (if applicable):							

"Access to Safe. Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"









5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development PREAPPROVAL OF MANUFACTURER'S CERTIFICATION

OPM-0618

THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE

MANUFACTURER: STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

Sheet: 1 of 12 Date: 7/15/21

GENERAL NOTES

EQUIPMENT NAME:

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2019 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 1.60 & 2.20. SEE DETAIL FOR APPLICABILITY
- 4. FORCES PER ASCE 7-16 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3,

WHERE SDS = 1.60, a_p = 1.0, I_p = 1.5, R_p = 1.5, z/h = 0 AT CONCRETE SLAB

WHERE SDS = 2.20, \mathbf{a}_p = 1.0, \mathbf{I}_p = 1.5, \mathbf{R}_p = 1.5, \mathbf{z}/h = 0 AT CONCRETE SLAB & $\mathbf{z}/h \le 1$ AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω_0

- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 6. ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB DETAIL VALID FOR DEMANDS SHOWN AT OR BELOW GRADE. (i.e. z/h = 0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS TH REQUIREMENTS OF THE APPLICABLE ICC ESR REPORT AND THIS OPM.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



www.EquipmentAnchorage.com

STERIS CORPORATION

DES. J. ROBERSON

DATE

JOB NO. 14-2014

7/15/21

.. 1

40

SHEETS

SHEET

AMSCO 600-37 STEAM STERILIZER


10. EXPANSION ANCHORS:

A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension Test
5/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	4"	12"	18"	See Detail "A"	40 FT-LB	N/A
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	2"	8"	16"	3.25"	30 FT-LB	1190 lb
3/4"	Normal Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	3.75"	4"	30"	6"	110 FT-LB	3450 lb
3/4"	Normal Weight	3000	Hilti Kwik Bolt TZ2	ESR-4266	4.75"	4"	30"	6"	110 FT-LB	4646 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 16 & 28" AWAY MINIMUM (i.e. CORNER).

 SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING AND SPECIAL INSPECTION OF EXPANSION ANCHORS SHALL BE PERFORMED BY AN APPROVED INDEPENDENT AGENCY EMPLOYED BY THE FACILITY OWNER PER CBC 1704A & 1910A.5 AND CAC 7-149. ALL REPORTS SHALL BE SENT TO THE INSPECTOR OF RECORD, OWNER AND THE ARCHITECT OR ENGINEER IN RESPONSIBLE CHARGE.
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:
 - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
 - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
 - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
 - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
 - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL



7/15/21

DATE

www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER


DES. J. ROBERSON

JOB NO. 14-2014

OF 12 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB



NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16 STRENGTH DESIGN IS USED. ($a_p = 1.0$, $|_p = 1.5$, $R_p = 1.5$, $\Omega_0 = 2.0$, z/h = 0)

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.



www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

DES. J. ROBERSON

14-2014 JOB NO.

7/15/21 DATE


SHEET

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

MAX Sps ≤ 1.60

CONCRETE SLAB



(2)- 1/2" ASTM A354 (GR 5) BOLTS & NUTS (4 PLACES, & TOTAL) (BRACKET TO UNIT)

(2)- 3/4" PHLTI KB-TZ2 EXPANSION ANCHORS (MIN. EMBED. (hef) = 3.75") (4 PLACES, 8 TOTAL) (SEE DETAIL "B")

WEIGHT (LB)	DIM. "A" (IN.)	(LB/BOLT)	/** Yu (LB/BOLT)
3800	54,0	1920	889
4200	66.0	2058	9830
4700	78.0	2253	1100
3800	54.5	1917	889
4200	66.5	2056	983
4700	78.5	2252	1100
	(LB) 3800 4200 4700 3800 4200	(LB) (IN.) 3800 54.0 4200 66.0 4700 78.0 3800 54.5 4200 66.5	(LB) (IN.) (LB/BOLT) 3800 54.0 J920 4200 66.0 2058 4100 78.0 2253 3800 54.5 J917 4200 66.5 2056

^{**} VALUES INCLUDE Ω.



www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

DES. J. ROBERSON

14-2014 JOB NO.

7/15/21 DATE


SHEET

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

1.60 < MAX Sos ≤ 2.20

CONCRETE SLAB



(2)- 1/2" ASTM A354 (GR 5) BOLTS & NUTS (4 PLACES, & TOTAL) (BRACKET TO UNIT)

(2)- 3/4" PHLTI KB-TZ2 EXPANSION ANCHORS (MIN. EMBED. (hef) = 4.75") (4 PLACES, 8 TOTAL) (SEE DETAIL "B")

MODEL	WEIGHT	DIM. "A"	1 0	** Yu
TIODEE .	(LB)	(IN.)	(LB/BOLT)	(LB/BOLT)
AMSCO 600-37 42 SD	3800	54.0	1920	889
AMSCO 600-37 54 SD	4200	66.0	2058	9830
AMSCO 600-37 66 SD	4700	78.0	2253	1100
AMSCO 600-37 42 DD	3800	54.5	1917	889
AMSCO 600-37 54 DD	4200	66.5	2056	983
AMSCO 600-37 66 DD	4700	78.5	2252	1100

^{**} VALUES INCLUDE Ω_{o}



www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

DES. J. ROBERSON

JOB NO. 14-2014


DATE 7/15/21

SHEET

12 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK



NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 2.20, $\Delta p = 1.0$, Ip = 1.5, Rp = 1.5, $\Omega_0 = 2.0$, z/h < 1)

HORIZONTAL FORCE (En) = 2.64 Wp

HORIZONTAL FORCE (Emh) = 5.28 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (Ev) = 0.44 Wp

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THESE CALCULATIONS ENCOMPASS ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL VERIFY ALL CONDITIONS AND PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.



www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

DES. J. ROBERSON

JOB NO. 14-2014


DATE 7/15/21

SHEET 7

12 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK



MODEL	WEIGHT (LB)	NDIM, "A"n (IN.)	**S (Tajeh 1 : (LB/BOLT)	** Yu (LB/BOLT)
AMSCO 600-37 42 SD	3800	54.0	3614	1556 0
AMSCO 600-37 54 SD	4200	66.0	3882	1720
AMSCO 600-37 66 SD	4700	78.0	4257	1925/
AMSCO 600-37 42 DD	3800	54.5	3609	1556
AMSCO 600-37 54 DD	4200	66.5	3878	1720
AMSCO 600-37 66 DD	4700	78.5	4254	1925

^{**} VALUES INCLUDE Ω₀



www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

DES. J. ROBERSON

JOB NO. 14-2014


DATE 7/15/21

SHEET

OF 12 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB



FRONT ELEVATION
(STEAM GENERATOR)

SIDE ELEVATION (STEAM GENERATOR)

> AT $z/h \le 10$ Tu = 438 LB/BOLT (MAX) Vu = 462 LB/BOLT (MAX) (VALUES INCLUDE Ω_0)

NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 1.60, Δp = 1.0, lp = 1.5, Rp = 1.5, Ω_0 = 2.0)

HORIZONTAL FORCE (En) = 1.92 Wp

HORIZONTAL FORCE (Emh) = 3.84 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (Ev) = 0.32 Wp

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.



www.EquipmentAnchorage.com

STERIS CORPORATION

AMSCO 600-37 STEAM STERILIZER

DES. J. ROBERSON

JOB NO. 14-2014


DATE 7/15/21

SHEET

F 12 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK



NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 2.20, Δp = 1.0, Ip = 1.5, Rp = 1.5, Ω_0 = 2.0)

HORIZONTAL FORCE (Eh) = 2.64 Wp HORIZONTAL FORCE (Emh) = 5.28 Wp (FOR CONCRETE ANCHORAGE) VERTICAL FORCE (Ev) = 0.44 Wp

- 2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
- 3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.
- 4. SEE GENERAL NOTES: SHEETS 1 AND 2.



www.EquipmentAnchorage.com

STERIS CORPORATION

JOB NO. 14-2014

DES. J. ROBERSON

7/15/21


10

AMSCO 600-37 STEAM STERILIZER

JOB NO. 14-2014

DATE

12 SHEETS



NOTES:

1. FORCES ARE DETERMINED PER 2019 CALIFORNIA BUILDING CODE AND ASCE 7-16.

STRENGTH DESIGN IS USED. (SDS = 2.20, Δp = 1.0, |p| = 1.5, Rp = 1.5, Ω_0 = 2.0 z/h < 1)

HORIZONTAL FORCE (En) = 2.64 Wp

HORIZONTAL FORCE (Emh) = 5.28 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (Ev) = 0.44 Wp

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER LOADS THAT MAY BE PRESENT.



www.EquipmentAnchorage.com

STERIS CORPORATION

DES. J. ROBERSON

7/15/21


11

AMSCO 600-37 STEAM STERILIZER

JOB NO. 14-2014


DATE

F 12 SHEETS



DECK REQUIREMENTS AND STRUT DETAIL

OPNIA BUILDING



www.EquipmentAnchorage.com

STERIS CORPORATION

DES. J. ROBERSON
JOB NO. 14-2014

12

AMSCO 600-37 STEAM STERILIZER

DATE 7/15/21

12 SHEETS

BRACKET DETAILS SEISMIC SUPPORTS & ATTACHMENTS 5.625" BRACKET 1.25" (3/8" THK A36) 1/4 35. INSIDE BOTH SIDES 1.0625 FRONT SIDE 1.31" 3.5" TYP 5.625" - (2)- 0.813"¢ HOLES 7.625" PLAN

BRACKET DETAIL B

