



DEPARTMENT OF HEALTH CARE ACCESS AND INFORMATION
OFFICE OF STATEWIDE HOSPITAL PLANNING AND DEVELOPMENT

APPLICATION FOR HCAI PREAPPROVAL OF
MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY

APPLICATION #: OPM-0627

HCAI Preapproval of Manufacturer's Certification (OPM)

Type: New Renewal/Update

Manufacturer Information

Manufacturer: Hospital Systems, Inc.

Manufacturer's Technical Representative: Becca Teutle

Mailing Address: 750 Garcia Ave., Pittsburg, CA 94565

Telephone: (925) 427-7800

Email: bteutle@hsiheadwalls.com

Product Information

Product Name: HSI AXIOM & ELOQUENCE HEADWALLS

Product Type: Hospital Patient Headwalls

Product Model Number: AXIOM; ELOQUENCE

General Description: Patient Headwalls

Applicant Information

Applicant Company Name: Hospital Systems, Inc.

Contact Person: Becca Teutle

Mailing Address: 750 Garcia Ave., Pittsburg, CA 94565

Telephone: (925) 427-7800

Email: bteutle@hsiheadwalls.com

Title: President

"A healthier California where all receive equitable, affordable, and quality health care"

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY





**DEPARTMENT OF HEALTH CARE ACCESS AND INFORMATION
OFFICE OF STATEWIDE HOSPITAL PLANNING AND DEVELOPMENT**

Registered Design Professional Preparing Engineering Recommendations

Company Name: CYS STRUCTURAL ENGINEERS, INC.
 Name: Dieter Siebald California License Number: S4346
 Mailing Address: 2710 Gateway Oaks Drive, Suite 190N, Sacramento, CA 95833
 Telephone: (916) 920-2020 Email: dieters@cyseng.com

HCAI Special Seismic Certification Preapproval (OSP)

Special Seismic Certification is preapproved under OSP OSP Number: _____

Certification Method

Testing in accordance with: ICC-ES AC156 FM 1950-16
 Other(s) (Please Specify): _____

*Use of criteria other than those adopted by the California Building Standards Code, 2022 (CBSC 2022) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2022 may be used when approved by HCAI prior to testing.

- Analysis
- Experience Data
- Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

HCAI Approval

Date: 2/7/2025
 Name: Timothy Piland Title: Senior Structural Engineer
 Condition of Approval (if applicable): _____

"A healthier California where all receive equitable, affordable, and quality health care"

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY



AXIOM & ELOQUENCE HEADWALLS




TABLE OF CONTENTS OPM-0627

	PAGE
GENERAL NOTES	2
ABBREVIATIONS	3
DESIGN CRITERIA & LOAD COMBINATIONS	4
VERTICAL SURFACE MOUNTED HEADWALL ELEVATIONS & DETAILS	5
HORIZONTAL FLUSH MOUNTED HEADWALL ELEVATION & DETAILS	6
VERTICAL FLUSH MOUNTED HEADWALL ELEVATION & DETAILS	7
HORIZONTAL SURFACE MOUNTED HEADWALL ELEVATIONS & DETAILS	8
SURFACE MOUNTED HEADWALL ATTACHMENT TO STUD WALLS	9-10

- NOTES:**
1. THESE DRAWINGS ARE PREPARED FOR HOSPITAL SYSTEM, INC., PITTSBURG, CA
 2. THE CONTRACTOR AND INSPECTOR OF RECORD SHALL OBTAIN A COPY OF THIS PRE-APPROVAL FROM THE OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT (OSHPD) PRE-APPROVAL PROGRAMS WEBSITE.
 3. THIS PRE-APPROVAL COVERS THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE SUPPORTING STRUCTURE. THE EQUIPMENT IS SUPPLIED BY THE MANUFACTURER. THE SCREWS & BACKING PLATES SHOWN IN THIS OPM SHALL BE SUPPLIED AND INSTALLED BY THE CONTRACTOR.



SHEET TITLE: TABLE OF CONTENTS

 CYS STRUCTURAL ENGINEERS, INC. 2710 GATEWAY OAKS DRIVE, SUITE 190N SACRAMENTO, CA 95833	TEL (916) 920-2020 www.cyseng.com	Job No: 24070 Date: 01-08-2025 Page: 1 of 10
--	--------------------------------------	--


C:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025--02:37pm Login:FalkR DimScale:1 LTRScale:6

AXIOM & ELOQUENCE HEADWALLS



GENERAL NOTES:

1. THIS OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2022. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM MUST BE BASED ON THE CBC 2022.
2. IT IS THE RESPONSIBILITY OF THE STRUCTURAL ENGINEER OF RECORD FOR A SITE SPECIFIC PROJECT TO VERIFY:
 - A. THE ADEQUACY OF THE NEW OR (E) STRUCTURE TO RESIST THE FORCES & WT SPECIFIED FOR EA EQUIP IN ADDITION TO ALL OTHER LOADS. PROVIDE & DESIGN SUPPLEMENTARY MEMBERS AS REQ.
 - B. THAT THE INSTALLATION IS IN CONFORMANCE W/ THE CBC 2022 & W/ THE DETAILS SHOWN IN THIS PRE-APPROVAL.
 - C. THAT THE ACTUAL EQUIP'S WT, CENTER OF GRAVITY (CG) LOCATION, ATTACHMENT LOCATIONS, ATTACHMENT DETAILS, & THE MATERIAL & GA OF THE EQUIP WHERE ATTACHMENTS ARE MADE, AGREE W/ THE INFO SHOWN ON THE PRE-APPROVAL DOCUMENTS.
 - D. THAT THE PROJECT SPECIFIC VALUES OF S_{DS} & z/h RESULT IN SEISMIC FORCES THAT DO NOT EXCEED THE VALUES PROVIDED IN THE DESIGN CRITERIA.
3. ONE (1) CASE OF ATTACHMENT IS SPECIFIED & PRESENTED IN THIS PRE-APPROVAL:



CASE 1: ATTACHMENT DETAILS LOCATED AT FLRS AT OR ABV THE BASE OF A BLDG ($z/h \leq 1.0$), IT IS ASSUMED THAT THE WALLS ARE BUILT OF $\frac{5}{8}$ " THK GWB OVER 20 GA MIN STUD WALLS. MAY BE USED AT ANY GEOGRAPHICAL LOCATION IN THE STATE OF CALIFORNIA WHERE S_{DS} IS LESS THAN OR EQ TO 2.5.

4. SHEET METAL SCREWS SHALL BE HILTI SELF-DRILLING & SELF-PIERCING SCREWS PER ICC-ES ESR-2196.
5. SHEET METAL SCREWS SHOWN TO PROJECT THROUGH FRAMING MEMBERS SHALL PROJECT BEYOND THE MEMBER BY 3 FULL THREADS.
6. THE CONTRACTOR SHALL BE RESPONSIBLE FOR SELECTING THE PROPER SHEET METAL SCREW FEATURES UNLESS NOTED:
 - A. WAFER HEAD FASTENERS SHALL BE USED AT FRAMING CONNECTIONS COVERED WITH PLYWOOD, GYPSUM BOARD OR OTHER MATERIAL THAT MAY BE IMPEDED BY THE PROJECTION OF THE FASTENER HEAD.
 - B. HEX WASHER HEAD FASTENERS SHALL BE USED AT ALL OTHER CONDITIONS.
 - C. THREAD PITCH SHALL BE COMPATIBLE WITH THE THICKNESS OF THE PARTS BEING CONNECTED. THINNER GAUGE PARTS REQUIRE COARSER THREADS COMPARED TO THICKER GAUGE PARTS.



SHEET TITLE: GENERAL NOTES

CYS STRUCTURAL ENGINEERS, INC.

2710 GATEWAY OAKS DRIVE, SUITE 190N
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	24070
Date:	01-08-2025
Page:	2 of 10

AXIOM & ELOQUENCE HEADWALLS



ABBREVIATIONS:

Ω_o	SEISMIC OVERSTRENGTH FACTOR	INFO	INFORMATION
@	AT	JT	JOINT
ABV	ABOVE	KSI	KIPS PER SQUARE INCH
ALUM	ALUMINUM	LBS	POUNDS
ASCE	AMERICAN SOCIETY OF CIVIL ENGINEERS	LL	LIVE LOAD
ASD	ALLOWABLE STRESS DESIGN	LRFD	LOAD AND RESISTANCE FACTOR DESIGN
ASTM	AMERICAN SOCIETY FOR TESTING & MATERIALS	MAX	MAXIMUM
BLDG	BUILDING	MFR	MANUFACTURER
BLW	BELOW	MIN	MINIMUM
CBC	CALIFORNIA BUILDING CODE	MTL	METAL
CG	CENTER OF GRAVITY	NO. (#)	NUMBER OR POUNDS
\bar{C}	CENTERLINE	NWC	NORMAL WEIGHT CONCRETE
CLR	CLEAR	OPM	OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION
CONC	CONCRETE	OSHPD	OFFICE OF STATEWIDE HEALTH PLANNING & DEVELOPMENT
CONT	CONTINUOUS	PCF	POUNDS PER CUBIC FOOT
DF	DOUG FIR	PERP	PERPENDICULAR
DIA (ϕ)	DIAMETER	PG	PAGE
DL	DEAD LOAD	\bar{P}	PLATE
(E)	EXISTING	PSI	POUNDS PER SQUARE INCH
EA	EACH	PSF	POUNDS PER SQUARE FOOT
ELEV	ELEVATION	REQ	REQUIRED
EQ	EQUAL	SEOR	STRUCTURAL ENGINEER OF RECORD
EQUIP	EQUIPMENT	SLWC	SAND-LIGHTWEIGHT CONCRETE
ES	EACH SIDE	SMS	SHEET METAL SCREW
f'_c	MINIMUM ULTIMATE COMPRESSIVE STRENGTH OF CONCRETE	SPCG	SPACING
FLG	FLANGE	STL	STEEL
FLR	FLOOR	T	TENSION
FT (')	FOOT/FEET	THK	THICK/THICKNESS
F_p	HORIZONTAL SEISMIC FORCE PER ASCE 7-16 SEISMIC FORCE REQUIREMENTS	Tu	ANCHORAGE TENSION REACTION DUE TO SEISMIC FORCE
F_v	VERTICAL SEISMIC DESIGN FORCE PER ASCE 7-16 SECTION 12.4-4	TYP	TYPICAL
F_y	SEISMIC DESIGN FORCE REQUIREMENTS SPECIFIED MINIMUM YIELD STRESS OF STEEL	V	SHEAR
GA	GAUGE	Vu	ANCHORAGE SHEAR REACTION DUE TO SEISMIC FORCE
GR	GRADE	W/	WITH
GWB	GYPNUM WALLBOARD	Wp	OPERATING WEIGHT
HORIZ	HORIZONTAL	WS	WOOD SCREW
HT	HEIGHT	WT	WEIGHT
ICC	INTERNATIONAL CODE COUNCIL		
IN (")	INCH		



SHEET TITLE: ABBREVIATIONS

	CYS STRUCTURAL ENGINEERS, INC.		Job No: 24070
	2710 GATEWAY OAKS DRIVE, SUITE 190N SACRAMENTO, CA 95833		Date: 01-08-2025
	TEL (916) 920-2020 www.cyseng.com		Page: 3 of 10

c:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR Dimscale:1 LTScale:6

AXIOM & ELOQUENCE HEADWALLS

DESIGN CRITERIA

SUPPORT & ATTACHMENT DESIGN IS PER 2022 CBC AT ASD LEVEL FORCES. OTHER RIGID COMPONENTS LOW DEFORMABILITY ELEMENTS & ATTACHMENTS.. PER TABLE 13.5-1 OF ASCE 7-16 SUPPLEMENT #1

$$a_p = 1.0 \quad R_p = 1.5 \quad I_p = 1.5$$

MAX W_p AS SHOWN ON PG 5.

FOR CASE 1 – UPPER FLRS ABV THE BASE, $z/h \leq 1.0$

$$S_{DS} = 2.5$$

$$0.7 F_p = \frac{(0.7) 0.4 a_p S_{DS} W_p (1+2 z/h)}{(R_p/I_p)} = 2.1 W_p \quad \text{ASCE 7-16 (13.3-1)}$$

$$0.7 F_p \text{ (MAX)} = (0.7) 1.6 S_{DS} I_p W_p = 4.2 W_p \quad \text{ASCE 7-16 (13.3-2)}$$

$$0.7 F_p \text{ (MIN)} = (0.7) 0.3 S_{DS} I_p W_p = 0.7875 W_p \quad \text{ASCE 7-16 (13.3-3)}$$

$$0.7 (E_v+F_v) = (0.7) \pm 0.2 S_{DS} W_p = 0.35 W_p \quad \text{ASCE 7-16 (12.4-4)}$$

LOAD COMBINATIONS

$$(1.0+0.14 S_{DS}) D+0.7E \quad \text{ASD}$$

LOAD COMBINATIONS WERE RUN FOR 100% OF HORIZ FORCE IN ONE DIRECTION & 30% OF HORIZ FORCE IN THE PERP DIRECTION.



SHEET TITLE: DESIGN CRITERIA & LOAD COMBINATIONS



CYS STRUCTURAL ENGINEERS, INC.


2710 GATEWAY OAKS DRIVE, SUITE 190N
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No:	24070
Date:	01-08-2025
Page:	4 of 10

C:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR Dimscale:1 LTScale:6

AXIOM & ELOQUENCE HEADWALLS



FRONT ELEV
SURFACE MOUNTING
PANEL

ATTACH HEADWALL TO BACKING PLATE W/ #10 SMS AT MFR PROVIDED MOUNTING HOLES. TYPICAL OF 3 SCREWS MINIMUM TOP, BOTT, & MID-HT (9 MIN TOTAL)

SIDE ELEV

NOTE:
FOR HEADWALL ATTACHMENT TO BACKING PLATE & BACKING PLATE ATTACHMENT TO STUD WALLS, SEE PGS 9 & 10.



SHEET TITLE: VERTICAL SURFACE MOUNTED HEADWALL
ELEVATIONS & DETAILS



CYS STRUCTURAL ENGINEERS, INC.

2710 GATEWAY OAKS DRIVE, SUITE 190N
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com


Job No:	24070
Date:	01-08-2025
Page:	5 of 10

c:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR Dimscale:1 LTScale:6

AXIOM & ELOQUENCE HEADWALLS



NOTE:
COORDINATE OPENING SIZE W/
HEADWALL MFR PRIOR TO INSTALLATION



(E) MTL STUD SUPPORTING
WALL BY OTHERS
(16" MAX, 20 GA MIN,
Fy= 33 KSI MIN).
SEE GENERAL NOTE 2

#8 SMS HEADWALL TO HDR
AT MFR PROVIDED
MOUNTING HOLES (16" OC)

#8 SMS HEADWALL TO
STUD AT MFR PROVIDED
MOUNTING HOLES (16" OC)

HEADWALL FLG,
CONT AROUND
PERIMETER OF UNIT


12 GA ANODIZED ALUM
6063-T5 HEADWALL, TYP.
Fy= 21 KSI MIN (SHOWN
HATCHED)

#8 SMS HEADWALL TO
TRACK AT MFR PROVIDED
MOUNTING HOLES (16" OC)

1- LAYER 5/8" GWB
BOTH SIDES

FRONT ELEV
FLUSH MOUNTED PANEL

SECTION A-A



PLAN VIEW B-B



SHEET TITLE: AXIOM ELOQUENCE HORIZONTAL FLUSH MOUNTED HEADWALL
ELEVATION & DETAIL

CYS STRUCTURAL ENGINEERS, INC.


2710 GATEWAY OAKS DRIVE, SUITE 190N
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No: 24070
Date: 01-08-2025
Page: 6 of 10

c:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR DimScale:1 LTScale:6

AXIOM & ELOQUENCE HEADWALLS



NOTE:
COORDINATE OPENING SIZE
W/ HEADWALL MFR PRIOR
TO INSTALLATION

HEADWALL FLG, CONT
AROUND PERIMETER OF UNIT

#8 SMS HEADWALL TO HDR
AT MFR PROVIDED
MOUNTING HOLES (16" OC)

#8 SMS HEADWALL TO
STUD AT MFR PROVIDED
MOUNTING HOLES (16" OC)


12 GA ANODIZED ALUM
6063-T5 HEADWALL, TYP.
Fy= 21 KSI MIN
(HATCHED FOR CLARITY)

(E) STUD

#8 SMS HEADWALL TO
TRACK AT MFR PROVIDED
MOUNTING HOLES (16" OC)

16 GA MIN STUD BY
SEOR, ES OF HEADWALL

MTL STUD SUPPORTING WALL
BY OTHERS (16" MAX, 20 GA
MIN, Fy= 33 KSI MIN). SEE
GENERAL NOTE 2




FRONT ELEV
FLUSH MOUNTED
PANEL

16 GA MIN TRACK
BY OTHERS

1- LAYER 5/8" GWB
BOTH SIDES


SECTION A-A



PLAN VIEW B-B




SHEET TITLE: AXIOM ELOQUENCE VERTICAL FLUSH MOUNTED HEADWALL
ELEVATIONS & DETAIL

 CYS STRUCTURAL ENGINEERS, INC. 2710 GATEWAY OAKS DRIVE, SUITE 190N SACRAMENTO, CA 95833	TEL (916) 920-2020 www.cyseng.com	Job No: 24070 Date: 01-08-2025 Page: 7 of 10
--	--------------------------------------	--

C:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR Dimscale:1 LTScale:6

AXIOM & ELOQUENCE HEADWALLS



ATTACH HEADWALL TO BACKING PLATE W/ #10 SMS AT MFR PROVIDED MOUNTING HOLES. TYPICAL OF 3 SCREWS MINIMUM T&B (6 MIN TOTAL)

12 GA ANODIZED ALUM 6063-T5 HEADWALL, TYP. Fy = 21 KSI, MIN

SUPPORTING FLR

FRONT ELEV

SURFACE MOUNTING PANEL


MTL STUD SUPPORTING WALL BY OTHERS (16" MAX, 20 GA MIN, Fy = 33 KSI, MIN). SEE GENERAL NOTE 2

12 GA ANODIZED ALUM 6063-T5 HEADWALL, TYP. Fy = 21 KSI, MIN

NOTE:
FOR HEADWALL ATTACHMENT TO BACKING PLATE & BACKING PLATE ATTACHMENT TO STUD WALLS, SEE PGS 9 & 10.

TYP WALL ATTACHMENT. SEE NOTE

SIDE ELEV



SHEET TITLE: HORIZONTAL SURFACE MOUNTED HEADWALL ELEVATIONS & DETAIL



CYS STRUCTURAL ENGINEERS, INC.

2710 GATEWAY OAKS DRIVE, SUITE 190N
SACRAMENTO, CA 95833

TEL (916) 920-2020
www.cyseng.com

Job No: 24070
Date: 01-08-2025
Page: 8 of 10

C:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR DimScale:1 LTScale:6



AXIOM & ELOQUENCE HEADWALLS



MAX ASD FORCES
AT EA SCREW (LBS)
PLATE TO STUDS

16 GA MTL	
T	V
CASE 1	47#
	36#

MTL STUD SUPPORTING WALL
BY OTHERS (16" MAX, 20 GA
MIN, Fy= 33 KSI MIN). SEE
GENERAL NOTE 2



CONNECTION PATTERN

SCREW INSTALL VARIANCE



SHEET TITLE: SURFACE MOUNTED HEADWALL
ATTACHMENT TO STUD WALLS


	CYS STRUCTURAL ENGINEERS, INC.		Job No: 24070
	2710 GATEWAY OAKS DRIVE, SUITE 190N	TEL (916) 920-2020	Date: 01-08-2025
	SACRAMENTO, CA 95833	www.cyseng.com	Page: 9 of 10

c:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR Dimscale:1 LTScale:6

AXIOM & ELOQUENCE HEADWALLS



NOTE:
HEADWALL UNIT NOT SHOWN IN SECTION B-B FOR CLARITY.



SHEET TITLE: SURFACE MOUNTED HEADWALL ATTACHMENT TO STUD WALLS

<p>CYS STRUCTURAL ENGINEERS, INC. 2710 GATEWAY OAKS DRIVE, SUITE 190N SACRAMENTO, CA 95833</p>	TEL (916) 920-2020 www.cyseng.com	Job No: 24070 Date: 01-08-2025 Page: 10 of 10
---	--------------------------------------	---

c:\Users\FalkR\AppData\Local\Temp\AcPublish_10664\S1_Axiom Headwall.dwg Time:Jan27,2025-02:37pm Login:FalkR Dimscale:1 LTScale:6