

DEPARTMENT OF HEALTH CARE ACCESS AND INFORMATION FACILITIES DEVELOPMENT DIVISION

whith his and				
APPLICATION FOR HCAI PREAPPROVAL OF	OFFICE USE ONLY			
MANUFACTURER'S CERTIFICATION (OPM)	APPLICATION #: OPM-0654			
HCAI Preapproval of Manufacturer's Certification (OPM)				
Type: X New Renewal/Update				
Manufacturer Information				
Manufacturer: Skytron LLC				
Manufacturer's Technical Representative: BOB VREELAND				
Mailing Address: 5085 Corporate Exchange Blvd SE, Grand Rapids, MI 49512				
Telephone: (616) 656-1189 Email: bvreeland@skytron.c	om			
Product Information	Z			
Product Name: Skytron Freedom Series booms with optional Skytron lights	G			
Product Type: Hospital Ceiling Mounted Booms				
Product Model Number: F110, F120, F130, F10L, F20H, F200, F212, F221, F22 F440	2, F <mark>310,</mark> F320, F330, F340, F350, F410, F420,			
General Description: Skytron's Freedom Series ceiling mounted booms are mou include optional surgical, procedure or examination lights. for equipment, monitors and lights.				
Applicant Information				
Applicant Company Name: ISAT SEISMIC BRACING				

Telephone: (757) 817-1893

Contact Person: KRIS CLEMENTE

Email: kclemente@isatsb.com

Title: ENGINEER

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

Mailing Address: 1020 CREWS RD STE Q, MATTHEWS, ND 28105





STATE OF CALIFORNIA - HEALTH AND HUMAN SERVICES AGENCY



DEPARTMENT OF HEALTH CARE ACCESS AND INFORMATION FACILITIES DEVELOPMENT DIVISION

Registered Design Professonal Preparing Engineering Recommendations				
Company Name: ISAT SEISMIC BRACING				
Name: Anthony Rubalcava California License Number: S4710				
Mailing Address: 14848 Northam Street, La Mirada, CA 90638				
Telephone: (714) 356-3286 Email: arubalcava@isatsb.com				
HCAI Special Seismic Certification Preapproval (OSP)				
Special Seismic Certification is preapproved under OSP OSP Number:				
an CODE a				
Certification Method				
Testing in accordance with: ICC-ES AC156 FM 1950-16				
Other(s) (Please Specify): Equipment is considered to be rugged. OPM is for anchorage to concrete slabs				
*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by HCAI prior to testing.				
X Analysis				
Experience Data DATE: 02/18/2022				
Combination of Testing, Analysis, and/or Experience Data (Please Specify):				
COE:				
HCAI Approval				
Date: <u>2/18/2022</u>				
Name: William Staehlin Title: Senior Structural Engineer				
Condition of Approval (if applicable):				

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"







Innovation • Engineering • BIM • Fabrication

A Division of Tomarco Contractor Specialties

Submittal Documents

OPM-0654

HCAI OPM-0654

SUPPORT AND ATTACHMENT OPM CONSTRUCTION DRAWINGS FOR SKYTRON FREEDOM SERIES BOOMS WITHOPTIONAL SKYTRON LIGHTS



ISAT 1020 Crews Road Suite Q Matthews, N.C. 28105 704-841-4080

FILE NO.: NC32527.200690 "Empowered by Experience"

REV 1

HCAI OPM-0654 CONST DWG - i

02/01/2022

S 5867



HCAI OPM-0654 OPM DRAWING INDEX

INDEX CODE CO				
Cover Page	рi			
Index Page	p ii			
Cover Page Index Page General Notes OPM-0654 Attachment Notes	p 1			
Attachment Notes	p 2			
Skytron Freedom Series at Metal Deck Slabs	р3			
Skytron Freedom Series at Solid Concrete Slabs				
Skytron Freedom Series Miscellaneous Steel Page 1				
Skytron Freedom Series Miscellaneous Steel Page 2				
Skytron Freedom Series Miscellaneous Steel Page 3				
Skytron Freedom Series Top Plate Details				
Skytron Freedom Series Mounting Plate Details				
Skytron Freedom Series Forces and Moments				
Skytron Freedom Series Anchor Forces	p 11			

HCAI OPM-0654

MANUFACTURE: SKYTRON
EQUIPMENT TYPE: CEILING MOUNTED BOOM
AND SURGICAL LIGHTS (NON-MOTORIZED)

GENERAL NOTES:

- 1. THIS HCAI PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2019. THE DESIGN FORCES FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2019 AND ASCE 7-16.
- 2. SEISMIC CRITERIA USED: $S_{DS} = 2.5 I_P = 1.5 ap = 1.0 Rp = 1.5$ (OTHER MECHANICAL AND ELECTRICAL COMPONENTS) $z/h \le 1.0 FpH = 3.00 AND FpV = 0.50$.
- 3. SUPPORT AND ATTACHMENT FORCES ARE DETERMINED USING ASCE 7-16 CHAPTER 13 "SEISMIC DESIGN REQUIREMENTS FOR NONSTRUCTURAL COMPONENTS". AN OVERSTRENGTH FACTOR Ω_0 = 2.0 IS USED FOR CONCRETE MATERIALS PER ASCE 7-16 TABLE 13.6-1. LOADS SHOWN ARE STRENGTH DESIGN LOADS PER CBC 2019.
- 4. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- 5. STEEL MATERIALS: ANGLES AND PLATE ASTM A36, ROUND HSS ASTM A500 GRADE B (FY = 42 KSI), RECTANGULAR HSS ASTM A500 GRADE B (FY = 46 KSI), WELDING ELECTRODES 70 KSI TENSILE STRENGTH, ALL THREAD ROD ASTM A193 GR. B7, NUTS ASTM A194 GR.2H, WASHERS ASTM F436. ALL THREAD ROD SUPPLIED BY SKYTRON ARE ASTM A193 GR.B7.
- 6. CONCRETE SLABS: //
 - a. FOR ELEVATED SOLID CONCRETE SLABS: 6" THICKNESS OF NORMAL WEIGHT CONCRETE WITH 4000 PSI MINIMUM STRENGTH.
 - b. METAL DECK: 3" DEEP COMPOSITE STEEL DECK, 20 GAGE MINIMUM, 4 1/2 INCH MINIMUM BOTTOM FLUTE WIDTH AND FLUTE SPACING IS 12", WITH 3 1/4 INCH SAND LIGHT WEIGHT CONCRETE CONCRETE COVER AT 4000 PSLMINIMUM STRENGTH.
- 7. POST-INSTALLED CONCRETE ANCHORS: HILT! KWIK BOLT TZ2 (ESR-4266) 5/8" DIAMETER x 4 3/4" MIN. HOLE DEPTH (4" EFFECTIVE EMBEDMENT) AND 40 FT-LBS INSTALLATION TORQUE; 12" MINIMUM EDGE DISTANCE. ANCHOR SPACING IS SHOWN ON THE OPM DRAWINGS.
- 8. FREEDOM SURGICAL LIGHTS THAT ARE PART OF THIS OPM ARE ATTACHED TO THE FREEDOM BOOMS AND HUB.

OPM-0654 SKYTRON FREEDOM SERIES GENERAL NOTES



HCAI OPM-0654

MANUFACTURE: SKYTRON

EQUIPMENT TYPE: CEILING MOUNTED BOOM AND SURGICAL LIGHTS (NON-MOTORIZED)

ATTACHMENT NOTES:


- 1. THIS HCAI PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2019. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2019.
- 2. BRACE ARM INCLINATION MAY VARY FROM 30° TO 45° FROM HORIZONTAL.
- 3. PERIODIC SPECIAL INSPECTION PER CBC 2019 SECTION 1705 AND TABLE 1705.3 INCLUDING VERIFICATION OF ANCHOR TYPE, ANCHOR DIMENSIONS, CONCRETE TYPE, CONCRETE COMPRESSIVE STRENGTH, ANCHOR SPACING, EDGE DISTANCES, CONCRETE MEMBER THICKNESS, TIGHTENING TORQUE, HOLE DIMENSIONS, ANCHOR EMBEDMENT AND ADHERENCE TO THE MANUFACTURER'S PRINTED INSTALLATION INSTRUCTIONS. IN ADDITION, FOLLOW THE PROVISIONS OF THE 2019 CALIFONIA BUILDING CODE SECTION 1910A.5 BY CONFIRMING THE INSTALLATION TORQUE SPECIFIED BY THE MANUFACTURER. TESTING IS NOT TO OCCUR UNTIL A MINIMUM OF 24 HOURS HAS ELAPSED AFTER THE INSTALLATION OF THE SUBJECT ANCHORS. TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR. TEST 50% OF THE ANCHORS FOR EACH PIECE OF EQUIPMENT. USING A CALIBRATED TORQUE WRENCH VERIFY THE INSTALLATION TORQUE IS OBTAINED WITHIN 1/2 TURN OF THE NUT. REPORT OF TEST RESULTS ARE TO BE SUBMITTED TO THE ENFORCEMENT AGENGY. THE SEOR SHALL POVIDE REMEDIAL ANCHORAGE DETAILS IN THE EVENT THAT AN ANCHOR FAILS TO MEET THE TEST REQUIREMENTS.
- 4. WELDS ARE TO BE VISUALLY INSPECTED BY A QUALIFIED WELDING INSPECTOR.
- 5. STRENGTH DESIGN WAS USED FOR ANCHOR FORCE CALCULATIONS INCLUDING Ω_0 PER ACI 318-14 WHERE REQUIRED FOR ATTACHMENT TO CONCRETE.
- 6. EXCERCISE DUE CARE WHEN DRILLING POST-INSTALLED ANCHORS TO AVOID DAMAGING CONCRETE REINFORCEMENT OR TENDONS.
- 7. PROVIDE FULL THREAD ENGAGEMENT OF NUT AND WASHER.


RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD


- 1. CONFIRM THE MATERIAL PROPERTIES AND THICKNESS OF THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ATTACHED MEETS THE REQUIREMENTS OF THIS OPM.
- 2. PROVIDE A PLAN FOR INSPECTION OF SUPPORTS AND ATTACHMENTS AND VERIFY ITS IMPLEMENTATION.
- 3. CONFIRM THE SPECIFIED MINIMUM CONCRETE EDGE DISTANCES ARE MAINTAINED BASED ON THE ACTUAL EQUIPMENT LOCATION. VERIFY THAT EXISTING OR NEW ANCHORS ARE AN ADEQUATE DISTANCE FROM THIS UNIT'S ATTACHMENT.
- 4. VERIFY THAT THE EXISTING STRUCTURE IS ADEQUATE FOR THE IMPOSED DEAD, LATERAL AND TENSION FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- 5. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH CBC 2019 AND WITH DETAILS IN THIS OPM INCLUDING MATERIALS AND DIMENSIONS OF THE SUPPORT WHERE THE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN.
- 6. VERIFY THAT THE PROJECT SPECIFIC S_{DS} AND z/h VALUES RESULT IN SEISMIC FORCES (Eh AND Ev) THAT DO NOT EXCEED THE VALUES SHOWN IN THESE DETAILS.


OPM-0654 SKYTRON FREEDOM SERIES ATTACHMENT NOTES























		Seismic Design Forces at Ultimate Loads							
Model Number⁵	Weight at 1g (lbs.)	Vertical Force (lbs.)	Lateral force (lbs.)	Mx (ftlbs.)	My (ftlbs.)	Mz (ftlbs.)			
F130 SERIES	761	1,293	2,282	693	589	1,039			
F200 SERIES	1050	1,785	3,150	761	766	1,352			
F222 SERIES	1050	1,785	3,150	780	766	1,352			
F350 SERIES	1454	2,473	4,363	1,159	1,221	2,155			
F440 SERIES	1442	2,452 - 065	4,326	1,068	1,142	2,015			
DESIGN CASE	1500	2,500	4,500	1,250	1,250	2,200			

- 1. WEIGHTS AND MOMENTS ARE FACTORED LOADS USING STRENGTH DESIGN AND INCLUDE THE FOLLOWING FACTORS: DL = 1.2, FpV = 0.5 AND FpH = 3.0. MAXIMUM LOADING AFFECT IS USED FOR DESIGN.
- 2. Mz are moments about the vertical axis due to vertical seismic acceleration times the horizontal eccentricity.
- 3. Mx are moments about the horizontal axis due to horizontal seismic acceleration times the vertical eccentricity.
- My ARE MOMENTS ABOUT THE HORIZONTAL AXIS DUE TO VERTICAL SEISMIC ACCELERATION TIMES THE HORIZONTAL ECCENTRICITY.
- 5. THE MODELS LISTED ARE THE MAXIMUM OF EACH SERIES, SEE PAGE "ANCHOR FORCES" FOR ALL MODELS IN THIS OPM

OPM-0654 SKYTRON FREEDOM SERIES FORCES AND MOMENTS



EQUIPMENT ATTACHMENT DATA AND ANCHORAGE FORCES ³							
	Bolt Diameter	Install	All Concrete	Top Plate Max. 2		Brace Anchoragge 2	
	at Mounting	Torque	Anchorage ¹	Tension	Shear	Tension	Shear
Model Series	Plate ⁴	(ft-lbs)		(lbs)	(lbs)	(lbs)	(lbs)
F110 Low duty single	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F120 Low duty double	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F130 Low duty triple	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F10L Stackrotation unit Low Duty	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F20H Stackrotation unit Heavy Duty	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F200 Single Q1	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F212 Twin TL+Q1	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F221 Twin Q1+TL	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F222 Twin Q1+Q1	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F310 Heavy duty single	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F320 Heavy duty + Single Low duty	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F330 Heavy duty + Double Low duty	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F340 Heavy duty + Heavy tandem	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F350 Heavy duty+Heavy tandem+Low duty	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F410 Heavy duty double	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F420 Heavy duty double + Low duty	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771
F440 Heavy duty double + Heavy tandem	1-1/4"	660	Hilti KB TZ2 5/8"	2728	1442	1661	1771

- 1. SEE PAGE "GEN NOTE" GENERAL NOTE 7 HILTI KWIK BOLT TZ2 WITH 4" EFFECTIVE EMBEDMENT AND 40 FT-LBS INSTALLATION BY: William Staehlin TORQUE PER ESR-4266.
- 2. INCLUDES AN OVERSTRENGTH FACTOR OF $\Omega_0 = 2.0$.
- 3. FORCES SHOWN ARE FOR EACH ANCHOR. TE: 02/18/2022
- 4. SEE PAGE "MOUNTING PL" FOR BOLT INFORMATION.

OPM-0654 SKYTRON FREEDOM SERIES WITH STELLAR XL SURGICAL LIGHTS **EQUIPMENT MOUNTING AND ANCHOR FORCES**



ISAT TOTAL SUPPORT

A Division of Tomarco Contractor Specialties 1020 Crews Road, Suite Q, Matthews, NC 28105 704-841-4080 www.isatsb.com

DRAWN BY: KJC

DATE: 12/06/2021

REVISED BY: KJC

DATE: 01/24/2022 REV NO: 1

SCALE

PAGE N.T.S. ANCH FORCES