### APPLICATION FOR OSHPD SPECIAL SEISMIC **CERTIFICATION PREAPPROVAL (OSP) APPLICATION #:** OSP - 0397 - 10 **OSHPD Special Seismic Certification Preapproval (OSP) Manufacturer Information** GE HEALTHCARE Manufacturer: Manufacturer's Technical Representative: Tom Farnow Mailing Address: 3000 N. Grandview Blvd., Waukesha, WI 53188-1696 Telephone: 888-406-1101 Email: Tom.Farnow@gehcseismic.com **Product Information** Product Name: LIGHTSPEED VCT System Product Type: MEDICAL IMAGING EQUIPMENT Product Model Number: See Attachment 1 (List all unique product identification numbers and/or part numbers) General Description: System components of multiple-component Computed Tomography (CT) medical imaging system. Equipment shall include all enhancements present in test specimens. Mounting Description: Rigid base mounted **Applicant Information** Applicant Company Name: EASE LLC Contact Person: JONATHAN ROBERSON, S.E. Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA 91709 Telephone: (406) 541-EASE (3273) Email: j.roberson@easeco.com I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2016. Signature of Applicant: Date: Title: Principal Engineer Company Name: EASE LLC

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs'





Page 1 of 3

OFFICE USE ONLY



Page 2 of 9

02/27/2017



# OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Seismic Parameters
Design in accordance with ASCE 7-10 Chapter 13: ⊠ Yes □ No
Design Basis of Equipment or Components (F <sub>p</sub> /W <sub>p</sub> ) = See Attachment 1, Table 2
S <sub>DS</sub> (Design spectral response acceleration at short period, g) = See Attachment 1, Table 2
a <sub>p</sub> (In-structure equipment or component amplification factor) = See Attachment 1, Table 2
R <sub>p</sub> (Equipment or component response modification factor) = See Attachment 1, Table 2
$Ω_0$ (System overstrength factor) = See Attachment 1, Table 2
I <sub>P</sub> (Importance factor) = <b>1.5</b>
z/h (Height factor ratio) = See Attachment 1, Table 2
Equipment or Component Natural Frequencies (Hz) = See Attachment 2
Overall dimensions and weight (or range thereof) = See Attachment 1, Table 1
Equipment or Components @ grade designed in accordance with ASCE 7-10 Chapter 15:   Yes   No
Design Basis of Equipment or Components (V/W) =
S <sub>DS</sub> (Design spectral response acceleration at short period, g) =
S <sub>D1</sub> (Design spectral response acceleration at 1 second period, g) =
R (Response modification coefficient ) =
$\Omega_0$ (System overstrength factor) =
C <sub>d</sub> (Deflection amplification factor) =
I <sub>p</sub> (Importance factor) = 1.5
Height to Center of Gravity above base =
Equipment or Component Natural Frequencies (Hz) =
Overall dimensions and weight (or range thereof) =
Tank(s) designed in accordance with ASME BPVC, 2015: ☐ Yes ☒ No
List of Attachments Supporting Special Seismic Certification
☐ Test Report(s) ☐ Drawings ☐ Calculations ☐ Manufacturer's Catalog
Other(s) (Please Specify): Attachments 1 & 2
OSHPD Approval (For Office Use Only) – Approval Expires on December 31, 2022
1// 1 / 1 / 1 / 1 / 1 / 1 / 2 / 2047
Signature: Date: February 27, 2017
Print Name: Timothy J. Piland Title: SSE
Special Seismic Certification Valid Up to : Sps (g) = See Above z/h = See Above
Condition of Approval (if applicable):

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





Page 3 of 3



#### **ATTACHMENT 1: SEISMIC CERTIFIED COMPONENTS**

ATTACHMENT PAGE | 1 OF 2

#### TABLE 1: SEISMIC CERTIFIED SYSTEMS & COMPONENTS

Manufacturer	GE HEALTHCARE	<b>■</b> /						
System	LIGHTSPEED VC	Γ /						
		MODEL	APPI	ROX. DIMENSION	IS (IN.)	MAX. WT.		
COMPONENT		NUMBER	W	D	Н	(LB.)	MOUNT	BASIS [1]
GANTRIES								
LightSpeed	VCT Gantry	5124069-11	89.25	39.65	74.6	4085	Floor	UUT1
LightSpeed	VCT Gantry	5124069-5	89.25	39.65	74.6	4085	Floor	SAME (UUT1)
LightSpeed	VCT Gantry	5124069-6	89.25	39.65	74.6	4085	Floor	SAME (UUT1)
LightSpeed	VCT Gantry	5124069-2	89.25	39.65	74.6	4130	Floor	UUT2
LightSpeed	VCT Gantry	5124069	89.25	39.65	74.6	4130	Floor	SAME (UUT2)
PATIENT TA	ABLES	•		1		<u> </u>	<u> </u>	
GT2000		5121647-3	25.6	114.5	19.2 / 41.3	1146 [2]	Floor	UUT4
GT2000		5121647	25.6	114.5	19.2 / 41.3	1146	Floor	SAME (UUT4)
GT2000		5121647-4	25.6	114.5	19.2 / 41.3	1146	Floor	SAME (UUT4)
GT1700		5122080-3	25.6	93.3	19.2 / 41.3	1047	Floor	INT
GT1700		5122080-4	25.6	93.3	19.2 / 41.3	1047	Floor	INT
GT1700		5122080	25.6	93.3	19.2 / 41.3	1047	Floor	INT
GT1700V		5122080-11	25.6	93.3	19.2 / 41.3	1059 <sup>[3]</sup>	Floor	UUT6
CONSOLES	[4	]				•	•	•
GOC 6.5		5212920-310	48.74	46.54 / 54.74	26.7 / 34.7	512	Floor	UUT5
GOC 6.6		5212920-350	48.74	46.54 / 54.74	26.7 / 34.7	405	Floor	INT
GOC 6.6		5212920-355	48.74	46.54 / 54.74	26.7 / 34.7	405	Floor	INT
GOC 6.6		5212920-360	48.74	46.54 / 54.74	26.7 / 34.7	405	Floor	INT
GOC 6.6		5212920-365	48.74	46.54 / 54.74	26.7 / 34.7	405	Floor	INT
GOC 6.6		5212920-150	48.74	46.54 / 54.74	26.7 / 34.7	405	Floor	UUT7
POWER DIS	TRIBUTION (PDU)	.1						
Power Distr		2326492-4	27.6	21.7	41.8	806	Floor	UUT3
Power Distr		2326492-60	27.6	21.7	41.8	818	Floor	SAME (UUT8)
Power Distr		2326492-61	27.6	21.7	41.8	818	Floor	UUT8
Mount	Floor (Rigid Base): fr support above the ba	ee-standing, base-mo				1		
Notes	SAME: Mod color, softwa     INT (Interpo established     Patient Table w     Patient Table w     Weight includes	ates that a test specin lel is physically, mechalare and/or GE manufalated or Extrapolated) through evaluation of veight does not include reight does not includes mounting brackets us listed above are manually in the series of the series and the series and the series are manufacted in the series are manufacted and series are	anically & electring location: indicates a material testing of other than 150 lb. since the 350 lb. since the 350 lb. sinced in testing.	trically the same a n. nodel that was not er, similar models i mulated patient lo mulated patient lo	s test specimen specifically testen the product lin ad present durin ad present durin	i. Difference is led, and by which e. ng testing. ng testing.	limited to mo	odel number,

#### **ATTACHMENT 1: SEISMIC CERTIFIED COMPONENTS**

ATTACHMENT PAGE | 2 OF 2

#### TABLE 2: ASCE 7-10 DESIGN BASIS FOR EQUIPMENT

COMPONENT	MODEL No.	F <sub>P</sub> /W <sub>P</sub>	S <sub>DS</sub>	z/h	a <sub>P</sub>	$R_P$	$\Omega_0$	
LightSpeed VCT Gantry	5124069-11 5124069-5 5124069-6	2.40	2.0	1.0	1	1 ½	1 ½	
Lightopeed von Gantry	5124069-2 5124069	1.13	2.5	0	,	1 /2	1 /2	
GT2000 Patient Table GT1700 Patient Table	5121647-3 5121647 5121647-4	2.40	2.0	1.0	1	1 ½	1 ½	
	5122080-3 5122080-4 5122080	1.13	2.5	0				
GT1700V Patient Table	5122080-11	2.4	2.0	1.0	1	1 ½	1 ½	
G11700V Falletil Table	3122000-11	1.17	2.6	0	Į.	1 /2	1 /2	
GOC6	5212920-310	1.50	2.0	1.0			_	
GOC6.6 GOC6.6	5212920-350 5212920-355	1.13	2.5	0	2 ½	6	2	
GOC6.6	5212920-150	1.50	2.0	1.0	2 ½	6	2	
GOC0.0	3212920-130	1.17	2.6	0	2 /2	0		
Power Distribution Unit	2326492-4	1.44	2.0	1.0	1	2 ½	2	
Fower Distribution Offic	2320492-4	1.13	2.5	0	ļ	Z /2		
Power Distribution Unit	2326492-60	1.44	2.0	1.0	1	2 ½	2	
Power Distribution Offic	2326492-61	1.17	2.6	0	]	Z /2	2	



### ATTACHMENT 2: TEST SPECIMEN SUMMARY

ATTACHMENT PAGE | 1 OF 4

UUT- 1	LightSpeed VCT Gantry
MANUFACTURER:	General Electric Company
IDENTIFICATION:	Model No.: 5124069-11
DESCRIPTION:	System Component of the LightSpeed VCT System. Gantry includes "Seismic Kit"
MOUNTING:	Floor mounted using (4) – 5/8" dia GR 8 hex head bolts(torqued to 60 lb-ft) w/ GEHC supplied Gantry foot assembly



PROPERTIES:										
	DIMENSIONS (in.)		LOWEST RESONANT FREQUENCY (Hz.)							
Width	Depth	Height		Weight (lb.)	X-Axis	Y	-Axis	Z-Axis		
89.25	39.65	74.6		4085	8.6		6.2	20.0		
SHAKE TABLE TEST PARAMETERS										
CODE	TEST CRITERIA	$S_{DS}$	z/h	I <sub>P</sub>	$A_{FLX-H}$	$A_{RIG-H}$	A <sub>FLX-\</sub>	, A <sub>RIG-V</sub>		
CBC 2013	ICC-ES AC156-12	2.0 2.5	1.0 0.0	1 1 5	3.2 2.5	2.4 1.0	1.34 1.68			
Unit maintained structural integrity and functionality after the ICC-ES AC 156 test										

UUT- 2	LightSpeed VCT Gantry
MANUFACTURER:	General Electric Company
IDENTIFICATION:	Model No.: 5124069-2
DESCRIPTION:	System Component of the LightSpeed VCT System. Gantry includes "Seismic Kit"
MOUNTING:	Floor mounted using (4) – 5/8" dia GR 8 hex head bolts(torqued to 60 lb-ft) w/ GEHC supplied Gantry foot assembly



PROPERTIES:										
	DIMENSIONS (in.)			LOW	ANT FREQUE	EQUENCY (Hz.)				
Width	Depth	Height		Weight (lb.)	Front-Axis	s Si	de-Axis	Vertical-Axis		
89.25	39.65	74.6		4130	8.8		7.1	21.7		
SHAKE TABLE TEST PARAMETERS										
CODE	TEST CRITERIA	$S_{DS}$	z/h	I <sub>P</sub>	$A_{FLX-H}$	$A_{RIG-H}$	A <sub>FLX-V</sub>	A <sub>RIG-V</sub>		
CBC 2013	ICC-ES AC156-12	2.0 2.5	1.0 0.0	1.5	3.2 2.5	2.4 1.0	1.34 1.68	0.54 0.68		
Unit maintained structural integrity and functionality after the ICC-ES AC 156 test										



### ATTACHMENT 2: TEST SPECIMEN SUMMARY

ATTACHMENT PAGE | 2 OF 4

UUT- 3	Power Distribution Unit
MANUFACTURER:	GE Hangwei Medical Systems CO. LTD.
IDENTIFICATION:	Model No.: 2326492-4
DESCRIPTION:	System Component of the LightSpeed VCT System.
MOUNTING:	Floor mounted using (4) – 3/8" dia. socket head bolts through GEHC supplied mounting brackets to aluminum interface plate.



PROPERTIES:											
DIMENSIONS (in.)						LOWEST RESONANT FREQUENCY (Hz.)					
Width	Depth	Height		Weight (lb.)		Front-Axis		Side-Axis		Vertical-Axis	
27.6	21.7	41.8		806-	+22 brackets	26.7		18.9		>50	
SHAKE TABLE TI	EST PARAMETERS										
CODE	TEST CRITERIA	S <sub>DS</sub>	z/ł	h	I <sub>P</sub>	A <sub>FLX-H</sub>	I	A <sub>RIG-H</sub>	RIG-H A <sub>FLX-\</sub>		$A_{RIG-V}$
CBC 2013	ICC-ES AC156-12	2.0 2.5	1.0 0.0	-	1.5	3.2 2.5		2.4 1.34 1.0 1.68			0.54 0.68
Unit maintained structural integrity and functionality after the ICC-ES AC 156 test											

UUT- 4	GT2000 Patient Table
MANUFACTURER:	GE Hangwei Medical Systems CO. LTD.
IDENTIFICATION:	Model No.: 5121647-3
DESCRIPTION:	System Component of the LightSpeed VCT System GT2000 N9 Patient Table Test specimen included a simulated patient load of 550 lb. Seismic Kit
MOUNTING:	Floor: (4) – 5/8" dia GR 8 hex head bolts (torqued to 60 lb-ft) w/ GEHC supplied patient table foot assembly.



PROPERTIES:											
DIMENSIONS (in.)						LOWEST RESONANT FREQUENCY (Hz					(Hz.)
Width	Depth	Height		W	eight (lb.)	Transverse-Axis		Longitudinal-Axis		Vertical-Axis	
25.6	114.5	19.2 / 41.3		1146	+550 Patient	2.7		7.1		5.7	
SHAKE TABLE TEST PARAMETERS											
CODE	TEST CRITERIA	S <sub>DS</sub>	z/ł	า	I <sub>P</sub>	A <sub>FLX-H</sub>		A <sub>RIG-H</sub> A <sub>FLX-1</sub>			$A_{RIG-V}$
CBC 2013	ICC-ES AC156-12	2.0 2.5	1.0 0.0		1.5	3.2 2.5		2.4 1.0	1.34 1.68		0.54 0.68
Unit maintained structural integrity and functionality after the ICC-ES AC 156 test											



### ATTACHMENT 2: TEST SPECIMEN SUMMARY

ATTACHMENT PAGE | 3 OF 4

UUT- 5	GOC6 VCT Console
MANUFACTURER:	GE Healthcare Japan Corporation
IDENTIFICATION:	Model No.: 5212920-310
DESCRIPTION:	System Component of the LightSpeed VCT System. GOC6.5 Console with Xtream HD upgrade
MOUNTING:	Floor mounted using (4) – 3/8" dia. socket head bolts through GEHC supplied mounting brackets to aluminum interface plate.



PROPERTIES:											
DIMENSIONS (in.)						LOWEST RESONANT FREQUENCY (Hz.)					
Width Depth Height					Veight (lb.)	Front-Axis		Side-Axis		Vertical-Axis	
48.74	46.54 to 54.74	26.7 to 34	4.7	501	+11 brackets	14.2		6.5		14.6	
SHAKE TABLE TEST PARAMETERS											
CODE	TEST CRITERIA	z/ł	h	I <sub>P</sub>	A <sub>FLX-H</sub>		A <sub>RIG-H</sub>	A <sub>FLX-V</sub>	,	$A_{RIG-V}$	
CBC 2013	ICC-ES AC156-12	2.0 2.5	1.0 0.0	-	1.5	3.2 2.5		2.4 1.0	1.34 1.68		0.54 0.68
Unit maintained structural integrity and functionality after the ICC-ES AC 156 test											

UUT- 6	GT1700V Table
MANUFACTURER:	GE Hangwei Medical Systems Co., LTD
IDENTIFICATION:	Model No.: 5122080-11
DESCRIPTION:	System component of the <b>Optima CT660 System</b> Test specimen included a simulated patient load of 350 lb.
MOUNTING:	Floor mounted using (4) – 5/8" dia GR 8 hex head bolts (torqued to 60 lb-ft) w/ GEHC supplied patient table foot assembly. Leveling foot bushing excluded.



PROPERTIES:											
DIMENSIONS (in.)						LOWEST RESONANT FREQUENCY (Hz.)					Hz.)
Width	Depth	Weight (lb.)		Transverse-Axis		Longitudinal-Axis		Vertical-Axis			
25.6	93.3	19.2 / 41	.3	1059+	350 Patient	3.9		15.2			14.2
SHAKE TABLE T	SHAKE TABLE TEST PARAMETERS										
CODE	TEST CRITERIA	h	I <sub>P</sub>	$A_{FLX-H}$		A <sub>RIG-H</sub>	A <sub>FLX-V</sub>		$A_{\text{RIG-V}}$		
CBC 2013 ICC-ES AC156-12 2.0 1.0 1.5 3.2 2.40 1.34 0.54 0.70											
Unit maintained structural integrity and functionality after the ICC-ES AC 156 test											



#### ATTACHMENT 2: TEST SPECIMEN SUMMARY

ATTACHMENT PAGE | 4 OF 4

UUT- 7	GOC6.6 Console
MANUFACTURER:	General Electric Company
IDENTIFICATION:	Model No.: 5212920-150
DESCRIPTION:	System component of the Discovery CT750 HD System
MOUNTING:	Floor mounted using (4) – 3/8" dia. socket head bolts through GEHC supplied mounting brackets to aluminum interface plate.



#### PROPERTIES:

DIMENSIONS (in.)					LOWEST RESONANT FREQUENCY (Hz.)					
Width	Depth	Height		Weight (lb.)	Side-Axis	Fror	nt-Axis	Vertical-Axis		
48.74	46.54 / 54.74	26.7 / 34	.7	405	6.2	1	7.3	22.9		
SHAKE TABLE T	EST PARAMETERS									
CODE	TEST CRITERIA	$S_{DS}$	z/h	I <sub>P</sub>	A <sub>FLX-H</sub>	A <sub>RIG-H</sub>	$A_{\text{FLX-V}}$	$A_{RIG-V}$		
CBC 2013	ICC-ES AC156-12	2.0 2.6	1.0 0.0	1.5	3.2 2.6	2.40 1.04	1.34 1.74	0.54 0.70		
Unit maintained s	tructural integrity and fund			S AC 156 test	2.0	1.04	1.74	0.70		

UUT- 8	Power Distribution Unit
MANUFACTURER:	GE Hangwei Medical Systems CO., LTD
IDENTIFICATION:	Model No.: 2326492-61
DESCRIPTION:	System component of the Optima CT580 System
MOUNTING:	Floor mounted using (4) – 3/8" dia. socket head bolts through GEHC supplied mounting brackets to aluminum interface plate.



#### PROPERTIES:

DIMENSIONS (in.)					LOWEST RESONANT FREQUENCY (Hz.)					
Width	Depth	Height		Weight (lb.)	X-Axis	Y	-Axis	Z-Axis		
27.6	21.7	41.8		796+22 brackets	15.7		19.9	45		
SHAKE TABLE TEST PARAMETERS										
CODE	TEST CRITERIA	$S_{DS}$	z/h	I <sub>P</sub>	$A_{FLX-H}$	$A_{RIG-H}$	$A_{\text{FLX-V}}$	$A_{RIG-V}$		
T CBC 2013 T ICC-ES ΔC156-12 T T T T T T T T T T T T T T T T T T T								0.54 0.70		
Unit maintained s	Unit maintained structural integrity and functionality after the ICC-ES AC 156 test									