



**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

**APPLICATION FOR OSHPD SPECIAL SEISMIC
CERTIFICATION PREAPPROVAL (OSP)**

OFFICE USE ONLY	
APPLICATION #:	OSP – 0614

OSHPD Special Seismic Certification Preapproval (OSP)

Type: New Renewal

Manufacturer Information

Manufacturer: Philips Medical Systems DMC GmbH

Manufacturer's Technical Representative: Ming Xiao

Mailing Address: Roentgenstrasse 24, 22335 Hamburg

Telephone: +49 40 349712306

Email: Ming.Xiao@Philips.com

Product Information

Product Name: DigitalDiagnost

Product Type: Radiography and Fluoroscopy Medical Imaging System

Product Model Number: See Attachment

(List all unique product identification numbers and/or part numbers)

General Description: Multiple component digital radiography and fluoroscopy medical diagnostic imaging system.

Seismic enhancements incorporated into the test units and enhancements required to address anomalies observed during the tests shall be incorporated into the certified units.

Mounting Description: Rigid floor mounted, rigid wall mounted, ceiling mounted, and combined floor / wall mounted, see attachment.

Applicant Information

Applicant Company Name: W.E. Gundy & Associates, Inc.


Contact Person: Travis Soppe, SE

Mailing Address: 1199 Shoreline Drive, Suite 310, Boise, ID 83702

Telephone: (208) 342-5898 Ext. 115

Email: tsoppe@wegai.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2016.


Signature of Applicant: 

Date: 06-13-2019

Title: President

Company Name: W.E. Gundy & Associates, Inc.

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

California Licensed Structural Engineer Responsible for the Engineering and Test Report(s)

Company Name: W.E. Gundy & Associates, Inc.

Name: Travis Soppe, SE California License Number: S6115

Mailing Address: 1199 Shoreline Drive, Suite 310, Boise, ID 83702

Telephone: (208) 342-5898 Ext. 115 Email: tsoppe@wegai.com

Supports and Attachments Preapproval


Supports and attachments are preapproved under OPM-
(Separate application for OSHPD Preapproval of Manufacturer's Certification (OPM) of Supports and attachments is required)

Supports and attachments are not preapproved

Certification Method

Testing in accordance with: ICC-ES AC156

Other (Please Specify): _____




Testing Laboratory

Company Name: IABG mbH

Contact Name: Dr. Steffen Roedling

Mailing Address: Einsteinstrasse 20, Ottobrunn, Germany D-85521

Telephone: +49 (0) 89 / 6088-2052 Email: roedling@iabg.de





OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

Seismic Parameters

Design in accordance with ASCE 7-10 Chapter 13: [X] Yes [] No

Design Basis of Equipment or Components (Fp/Wp) = See Attachment

SDS (Design spectral response acceleration at short period, g) = 2.00, z/h = 1.0 and 2.50, z/h = 0

ap (In-structure equipment or component amplification factor) = Multiple, see attachment

Rp (Equipment or component response modification factor) = Multiple, see attachment

Omega_0 (System overstrength factor) = Multiple, see attachment

Ip (Importance factor) = 1.5

z/h (Height factor ratio) = 1, SDS = 2.00 and 0 at SDS = 2.50

Equipment or Component Natural Frequencies (Hz) = Multiple, see attachment

Overall dimensions and weight (or range thereof) = Multiple, see attachment

Equipment or Components @ grade designed in accordance with ASCE 7-10 Chapter 15: [] Yes [X] No

Design Basis of Equipment or Components (V/W) =

SDS (Design spectral response acceleration at short period, g) =

SD1 (Design spectral response acceleration at 1 second period, g) =

R (Response modification coefficient) =

Omega_0 (System overstrength factor) =

Cd (Deflection amplification factor) =

Ip (Importance factor) = 1.5

Height to Center of Gravity above base =

Equipment or Component Natural Frequencies (Hz) =

Overall dimensions and weight (or range thereof) =

Tank(s) designed in accordance with ASME BPVC, 2015: [] Yes [X] No

List of Attachments Supporting Special Seismic Certification

[X] Test Report(s) [] Drawings [X] Calculations [X] Manufacturer's Catalog

[X] Other(s) (Please Specify): Certified System Matrix, UUT Summary Sheets, Subcomponent Certification Letter

OSHPD Approval (For Office Use Only) - Approval Expires on December 31, 2025

Signature: [Signature] Date: February 18, 2021

Print Name: Timothy J. Piland Title: SSE

Special Seismic Certification Valid Up to: SDS (g) = See Above z/h = See Above

Condition of Approval (if applicable):




Table 1	PHILIPS MEDICAL SYSTEMS DMC GmbH SPECIAL SEISMIC CERTIFICATION CERTIFIED SYSTEM AND COMPONENTS					 WEGAI <small>W.E. GUNDY & ASSOCIATES, INC. STRUCTURAL & EARTHQUAKE ENGINEERING</small>		

System: DigitalDiagnost Radiography System				Manufacturer: Philips Medical Systems DMC GmbH			
---	--	--	--	---	--	--	--

System Component ¹	Philips Part Number	Dimensions (in)			Weight (lb)	Mounting	UUT ⁴
		Width	Length	Height			

Ceiling Suspension CSM Product Family

Ceiling Suspension CSM3 (Long) Version: Comfort Move	9890-010-87703	31.0	159.0	45.5-110.7	577	ceiling	UUT _w -1
---	----------------	------	-------	------------	-----	---------	---------------------

Ceiling Suspension CSM3 (Short) Version: Comfort Move	9890-010-87703	31.0	92.8	45.5-110.7	532	ceiling	UUT _w -2
--	----------------	------	------	------------	-----	---------	---------------------

Digital Table TH Product Family

Table TH2 Wifi Detector	9848-600-02222	33.5	94.5	20.8-36.4	612 ²	floor	UUT _y -1
----------------------------	----------------	------	------	-----------	------------------	-------	---------------------

Table TH2 Fixed Detector	9848-600-02222	33.5	94.5	20.8-36.4	615 ²	floor	UUT _y -2
-----------------------------	----------------	------	------	-----------	------------------	-------	---------------------

Table TH2 Wifi Detector	9848-600-02222	33.5	94.5	20.8-36.4	603 ²	floor	UUT _w -3
----------------------------	----------------	------	------	-----------	------------------	-------	---------------------

Vertical Stand VS Product Family

Vertical Stand 2 VS2 Non-Tilting - Fixed Detector	9897-010-01411	28.3	60.6	82.7	326	floor/wall	UUT _x -2A
--	----------------	------	------	------	-----	------------	----------------------

Vertical Stand 2 VS2 Non-Tilting - Wifi Detector	9897-010-01411	28.3	60.6	82.7	326	floor/wall	UUT _x -2B
---	----------------	------	------	------	-----	------------	----------------------

Vertical Stand 2 VS2 Tilting - Fixed Detector	9897-010-01411	23.5	35.5	81.9	562	floor/wall	interpolated
--	----------------	------	------	------	-----	------------	--------------

Vertical Stand 2 VS2 Tilting - Wifi Detector	9897-010-01411	23.5	35.5	81.9	562	floor/wall	UUT _z -8
---	----------------	------	------	------	-----	------------	---------------------

¹ All components are manufactured by Philips Medical Systems DMC GmbH unless noted in (). The part numbers listed uniquely identify the type of component, manufacturer, and material of construction for each sub-component within the tested units.

² Table TH2 weights listed do not include 560lb UUT_y-1/2 and 660lb UUT_w-3 simulated patient weight used during the seismic tests.

⁴ The units were tested at different times and the subscripts on the UUT's reference the following seismic certification test reports:
w = TAB3-PB-19-020-V1 x = TAF4-PB-18-082-V1 y = SQ13-1703-01 z = SQ13-1701-01

SEISMIC CERTIFICATION LIMITS

System Component	Code	S _{DS} (g)	z / h	I _P	a _P	R _P	Ω ₀	F _P / W _P
Ceiling Suspension CSM	CBC 2016	2.0	1.0	1.50	2.5	2.5	2.0	3.60
		2.5	0					1.50
Digital Table TH	CBC 2016	2.0	1.0	1.50	1.0	1.5	1.5	2.40
		2.5	0					1.13
Vertical Stand VS	CBC 2016	2.0	1.0	1.50	1.0	1.5	1.5	2.40
		2.5	0					1.13

Table 1	PHILIPS MEDICAL SYSTEMS DMC GmbH SPECIAL SEISMIC CERTIFICATION CERTIFIED SYSTEM AND COMPONENTS					 WEGAI <small>W.E. GUNDY & ASSOCIATES, INC. STRUCTURAL & EARTHQUAKE ENGINEERING</small>	
	System: DigitalDiagnost Radiography System			Manufacturer: Philips Medical Systems DMC GmbH			

System Component ¹	Philips Part Number	Dimensions (in)			Weight (lb)	Mounting	UUT ⁴
		Width	Length	Height			
Power Distribution							
M-Cabinet CXA	9890-010-88393	19.2	27.6	39.0	236	floor	UUT _y -4

User Interface and Miscellaneous							
AWS-2X PC (ADVANTECH)	4512-201-11781	8.5	10.9	17.6	33	floor	UUT _w -4
WIFI Access Point (ARUBA)	4512 201 12191	5.9	5.9	2.8	1	wall	UUT _x -4
WIFI Access Point (ARUBA)	4512 201 1219a ³ 4512 201 1220a ³	5.9	5.9	2.8	1	wall	same

¹ All components are manufactured by Philips Medical Systems DMC GmbH unless noted in (). The part numbers listed uniquely identify the type of component, manufacturer, and material of construction for each sub-component within the tested units.

³ The "a" designation indicates differences in color and software for the identified access point configuration.

⁴ The units were tested at different times and the subscripts on the UUT's reference the following seismic certification test reports:
w = TAB3-PB-19-020-V1 x = TAF4-PB-18-082-V1 y = SQ13-1703-01 z = SQ13-1701-01

DATE: 02/18/2021

SEISMIC CERTIFICATION LIMITS								
System Component	Code	S _{DS} (g)	z / h	I _P	a _P	R _P	Ω ₀	F _P / W _P
M-Cabinet CXA	CBC 2016	2.0	1.0	1.50	2.5	6.0	2.0	1.50
		2.5	0					1.13
AWS-2X PC	CBC 2016	2.0	1.0	1.50	1.0	2.5	2.0	1.44
		2.5	0					1.13
WIFI Access Point	CBC 2016	2.0	1.0	1.50	1.0	2.5	2.0	1.44
		2.5	0					1.13

Table 2	PHILIPS MEDICAL SYSTEMS DMC GmbH SPECIAL SEISMIC CERTIFICATION CERTIFIED SUBCOMPONENTS	 WEGAI <small>W.E. GUNDY & ASSOCIATES, INC. STRUCTURAL & EARTHQUAKE ENGINEERING</small>
----------------	---	---

System: DigitalDiagnost Radiography System	Manufacturer: Philips Medical Systems DMC GmbH
---	---

Subcomponent	Philips Part Number ¹	Dimensions (in)			Weight (lb)	UUT ²
		Width	Length	Height		

Vertical Stand VS Subcomponent Matrix

Portable Detector Skyplate Small (pixium 2430 EZ)	9897 010 00843	9.5	1.2	1.0	4	UUT _x -2B
	9897 010 0084a	9.5	1.2	1.0	4	interpolated
	4598 013 7483a	9.5	1.2	1.0	4	interpolated
	4512 134 7330a	9.5	1.2	1.0	4	interpolated
Portable Detector Skyplate Large (pixium 3543 EZ)	4598 012 4298a	13.8	16.9	1.0	7	interpolated
	4512 134 7214a	13.8	16.9	1.0	7	interpolated
	9897 010 0083a	13.8	16.9	1.0	7	interpolated
	9897 010 0268a	13.8	16.9	1.0	7	interpolated
	9897 010 02681	13.8	16.9	1.0	7	UUT _z -8
Fixed Detector (pixium 4343RCE)	4598 009 40141	20.0	19.4	2.0	44	UUT _x -2A
	4598 009 4014a	20.0	19.4	2.0	44	same
	4598 009 4015a	20.0	19.4	2.0	44	same

Digital Table TH / ProxiDiagnost Stand Subcomponent Matrix

Portable Detector Skyplate Large (pixium 3543 EZ)	4598-012-42982	13.8	16.9	1.0	7	UUT _w -3
	4598-012-4298a	13.8	16.9	1.0	7	interpolated
	4598-134-7214a	13.8	16.9	1.0	7	interpolated
	9897 010 0268a	13.8	16.9	1.0	7	interpolated
	9897 010 02681	13.8	16.9	1.0	7	UUT _y -1
Fixed Detector (pixium 4343RCE)	4598-009-40141	20.0	19.4	2.0	44	UUT _y -2
	4598-009-4014a	20.0	19.4	2.0	44	same
	4598-009-4015a	20.0	19.4	2.0	44	same

¹ "a" designation in Philips part number indicates differences in color, software, ect. for the identified subcomponent.

² The units were tested at different times and the subscripts on the UUT's reference the following seismic certification test reports:
w = TAB3-PB-19-020-V1 x = TAF4-PB-18-082-V1 y = SQ13-1703-01 z = SQ13-1701-01

UUT_w-1

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rails and connecting parts of the component itself (detailed in Phillips 4512 983 02641) bolt to unistrut grid spaced at 25.6" on center. The unistrut grid consisted of MURPO#150969 MPR-41/82/2.0 H-Profiles (Unistrut P1001 equivalent) anchored with 2 - M10 bolts with clamping claws (MURPO 157219) at each intersection to the ceiling fixture framing spaced at 23.6" on center.



Manufacturer: Philips Medical Systems DMC GmbH		
Component: Ceiling Suspension CSM3 (Long) Version: Comfort Move		
Model / Serial Number: 9890-010-87703 / 19 000005		
UUT Function: Digital system for making X-ray exposures of the body.		
UUT Description: DigitalDiagnost ceiling suspension system - long rail configuration		
Test Location: IABG mbH, Germany	Test Date: February 2019	Report: TAB3-PB-19-020-V1

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
577	31.0"	159.0"	45.5" - 110.7"	NA	NA	NA

The CS system moves laterally, rotates, and extends up and down to accommodate different patients and procedures. The system was tested in the normal operating position with the system horizontally centered, no rotation, and with an extension of 26.4" making the tested height = 71.9".

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _P	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_w-2

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rails and connecting parts of the component itself (detailed in Phillips 4512 983 02641) bolt to unistrut grid spaced at 25.6" on center. The unistrut grid consisted of MURPO#150969 MPR-41/82/2.0 H-Profiles (Unistrut P1001 equivalent) anchored with 2 - M10 bolts with clamping claws (MURPO 157219) at each intersection to the ceiling fixture framing spaced at 23.6" on center.



Manufacturer: Philips Medical Systems DMC GmbH

Component: Ceiling Suspension CSM3 (Short) Version: Comfort Move

Model / Serial Number: 9890-010-87703 / 19 000001

UUT Function: Digital system for making X-ray exposures of the body.

UUT Description: DigitalDiagnost ceiling suspension system - short rail configuration

Test Location: IABG mbH, Germany **Test Date:** February 2019 **Report:** TAB3-PB-19-020-V1

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
532	31.0"	92.8"	45.5" - 110.7"	NA	NA	NA

The CS system moves laterally, rotates, and extends up and down to accommodate different patients and procedures. The system was tested in the normal operating position with the system horizontally centered, no rotation, and with an extension of 26.4" making the tested height = 71.9".

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _p	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_w-3

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rigid floor mounted with 4 - 1/2" grade 8 bolts torqued to 110ft-lbs with 2 x 1 x 0.625" plate washers.



REVIEWED FOR CODE COMPLIANCE
DATE: 02/18/2021
BY: Timothy J Piland
OSP-0614
IABG

Manufacturer: Philips Medical Systems DMC GmbH		
Component: Table TH2 - Wifi Detector	Model / Serial Number: 9848-600-02222 / 19020008	
UUT Function: Motorized patient table and detector system used for X-Ray imaging.		
UUT Description: DigitalDiagnost patient table with installed detector: Wifi Pixium 3543 EZ - 4598-012-42982		
Test Location: IABG mbH, Germany	Test Date: February 2019	Report: TAB3-PB-19-020-V1

UUT PROPERTIES

Weight (lb) with Patient	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
1,263	33.5"	94.5"	20.8 - 36.4"	4.0	3.0	7.3

The patient table moves laterally both ways and vertically to accommodate different patients and procedures. The system was tested in the normal operating position with the table horizontally centered, a table top height of 29.3", and a total simulated patient weight of 660lbs.

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _P	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_w-4

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rigid floor mounting using Phillips provided seismic restraint kit 4512-134-74741. Seismic restraint kit include a hand tightend cam buckle strap (450lb WLL) looped thru angle brackets on each side of the unit. The angle brackets are attached to the table with individual 3/8" grade 8 bolts.



Manufacturer: ADVANTECH

Component: AWS-2X PC User Interface | **Model / Serial Number:** 4512-201-11781 / KSE0149757

UUT Function: Computational processing for image system.

UUT Description: DigitalDiagnost user interface PC

Test Location: IABG mbH, Germany | **Test Date:** February 2019 | **Report:** TAB3-PB-19-020-V1

UUT PROPERTIES

Weight (lb) with Patient	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
33	8.5"	10.9"	17.6"	> 33	28.1	>33

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _P	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67

Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_x-2A

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Floor/wall mount with 2 - M10 grade 8.8 bolts / washers torqued at 36ft-lb and 1 - M12 grade 8.8 bolt / washer torqued at 62ft-lb to the floor and 2 - M10 grade 10.9 bolts / washers torqued at 53ft-lb to the wall fixture. Electrical box mounted to wall fixture with 4 - M6 grade 8.8 bolts / washer torqued to 10ft-lb.



Manufacturer: Philips Medical Systems DMC GmbH

Component: Vertical Stand VS2 - Fixed Detector | **Model / Serial Number:** 9897 010 01411

UUT Function: Vertical X-Ray imaging system

UUT Description: Component of DigitalDiagnost Radiography and Fluoroscopy System with fixed detector (Pixium TX4343RCE - 4598 009 40141) and electrical box (4598-012-06871).

Test Location: IABG mbH, Germany | **Test Date:** February 2018 | **Report:** TAF4-PB-18-082-V1

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
326	28.3"	60.6"	82.7"	> 33	28.9	> 33

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _P	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_x-2B

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Floor/wall mount with 2 - M10 grade 8.8 bolts / washers torqued at 36ft-lb and 1 - M12 grade 8.8 bolt / washer torqued at 62ft-lb to the floor and 2 - M10 grade 10.9 bolts / washers torqued at 53ft-lb to the wall fixture. Electrical box mounted to wall fixture with 4 - M6 grade 8.8 bolts / washer torqued to 10ft-lb.



Manufacturer: Philips Medical Systems DMC GmbH

Component: Vertical Stand VS2 - Wifi Detector

Model / Serial Number: 9897 010 01411

UUT Function: Vertical X-Ray imaging system

UUT Description: Component of DigitalDiagnost Radiography and Fluoroscopy System with wifi detector (Pixium 2430EZ - 9897 010 00843) and electrical box (4598-012-06871).

Test Location: IABG mbH, Germany

Test Date: February 2018

Report: TAF4-PB-18-082-V1

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
326	28.3"	60.6"	82.7"	> 33	29.0	> 33

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _P	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_x-4

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Wall mount with 4 - size 3 wood screws.



Manufacturer: ARUBA

Component: WIFI Access Point

Model / Serial Number: 4512 201 12191

UUT Function: WIFI Module

UUT Description: Component of ProxiDiagnost Radiography and Fluoroscopy System

Test Location: IABG mbH, Germany

Test Date: February 2018

Report: TAF4-PB-18-082-V1

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
1	5.9"	5.9"	2.8"	NA	NA	NA

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{Ds} (g)	z / h	I _p	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_y-1

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rigid floor mounted with 4 - 1/2" grade 8 bolts torqued at 40ft-lb with 3.5" x 1.5" x 1/2" plate washers.



Manufacturer: Philips Medical Systems DMC GmbH

Component: Table TH2 - Wifi Detector | **Model / Serial Number:** 9848-600-02222 / 16000193

UUT Function: Motorized patient table and detector system used for X-Ray imaging.

UUT Description: DigitalDiagnost patient table with installed detector: Wifi Pixium 3543 EZ - 9897-010-02681.

Test Location: Environmental Testing Laboratory, Dallas, Texas

Test Date: January 2017

Report: SQ13-1703-01

UUT PROPERTIES

Weight (lb) with Patient	Dimensions (inches)			Natural Fequency (Hz)		
	Width	Depth	Height	FB	SS	V
1,172	33.5"	94.5"	20.8 - 36.4"	3.7	20.4	7.2

The patient table moves laterally both ways and vertically to accommodate different patients and procedures. The system was tested in the normal operating position with the table horizontally centered, a table top height of 30", and a total simulated patient weight of 560lbs.

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _p	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67

Note: The unit was full of contents during testing and remained fuctional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_y-2

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rigid floor mounted with 4 - 1/2" grade 8 bolts torqued at 40ft-lb with 3.5" x 1.5" x 1/2" plate washers.



Manufacturer: Philips Medical Systems DMC GmbH

Component: Table TH2 - Fixed Detector | **Model / Serial Number:** 9848-600-02222 / 16020429

UUT Function: Motorized patient table and detector system used for X-Ray imaging.

UUT Description: DigitalDiagnost patient table with fixed detector: Pixium 4343RCE 4598-009-40141

Test Location: Environmental Testing Laboratory, Dallas, Texas

Test Date: January 2017

Report: SQ13-1703-01

UUT PROPERTIES

Weight (lb) with Patient	Dimensions (inches)			Natural Fequency (Hz)		
	Width	Depth	Height	FB	SS	V
1,175	33.5"	94.5"	20.8 - 36.4"	3.2	12	7.2

The patient table moves laterally both ways and vertically to accommodate different patients and procedures. The system was tested in the normal operating position with the table horizontally centered, a table top height of 30", and a total simulated patient weight of 560lbs.

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _P	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67


Note: The unit was full of contents during testing and remained fuctional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_y-4

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Mounting Details: Rigid floor mounted with 4 - 3/8" A574 bolts with washers



Manufacturer: Philips Medical Systems DMC GmbH

Component: M-Cabinet CXA

Model / Serial Number: 9890-010-88393 / 16230001

UUT Function: X-Ray generator for Radiography System

UUT Description: Component of DigitalDiagnost Radiography and Fluoroscopy System

Test Location: Environmental Testing Laboratory, Dallas, Texas

Test Date: January 2017

Report: SQ13-1703-01

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Frequency (Hz)		
	Width	Depth	Height	FB	SS	V
236	19.2"	27.6"	39"	10.8	9.0	23.2

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _p	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67

Note: The unit was full of contents during testing and remained functional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

UUT_z-8

**UNIT UNDER TEST (UUT)
SUMMARY SHEET**



Note: Modifications required for unit to pass seismic test. All modifications incorporated in Phillips seismic kit 4512-201-12241.

Mounting Details: Rigid floor / wall mounted with 3 - 1/2" grade 5 bolts torqued at 40ft-lbs to the floor and 4 - 3/8" grade 2 bolts torqued at 17ft-lbs to 16 ga. steel backing plate on the wall



Manufacturer: Philips Medical Systems DMC GmbH

Component: Vertical Stand 2 VS2 - Wifi Detector **Model / Serial Number:** 9897-010-01411 / 16000012

UUT Function: Vertical X-Ray imaging system.

UUT Description: Component of DigitalDiagnost Radiography and Fluoroscopy System with wifi detector (Pixium 3543EZ - 9897 010 02681) and electrical box (4598-012-06641).

Test Location: Environmental Testing Laboratory, Dallas, Texas

Test Date: January 2017

Report: SQ13-1701-01

UUT PROPERTIES

Weight (lb)	Dimensions (inches)			Natural Fequency (Hz)		
	Width	Depth	Height	FB	SS	V
562	23.5"	35.5"	81.9"	10.2	4.2	6.1

SEISMIC TEST PARAMETERS

Building Code / Test Criteria	S _{DS} (g)	z / h	I _p	A _{FLX-H} (g)	A _{RIG-H} (g)	A _{FLX-V} (g)	A _{RIG-V} (g)
CBC 2019 / ICC-ES AC156	2.00	1.0	1.5	3.20	2.40		
	2.50	0.0	1.5			1.67	0.67

Note: The unit was full of contents during testing and remained fuctional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.