



**OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION**

**APPLICATION FOR OSHPD PREAPPROVAL OF
MANUFACTURER'S CERTIFICATION (OPM)**

OFFICE USE ONLY

APPLICATION #: OPM-0369

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: New Renewal/Update

Manufacturer Information

Manufacturer: Omnicell, Inc.

Manufacturer's Technical Representative: Todd Kijowski

Mailing Address: 51 Pennwood Place, Suite 400, Warrendale, PA 15086

Telephone: (724) 741-7777 Email: Todd.Kijowski@omnicell.com

Product Information

Product Name: HALF-HEIGHT CABINETS

Product Type: Automated medication Dispensing Cabinets

Product Model Number: MED-FRM-102, -103, -104, -020, -021, -029, -039

General Description: Medication storage and dispensing cabinets

Applicant Information

Applicant Company Name: Omincell, Inc

Contact Person: Todd Kijowski

Mailing Address: 51 Pennwood Place, Suite 400, Warrendale, PA 15086

Telephone: (724) 741-7777 Email: todd.Kkjowski@omnicell.com

Title: _____

Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY





OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations

Company Name: DEGENKOLB ENGINEERS

Name: Chad Closs

California License Number: S5946

Mailing Address: 225 Broadway, Suite 1325, San Diego, CA 92101

Telephone: (858) 699-5412

Email: ccloss@degenkolb.com

OSHPD Special Seismic Certification Preapproval (OSP)

Special Seismic Certification is preapproved under OSP

OSP Number: _____

Certification Method

Testing in accordance with: ICC-ES AC156 FM 1950-16

Other(s) (Please Specify): _____

*Use of criteria other than those adopted by the California Building Standards Code, 2019 (CBSC 2019) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2019 may be used when approved by OSHPD prior to testing.

Analysis

Experience Data

Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

OSHPD Approval

Date: 10/1/2021

Name: Kamalpreet Kalsi

Title: Senior Structural Engineer

Condition of Approval (if applicable): _____





OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION

OPM-0369-19

OMNICELL HALF-HEIGHT CABINET

DEGENKOLB ENGINEERS
225 Broadway, Suite 1325
San Diego, CA 92101
619.515.0299 PHONE
619.515.0298 FAX



MODEL NUMBERS

MED-FRM-102, MED-FRM-103, MED-FRM-104, MED-FRM-020, MED-FRM-021, MED-FRM-029, MED-FRM-039

GENERAL NOTES:

1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2019 CALIFORNIA BUILDING CODE (CBC). THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2019.
2. PRE-APPROVED DESIGN AND MATERIALS CONFORM WITH THE 2019 EDITION OF THE CALIFORNIA BUILDING CODE. DETAILS WITHIN THIS APPROVAL MAY BE USED ANYWHERE IN THE STATE OF CALIFORNIA WHERE $S_{DS} \leq 2.50$ FOR CASE 1 & 2, $S_{DS} \leq 1.50$ ($z/h \leq 1.0$) OR $S_{DS} \leq 2.25$ ($z/h \leq 0.5$) FOR CASE 3.
3. SEISMIC FORCES ON EQUIPMENT DETERMINED PER THE 2019 CBC & ASCE 7-16 SECTION 13.3. ALL LOADS IN THIS PRE-APPROVAL ARE AT STRENGTH LEVEL AND SHALL BE USED FOR STRENGTH DESIGN.
 - a. CASE 1 (EQUIPMENT ABOVE GRADE TO ROOF):
 $S_{DS}=2.50$, $a_p=1.0$, $R_p=1.5$, $I_p=1.5$, $\Omega_o=1.5$, $z/h \leq 1.0$
i. $F_p=3.00W_p$, $F_v=0.50W_p$
 - b. CASE 2 (EQUIPMENT AT OR BELOW GRADE, EXPANSION ANCHOR OPTION):
 $S_{DS}=2.50$, $a_p=1.0$, $R_p=1.5$, $I_p=1.5$, $\Omega_o=1.5$, $z/h \leq 0.0$
i. $F_p=1.13W_p$, $F_v=0.50W_p$
 - c. CASE 3 (EQUIPMENT ABOVE GRADE TO ROOF, EXPANSION ANCHOR OPTION):
 $S_{DS} \leq 1.50$, $a_p=1.0$, $R_p=1.5$, $I_p=1.5$, $\Omega_o=1.5$, $z/h \leq 1.0$
i. $F_p=1.80W_p$, $F_v=0.30W_p$
 $S_{DS} \leq 2.25$, $a_p=1.0$, $R_p=1.5$, $I_p=1.5$, $\Omega_o=1.5$, $z/h \leq 0.5$
i. $F_p=1.80W_p$, $F_v=0.45W_p$
4. THE STRUCTURAL ENGINEER-OF-RECORD (S.E.O.R.) IS RESPONSIBLE FOR THE FOLLOWING:
 - a. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB OPENINGS OR EDGES.
 - b. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS.
 - c. DESIGN ANY SUPPLEMENTARY MEMBERS AND THEIR ATTACHMENTS WHICH THE UNIT IS ANCHORED TO. VERIFY THE ADEQUACY OF ANY EXISTING MEMBERS AND THEIR ATTACHMENTS WHICH THE UNIT IS ANCHORED TO FOR THE FORCES EXERTED ON THEM BY THE UNIT IN ADDITION TO ALL OTHER LOADS AND FORCES.
 - d. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2019 CBC AND WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL. VERIFY THAT THE EQUIPMENT'S ACTUAL WEIGHT, CG LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS AND THE MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE FORMATION SHOWN IN THIS PRE-APPROVAL.

5. STRUCTURAL ENGINEER-OF-RECORD (S.E.O.R.) SHALL EVALUATE BRACKET ATTACHMENTS FOR CONDITIONS THAT VARY FROM THIS PRE-APPROVAL.
6. CONTRACTOR/INSPECTOR OF RECORD MUST VERIFY ANCHOR SPACING TO EXISTING ADJACENT ANCHORS IS TO BE GREATER THAN 8".
7. THIS OPM COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE UNIT TO THE STRUCTURE
8. EXPANSION OR WEDGE ANCHORS INTO CONCRETE: HILTI KB-TZ2 (ICC ESR 4266) AND DEWALT POWER-STUD+ SD2 (ICC ESR 2502). INSTALL ANCHORS IN ACCORDANCE WITH THE ICC REPORT AND MANUFACTURER'S RECOMMENDATIONS. TEST AT LEAST 50% OF ANCHORS NO SOONER THAN 24 HOURS AFTER INSTALLATIONS. TESTS SHALL BE CONDUCTED IN THE PRESENCE OF THE INSPECTOR OF RECORD (IOR) AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO OSHPD.

TEST PER THE FOLLOWING METHOD:

- a. TORQUE WRENCH TEST: TEST ANCHORS TO THE REQUIRED TORQUE LOAD GIVEN IN TABLE BELOW WITHIN THE LIMIT OF ONE-HALF TURN OF THE NUT.

ANCHOR TEST LOAD VALUES								
ANCHOR TYPE ^a	ANCHOR DIAMETER	EMBED h_{ef}	TORQUE LOAD (FT-LBS)	CONCRETE TYPE	f_c MIN (PSI)	MINIMUM SPACING	MINIMUM EDGE DIST. REQ.	ICC-ES ESR NO.
HILTI KB-TZ2	3/8"	2"	30	SAND LWC OR NWC	3,000	5"	36"	4266
DEWALT SD2	3/8"	2"	20	SAND LWC OR NWC	3,000	5"	36"	2502

- a. PROVIDE FOR FULL ENGAGEMENT OF NUT & WASHER
9. IF ANY ANCHOR FAILS DURING TESTING, UNIT MUST BE MOVED SO THAT NO ANCHOR IS WITHIN 8" OF AN ABANDONED ANCHOR.

10. A MANUFACTURER PROVIDED PERMANENT PLAQUE MUST BE AFFIXED ON THE UNIT STATING THE FOLLOWING: "WEIGHT OF CONTENTS SHALL NOT EXCEED 10 PCF". DESIGNED WEIGHT OF CONTENTS IS 10 PCF. VERIFY IN FIELD BEFORE INSTALLATION.
11. FOR BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT CONDITION (SNUG TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) IS ACHIEVED.
 - B. THROUGH BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.
12. INSTALLATION PROCEDURE:
 - a. MOUNT BASE ANGLE PROVIDED BY OMNICELL TO FLOOR WITH AS SHOWN IN THIS OPM.
 - b. POSITION UNIT WITH RESPECT TO BASE ANGLES. DOWEL INTO UNIT AS SHOWN. REFERENCE MFR'S INSTALLATION GUIDE FOR PROPER SPACING



OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION

OPM-0369-19


OMNICELL HALF-HEIGHT CABINET

DEGENKOLB ENGINEERS
225 Broadway, Suite 1325
San Diego, CA 92101
619.515.0299 PHONE
619.515.0298 FAX



MODEL NUMBERS
MED-FRM-102, MED-FRM-103, MED-FRM-104, MED-FRM-020, MED-FRM-021, MED-FRM-029, MED-FRM-039

CASE 1 - HALF-HEIGHT CABINETS ABOVE GRADE




MAX. DRILLED HOLE SIZE = BOLT DIAM. PLUS 1/16"

FRONT ELEVATION


SEE "MINIMUM STEEL DECK REQUIREMENTS" DETAIL FOR SLAB PROPERTIES

3/8" Ø A325 THROUGH BOLTS IN STANDARD HOLE AT BASE BRACKET AND ANGLE BELOW SLAB, (8) TOTAL PER BRACKET



PLAN SECTION C-C

BASE BRACKET SEE PAGE 5 FOR PROPERTIES



SIDE ELEVATION

L4x4x1/4
ASTM A36,
TYP.

	FORCES										CABINET PROPERTIES		
	SW _c (LBS)	W _p (LBS)	Rult 1 (LBS)	Vu PIN 1 (LBS)	Rult 2 (LBS)	Rult 3 (LBS)	Vult (LBS/BOLT)	Tult (LBS/BOLT)	Ω _o Vult (LBS/BOLT)	Ω _o Tult (LBS/BOLT)	L (in)	W (in)	H (in)
HALF-HEIGHT CABINETS	800	1,015	1,920	1,523	1,280	1,008	406	809	609	1,213	26 1/2	26 1/8	53

$F_p = 3.00 W_p$ [$S_{DS} \leq 2.50, I_p = 1.5, R_p = 1.5, a_p = 1.0, \Omega_o = 1.5, z/h < 1.0$]

$F_v = 0.50 W_p$

SW_c = SELF-WEIGHT OF THE CABINET

W_p = TOTAL WEIGHT; INCLUDES SW_c AND 10 pcf CONTENTS PER NOTE 10 ON PAGE 1

Rult = MAXIMUM BRACKET PIN UPLIFT FORCE AT STRENGTH LEVEL

Vu PIN1 = MAXIMUM SHEAR ON PIN 1

Vult = MAXIMUM SHEAR PER EXPANSION ANCHOR OR THROUGH BOLT AT STRENGTH LEVEL


Tult = MAXIMUM TENSION FORCE PER EXPANSION ANCHOR OR THROUGH BOLT AT STRENGTH LEVEL

SEE PAGE 5 OF 6 FOR FORCE VECTORS


NOTES: 0369

1. THE DESIGN OF SUPPORTS AND ATTACHMENTS CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE.
2. Rult, Vult AND Tult GIVEN ARE FACTORED LOADS AT STRENGTH LEVEL. FINAL DEMAND FORCES FOR ANCHORAGE TO CONCRETE SHALL INCLUDE OVERSTRENGTH FACTOR Ω_o AS DEFINED BY ASCE 7-16.
3. SEE GENERAL NOTES SECTION ON PAGE 1.
4. FOR THE SUPPORT AND ATTACHMENT DESIGN, THE MOST CRITICAL LOAD COMBINATION IS (0.9 - 0.2S_{DS})xDL
5. SEE PAGE 5 FOR LOCATION OF APPLIED FORCES IN BASE BRACKET.
6. SEE PAGE 5 AND PAGE 6 FOR MANUFACTURER BRACKET INFORMATION.
7. S.E.O.R. MAY RECALCULATE MAX. ANCHOR FORCES Rult, Vult AND Tult, AT THEIR DISCRETION, BASED ON PROJECT SPECIFIC SEISMIC DEMANDS SUBJECT TO OSHPD REVIEW AND APPROVAL.
8. TOTAL WEIGHT (W_p) IS A MAXIMUM. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM SHOWN.
9. EQUIPMENT MANUFACTURER MUST DESIGN UNIT TO MAKE Hcg EQUAL OR LESS THAN THE HEIGHT DIMENSION SHOWN.

MINIMUM STEEL DECK REQUIREMENTS



PROVIDE (2) EXPANSION ANCHORS PER REQUIREMENTS IN NOTE 2



PROVIDE (2) EXPANSION ANCHORS PER REQUIREMENTS IN NOTE 2

MINIMUM STEEL DECK REQUIREMENTS NOTES:



1. PROVIDE 36" MINIMUM DISTANCE TO EDGE OF SLAB OR OPENINGS.
2. PROVIDE (2) 3/8" Ø HILTI KB-TZ2 W/ 2" EMBED EXPANSION ANCHORS TO SUPPORT ANGLE. INSTALL ON THE SLAB RIB INDEPENDENT FROM THROUGH BOLTS. EXTEND ANGLE AS REQUIRED. DO NOT INSTALL EXPANSION ANCHORS IN SLAB RIBS WHERE THROUGH BOLTS ARE PRESENT.
3. W- STEEL DECK TO BE 20 GAGE MIN.




MODEL NUMBERS
MED-FRM-102, MED-FRM-103, MED-FRM-104, MED-FRM-020, MED-FRM-021, MED-FRM-029, MED-FRM-039

CASE 2 - HALF-HEIGHT CABINETS ON GRADE

	FORCES										CABINET PROPERTIES		
	SWc (LBS)	Wp (LBS)	Rult 1 (LBS)	Vu PIN 1 (LBS)	Rult 2 (LBS)	Rult 3 (LBS)	Vult (LBS/BOLT)	Tult (LBS/BOLT)	Ω _o Vult (LBS/BOLT)	Ω _o Tult (LBS/BOLT)	L (in)	W (in)	H (in)
HALF-HEIGHT CABINETS	800	1,015	593	571	417	315	152	269	228	403	26 1/2	26 1/8	53



$F_p = 1.13 W_p$ [$S_{DS} \leq 2.50$, $I_p = 1.5$, $R_p = 1.5$, $a_p = 1.0$, $\Omega_o = 1.5$, $z/h = 0$]
 $F_v = 0.50 W_p$
 SWc = SELF-WEIGHT OF THE CABINET
 Wp = TOTAL WEIGHT; INCLUDES SWc AND 10 pcf CONTENTS PER NOTE 10 ON PAGE 1
 Rult = MAXIMUM BRACKET PIN UPLIFT FORCE AT STRENGTH LEVEL
 Vu PIN1 = MAXIMUM SHEAR ON PIN 1
 Vult = MAXIMUM SHEAR PER EXPANSION ANCHOR OR THROUGH BOLT AT STRENGTH LEVEL
 Tult = MAXIMUM TENSION FORCE PER EXPANSION ANCHOR OR THROUGH BOLT AT STRENGTH LEVEL
 SEE PAGE 5 OF 6 FOR FORCE VECTORS



PLAN SECTION C-C

NOTES:

1. THE DESIGN OF SUPPORTS AND ATTACHMENTS CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE.
2. Rult, Vult AND Tult GIVEN ARE FACTORED LOADS AT STRENGTH LEVEL. FINAL DEMAND FORCES FOR ANCHORAGE TO CONCRETE SHALL INCLUDE OVERSTRENGTH FACTOR Ω_o AS DEFINED BY ASCE 7-16.
3. FOR THE SUPPORT AND ATTACHMENT DESIGN, THE MOST CRITICAL LOAD COMBINATION IS $(0.9 - 0.2S_{DS}) \times DL$
4. SEE GENERAL NOTES SECTION ON PAGE 1.
5. SEE PAGE 5 FOR LOCATION OF APPLIED FORCES IN BASE BRACKET.
6. SEE PAGE 5 AND PAGE 6 FOR MANUFACTURER BRACKET INFORMATION.
7. S.E.O.R. MAY RECALCULATE MAX. ANCHOR FORCES Rult, Vult AND Tult, AT THEIR DISCRETION, BASED ON PROJECT SPECIFIC SEISMIC DEMANDS SUBJECT TO OSHPD REVIEW AND APPROVAL.
8. TOTAL WEIGHT (Wp) IS A MAXIMUM. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM SHOWN.
9. EQUIPMENT MANUFACTURER MUST DESIGN UNIT TO MAKE Hcg EQUAL OR LESS THAN THE HEIGHT DIMENSION SHOWN.



OSHPD PRE-APPROVAL OF MANUFACTURER'S CERTIFICATION

OPM-0369-19

OMNICELL HALF-HEIGHT CABINET



DEGENKOLB ENGINEERS
225 Broadway, Suite 1325
San Diego, CA 92101
619.515.0299 PHONE
619.515.0298 FAX



MODEL NUMBERS

MED-FRM-102, MED-FRM-103, MED-FRM-104, MED-FRM-020, MED-FRM-021, MED-FRM-029, MED-FRM-039

CASE 3 - HALF-HEIGHT CABINETS ABOVE GRADE




	FORCES										CABINET PROPERTIES				
	SWc (LBS)	Wp (LBS)	Rult 1 (LBS)	Vu PIN 1 (LBS)	Rult 2 (LBS)	Rult 3 (LBS)	Vult (LBS/BOLT)	Tult (LBS/BOLT)	Ω _o Vult (LBS/BOLT)	Ω _o Tult (LBS/BOLT)	L (in)	W (in)	H (in)	MAX z/h	MAX S _{DS}
HALF-HEIGHT CABINETS	800	1,015	969	914	677	514	236	434	354	651	26 1/2	26 1/8	53	1.00	1.50
	800	1,015	1,046	914	715	552	236	455	354	682	26 1/2	26 1/8	53	0.50	2.25

$F_p = 1.80 W_p$, $F_v = 0.30 W_p$, [$S_{DS} \leq 1.50$, $I_p = 1.5$, $R_p = 1.5$, $a_p = 1.0$, $\Omega_o = 1.5$, $z/h \leq 1.0$]
 $F_p = 1.80 W_p$, $F_v = 0.45 W_p$, [$S_{DS} \leq 2.25$, $I_p = 1.5$, $R_p = 1.5$, $a_p = 1.0$, $\Omega_o = 1.5$, $z/h \leq 0.5$]
 SW_c = SELF-WEIGHT OF THE CABINET
 W_p = TOTAL WEIGHT; INCLUDES SW_c AND 10 pcf CONTENTS PER NOTE 10 ON PAGE 1
 Rult = MAXIMUM BRACKET PIN UPLIFT FORCE AT STRENGTH LEVEL
 Vu PIN1 = MAXIMUM SHEAR ON PIN 1
 Vult = MAXIMUM SHEAR PER EXPANSION ANCHOR OR THROUGH BOLT AT STRENGTH LEVEL
 Tult = MAXIMUM TENSION FORCE PER EXPANSION ANCHOR OR THROUGH BOLT AT STRENGTH LEVEL
 SEE PG 5 OF 6 FOR FORCE VECTORS

NOTES:



1. THE DESIGN OF SUPPORTS AND ATTACHMENTS CONFORMS TO THE 2019 CALIFORNIA BUILDING CODE.
2. Rult, Vult AND Tult GIVEN ARE FACTORED LOADS AT STRENGTH LEVEL. FINAL DEMAND FORCES FOR ANCHORAGE TO CONCRETE SHALL INCLUDE OVERSTRENGTH FACTOR Ω_o AS DEFINED BY ASCE 7-16.
3. SEE GENERAL NOTES SECTION ON PAGE 1.
4. FOR THE SUPPORT AND ATTACHMENT DESIGN, THE MOST CRITICAL LOAD COMBINATION IS (0.9 - 0.2S_{DS})xDL
5. SEE PAGE 5 FOR LOCATION OF APPLIED FORCES IN BASE BRACKET.
6. SEE PAGE 5 AND PAGE 6 FOR MANUFACTURER BRACKET INFORMATION.
7. S.E.O.R. MAY RECALCULATE MAX. ANCHOR FORCES Rult, Vult AND Tult, AT THEIR DISCRETION, BASED ON PROJECT SPECIFIC SEISMIC DEMANDS SUBJECT TO OSHPD REVIEW AND APPROVAL.
8. TOTAL WEIGHT (W_p) IS A MAXIMUM. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM SHOWN.
9. EQUIPMENT MANUFACTURER MUST DESIGN UNIT TO MAKE H_{cg} EQUAL OR LESS THAN THE HEIGHT DIMENSION SHOWN.

3/8" Ø POWER-STUD+ SD2 W/ 2" MIN EMBED, (8) TOTAL PER BRACKET



BASE BRACKET SEE PAGE 5 FOR PROPERTIES

MINIMUM STEEL DECK REQUIREMENTS





MINIMUM STEEL DECK REQUIREMENTS NOTES:

1. PROVIDE 36" MINIMUM DISTANCE TO EDGE OF SLAB OR OPENINGS
2. W- STEEL DECK TO BE 20 GAGE MIN.



MODEL NUMBERS

MED-FRM-102, MED-FRM-103, MED-FRM-104, MED-FRM-020, MED-FRM-021, MED-FRM-029, MED-FRM-039



NOTE:
 1. FOR ITEMS NOT NOTED SEE PAGE 6

FORCE DISTRIBUTION IN BASE BRACKET

MODEL NUMBERS

MED-FRM-102, MED-FRM-103, MED-FRM-104, MED-FRM-020, MED-FRM-021, MED-FRM-029, MED-FRM-039

