



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL
OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY

APPLICATION #: OPM-0450-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: ☒ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number: _____

Manufacturer Information

Manufacturer: Germfree Laboratories, Inc.

Manufacturer's Technical Representative: Ronald Lavicott

Mailing Address: 4 Sunshine Blvd. Ormond Beach, FL 32174

Telephone: 386-265-4343 Email: rlavicott@germfree.com

Product Information

Product Name: Horizontal Laminar Airflow Workstation

Product Type: ISO 5 / Class 100 Horizontal Laminar Flow Hood

Product Model Number: BZ-3, BZ-4, BZ-6

General Description: ISO 5 / Class 100 Horizontal Laminar Flow Hood designed for the handling of materials in a sterile working environment.

Applicant Information

Applicant Company Name: Sun Structural Engineering, Inc.

Contact Person: Changhua Sun

Mailing Address: 2091 Las Palmas Dr. Suite D Carlsbad, CA. 92011

Telephone: (760)438-1188 Email: sunengineering@sbcglobal.net

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2016.

Signature of Applicant:  Date: 9-20-2017

Title: President Company Name: Sun Structural Engineering, Inc.

"Access to Safe, Quality Healthcare Environments that Meet California's



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations

Company Name: Sun Structural Engineering, Inc.

Name: Changhua Sun

California License Number: S4609

Mailing Address: 2091 Las Palmas Dr. Suite D Carlsbad, CA. 92011

Telephone: (760)438-1188

Email: sunengineering@sbcglobal.net

OSHPD Special Seismic Certification Preapproval (OSP)

- ☐ Special Seismic Certification is preapproved under OSP-
(Separate application for OSP is required)
- ☒ Special Seismic Certification is not preapproved

Certification Method(s)

- ☐ Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-16
- ☐ Other* (Please Specify): _____

*Use of criteria other than those adopted by the California Building Standards Code, 2016 (CBSC 2016) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2016 may be used when approved by OSHPD prior to testing.

- ☐ Analysis
- ☐ Experience Data
- ☐ Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

List of Attachments Supporting the Manufacturer's Certification

- ☐ Test Report ☒ Drawings ☒ Calculations ☐ Manufacturer's Catalog
- ☐ Other(s) (Please Specify): _____

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2016 & ALL PRE-2016 CODE BASED PROJECTS

Signature: Jeffrey Kikumoto Date: 09-28-2018


Print Name: Jeffrey Kikumoto

Title: SSE

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY
OSH-FD-700 (REV 12/16/15)



OSHPD

Page 2 of 2

OFFICE STATEWIDE HEALTH PLANNING AND DEVELOPMENT

PRE-APPROVAL MANUFACTURER'S CERTIFICATION

OPM - 0450 - 13

THIS PRE-APPROVAL CONFORMS TO THE 2016 CALIFORNIA BUILDING CODE

MANUFACTURER: GERMFREE LABORATORIES

EQUIPMENT TYPE: HORIZONTAL LAMINAR AIRFLOW WORKSTATION; BZ3, BZ4 & BZ6 MODELS

GENERAL NOTES

1. THIS OSHPD PRE APPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2016. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2016.
2. THE DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR A SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THE DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
3. DESIGN CRITERIA AND LIMITATIONS:
 - A. $S_{DS} \leq 1.5$ AND $1.5 < S_{DS} \leq 2.5$
 - B. $a_p = 1.0$
 - C. $R_p = 2.5$
 - D. $I_p = 1.5$
 - E. $z/h = 0$, FOR CONCRETE SLAB ON GRADE
 $z/h \leq 1.0$, FOR SLAB ABOVE GRADE
 - F. $\Omega_o = 2.0$, PER 2016 CBC 1616A.1.23 FOR CONCRETE ANCHORAGE
4. ALL ANCHOR FORCES SHOWN ON THIS DOCUMENT ARE STRENGTH LEVEL LOADS AND MAY BE USED FOR STRENGTH DESIGN.
5. THIS PRE APPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF EQUIPMENT TO THE STRUCTURE (CONCRETE SLAB OR CONCRETE TOPPING OVER METAL DECK WITH CONCRETE COMPRESSIVE STRENGTH, $f'_c \geq 3000$ PSI, PER NOTE 3E)
6. SHEET METAL SCREWS (SMS) SHALL BE 18-8 STAINLESS STEEL SCREWS, ASTM A593D ($F_y = 40$ KSI $F_u = 80$ KSI).

RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD

1. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2016 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
2. VERIFY THAT THE PROJECT SPECIFIC VALUES OF S_{DS} AND z/h FOUND ON THE PREAPPROVAL DOCUMENTS ARE NOT EXCEEDED.
3. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE POST-INSTALLED ANCHOR ICC-ES EVALUATION REPORT.
4. VERIFY THAT THE POST-INSTALLED ANCHORS COMPLY WITH THE MINIMUM SPACING AND EDGE DISTANCE REQUIREMENTS DEFINED IN THE PREAPPROVAL DOCUMENT.
5. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR $6 \times h_{ef}$ FROM THE UNIT ANCHORS.
6. PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS. VERIFY THE ADEQUACY OF THE STRUCTURE THAT SUPPORTS THE EQUIPMENT FOR THE LOADS IMPOSED ON THEM BY THE EQUIPMENT IN ADDITION TO ALL OTHER LOADS.



SUN Structural Engineering, Inc.
Consulting Structural Engineers



GERMFREE LABORATORIES
HORIZONTAL LAMINAR AIRFLOW WORKSTATION
BZ3, BZ4 AND BZ6

DATE: 09/14/2017

SHEET NO. 1

OF 8

POST-INSTALLED ANCHOR NOTES

1. POST-INSTALLED ANCHORS FOR ATTACHING THE EQUIPMENT TO SUPPORT STRUCTURE HILTI KWIK BOLT TZ STAINLESS STEEL (ICC-ES ESR-1917: REISSUED MAY 2017) AND SHALL COMPLY WITH THE FOLLOWING:

ANCHOR DIAMETER	CONCRETE TYPE	CONCRETE TYPE	MIN. EMBED. h_{ef}	MIN. SPACING	MIN. EDGE DIST.	MIN. CONC. THICKNESS	TORQUE TEST
1/2"	SAND LIGHT WT. OR NORMAL WT. ($f_c' = 3000$ PSI MIN.)	SLAB-ON-GRADE OR ELEVATED SLAB	2"	6"	12"	4"	40 FT-LB

ANCHOR DIAMETER	CONCRETE TYPE	SLAB TYPE	MIN. EMBED. h_{ef}	MIN. SPACING	MIN. EDGE DIST.	MIN. CONC. THICKNESS	TORQUE TEST
1/2"	SAND LIGHT WT. OR NORMAL WT. ($f_c' = 3000$ PSI MIN.)	SLAB-ON-METAL DECK	2"	12"	12"	3 1/4"	40 FT-LB

TESTING OF POST-INSTALLED ANCHORS

2. TESTING OF EXPANSION ANCHOR BOLTS SHALL BE PER 2016 CBC, 1910A.5; TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR AND REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO OSHPD,
- i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, TORQUE TEST AT LEAST 50% OF THE ANCHOR BOLTS.
 - ii) ACCEPTANCE CRITERIA:
 - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE 1/2 TURN OF THE NUT FOR WEDGE TYPE.
 - iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
3. AVOID DAMAGING EXISTING STEEL REINFORCEMENT IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
4. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.

BOLTS THROUGH CONCRETE ON METAL DECK

5. BOLTS SHALL BE TORQUED BY 3/4" TURN OF NUTS AFTER SNUG TIGHT CONDITION IS ACHIEVED: SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING CONNECTED ELEMENTS INTO FIRM CONTACT.
6. HOLES DRILLED IN CONCRETE SLAB SHALL BE 1/16" LARGER THAN BOLT SIZE.
7. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH-BOLTS WITH STEEL-TO-STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.

WELDED STUDS

8. AFTER CABINET HEIGHT IS ADJUSTED TO THE DESIRED LEVEL, WELDED STUD NUT SHALL BE TIGHTENED TO 60 FT-LB TORQUE.

ATTACHMENT MATERIALS

9. ALL BENT STAINLESS STEEL PLATES & ANGLES: ASTM A240, TYPE 304, $F_y = 34$ KSI $F_u = 85$ KSI
10. ALL STRUCTURAL TUBES: ASTM A554, TYPE 304, $F_y = 34$ KSI $F_u = 85$ KSI
11. BOLTS: ASTM F593D, TYPE 304, $F_y = 40$ KSI $F_u = 80$ KSI
12. WELDED STUDS: ASTM F593A, TYPE 304, $F_y = 90$ KSI $F_u = 115$ KSI
13. HARDENED PLAIN WASHER PLATES: ASTM A240 TYPE 304, $F_y = 34$ KSI $F_u = 85$ KSI
- ALL MEMBERS TO BE STAINLESS STEEL



SUN Structural Engineering, Inc.
Consulting Structural Engineers




GERMFREE LABORATORIES
HORIZONTAL LAMINAR AIRFLOW WORKSTATION
BZ3, BZ4 AND BZ6

DATE: 09/14/2017

SHEET NO. 2

OF 8

WORKSTATION DIMENSIONS AND SCHEDULES



WORKSTATION TYPE	WORKSTATION DIMENSIONS AND CENTER OF GRAVITY						WORKSTATION WEIGHT
	A	B	C	X	Y	Z	
BZ3	85.875"	38.75"	37.5"	18.5"	40.25"	13.0"	452 LBS
BZ4	85.875"	50.75"	37.5"	24.5"	40.25"	13.0"	502 LBS
BZ6	85.875"	74.75"	37.5"	36.5"	40.25"	13.0"	648 LBS



SUN Structural Engineering, Inc.
Consulting Structural Engineers



GERMFREE LABORATORIES
HORIZONTAL LAMINAR AIRFLOW WORKSTATION
BZ3, BZ4 AND BZ6

DATE: 09/14/2017

SHEET NO. 3

OF 8

CABINET ATTACHMENT
1/4" S.S. BRACKETS
(ASTM A240 TYPE 304)
WELDED TO CABINET
AND ATTACHED USING
x2 M10 HEX BOLTS
W/ M10 REVIT NUT
FOR MOR INFO., SEE
DETAIL-A

SHOP WELD TO
UNDER-SIDE OF
CABINET, TYP. $\frac{1}{8}$ "
2- PLACES

TYP. $\frac{1}{8}$ "

TYP. ALL
ANGLES $\frac{1}{8}$ " x $\frac{1}{8}$ "

LIFT CYLINDER
ATTACHED TO BRACKET
USING 4x
M5 S.S. MACHINE
SCREWS TYP. 4 PLCS

2-1/2" x 2-1/2" x $\frac{3}{8}$ "
S.S. TUBE FRAME
(ASTM A554
TYPE 304)

1 1/2" x 1 1/2" x $\frac{3}{8}$ " S.S. TUBE
(ASTM 554, TYPE 304);
TYP. OF (4)

CABINET
1/4" S.S. BRACKETS
(ASTM A240 TYPE 304)
WELDED TO CABINET
AND ATTACHED USING
x2 M10 HEX BOLTS
W/ M10 REVIT NUT
FOR MOR INFO., SEE
DETAIL-A

D

L

D'

L'

G'

36"

18"

1 1/8"

4.5"

'H'

2"x2"x 1/8"
S.S. TUBE,
(ASTM A554
TYPE 304)
TYP.

TUBE WELDED
TO SS BASE

3/4" SIZE
HOLE ON TOP
OF CHANNEL
TYP. @ 4
PLACES

1/8" TYP.

1/8"

SHOP WELD TO SIDE OF CABINET, TYP. @ 2- PLACES

1/4" S.S. BRACKET (ASTM A240 TYPE 304)

16 GA. S.S.
SHEET METAL,
ASTM A240
TYPE 304
(SIDES AND
UNDERSIDE
OF CABINET)

$1\frac{1}{2} \times 1\frac{1}{2} \times \frac{3}{8}$ S.S.
TUBE FRAME
(ASTM A554
TYPE 304)

LIFT CYLINDER
BRACKETS ATTACHED
USING 3x M6 S.S.
SCREWS INTO TUBE
FRAME 4 PLCS
FOR MORE INFO.,
SEE DETAIL-D

HYDRAULIC CYLINDER
MAKE-MOVOTEC
MODEL# MLS-00249
CAPACITY 1000 LBS

S.S. LIFT TUBE,
TYP., 4 PLACES

1" DIA. CASTER,
TYP. 4 PLACES
CAPACITY: 300 LBS EA.

3x M6 S.S. SCREWS

18"

5.75"

1" x $\frac{1}{2}$ " BOTTOM
S.S. TUBE FRAME
ASTM A554 TYPE 304
FOR MOR INFO., SEE
DETAIL-C

ROD ATTACHED USING
M10 HEX BOLT, WASHER,
AND LOCK WASHER
TYP. 4 PLCS

SECTION C-C

Technical drawing of a roof curb detail, showing dimensions and materials. The drawing includes a side elevation and a cross-section view.

Dimensions:

- Overall height: 18"
- Top section height: 8 5/8"
- Bottom section height: 8 5/8"
- Top flange width: 1 1/4"
- Bottom flange width: 1 3/4"
- Bottom flange thickness: 3/8"
- Bottom flange width (inner): 2 1/4"
- Bottom flange width (outer): 2 1/4"
- Bottom flange width (inner): 2 1/4"
- Bottom flange width (outer): 2 1/4"
- Bottom flange width (inner): 2 1/4"
- Bottom flange width (outer): 2 1/4"

Materials and Fasteners:

- 4x M5 S.S. SCREWS
- 3x M6 S.S. SCREWS
- 1 1/8" - 1/4" S.S. SHIM PL., LINE UP EDGES W/ BRACKET (ASTM A240 TYPE 304)
- 3- SIDES
- 1- LAYER OF 1/8" S.S. SHEET METAL ON TOP & BOTTOM OF TUBE

Notes:

- BY: Jeffrey Y. Kimoto
- DATE: 09/28/2018
- DETAIL - A

1- 1/8" S.S. SHIM PL., LINE UP EDGES W/ BRACKET (ASTM A240 TYPE 304)

3-SIDES

1 3/4"

2"

3 7/8"

3/4"

3 3/4"

DETAIL-A

3 3/4"

3/4"

1 5/8"

3"

DETAIL-B

1- LAYER OF 3/16" S.S. SHEET METAL ON TOP & BOTTOM OF TUBE

4 1/2"


1"

1/8" 1-6"


2- LAYERS OF 3/16" S.S. SHEET METAL ON SIDES OF TUBE

DETAIL-C

WORKSTATION TYPE	WORKSTATION DIMENSIONS				
	D	E	F	G	H
BZ3	6.125"	23.0"	11.125"	13.0"	37.125"
BZ4	6.125"	35.0"	11.125"	25.0"	49.125"
BZ6	6.125"	59.0"	11.125"	49.0"	73.125"



WORKSTATION SUPPORTS AND ATTACHMENTS FORCES AND DETAILS (FOR $S_{DS} \leq 1.5$)



WORKSTATION TYPE	ANCHORAGE VERTICAL (UPLIFT, T_u) AND HORIZONTAL (SHEAR, V_u) REACTIONS AT INTERFACE WITH CONCRETE SLAB (STRENGTH DESIGN)			
	$z/h=0$, $S_{DS} \leq 1.5$		$z/h \leq 1.0$, $S_{DS} \leq 1.5$	
	T_u /BOLT (LBS)	V_u /BOLT (LBS)	T_u /BOLT (LBS)	V_u /BOLT (LBS)
BZ3	354	100	273	159
BZ4	374	110	288	176
BZ6	459	143	351	228



SUN Structural Engineering, Inc.
Consulting Structural Engineers




GERMFREE LABORATORIES
HORIZONTAL LAMINAR AIRFLOW WORKSTATION
BZ3, BZ4 AND BZ6

DATE: 09/14/2017

SHEET NO. 5

OF 8

WORKSTATION SUPPORTS AND ATTACHMENTS FORCES AND DETAILS (FOR $1.5 < S_{DS} \leq 2.5$)



ANCHORAGE VERTICAL (UPLIFT, T_u) AND HORIZONTAL (SHEAR, V_u) REACTIONS AT INTERFACE WITH CONCRETE SLAB (STRENGTH DESIGN)

WORKSTATION TYPE	$z/h=0, 1.5 < S_{DS} \leq 2.5$		$z/h \leq 1.0, 1.5 < S_{DS} \leq 2.5$	
	T_u /BOLT (LBS)	V_u /BOLT (LBS)	T_u /BOLT (LBS)	V_u /BOLT (LBS)
BZ3	332	83	263	132
BZ4	347	92	275	147
BZ6	719	119	333	190



SUN Structural Engineering, Inc.
Consulting Structural Engineers



GERMFREE LABORATORIES
HORIZONTAL LAMINAR AIRFLOW WORKSTATION
BZ3, BZ4 AND BZ6

DATE: 09/14/2017

SHEET NO. 6

OF 8

4- 1/4"Ø S.S. SELF-TAPPING SCREWS


2"x2"x1/8" S.S. TUBE (ASTM A554 TYPE 304)

1/4"x8" WIDE S.S. BENT PLATE (ASTM A240 TYPE 304) FOR MORE INFO. SEE DETAIL-A

CONCRETE SLAB ON GRADE

(2) 1/2"Ø S.S. HILTI KB-TZ BOLTS EMBEDDED 2" INTO CONCRETE

1" 1" 6"



1

2"x2"x $\frac{1}{8}$ "
S.S. TUBE
(ASTM A554
TYPE 304)

1/8" ∇ AROUND NUT
TO $\frac{1}{2}$ " PLATE

4- $\frac{1}{4}$ " \varnothing S.S.
SELF-TAPPING SCREWS

CONCRETE SLAB
ON METAL DECK

1/4"x2"x2" STL. PL.
(ASTM A240 TYPE 304)

1/8"

3 $\frac{1}{4}$ " MIN.
3" MAX.

MIN. 20 GA.
STEEL W-DECK

(2) 1/2" \varnothing S.S. THRU
BOLTS

S.S. 3"x3"x1/4"
x15" LONG BENT PLATE
(ASTM A240 TYPE 304)
W/ 2- 1/2" \varnothing S.S. HILTI
KB-TZ EMBEDDED 2" INTO
CONCRETE (TOTAL OF 2)

12" MIN. 1" MAX.
OFFSET

4 $\frac{1}{2}$ " MIN.

1/4"x8" WIDE S.S.
BENT PLATE
FOR MORE INFO.
SEE DETAIL-A

HEX NUT TOP & BOTTOM
OF STEEL ANGLE FLANGE, TYP.
IF NUT CAN NOT BE INSTALLED
AT TOP SIDE, PROVIDE TAPPED
HOLE THRU ANGLE FLANGE

\varnothing OF FLUTE \varnothing OF FLUTE

1" MIN.

DATE: 07/28/2016

BY: Jeffery Y. Kiku

OPM 455-13

2nd COMPLIANCE

(4) 1/4" \varnothing HOLES
@ 2" O.C.

1/4"x8" WIDE S.S.
BENT PLATE (ASTM
A240 TYPE 304)

6.75"

1/2"

1 1/2"

3/4"

1 1/2"

1"

TYP. 8"


PLAN VIEW

A STEEL BENT PLATE DETAIL




OF 8


WORKSTATION ANCHORAGE DETAILS (FOR $1.5 < S_{DS} \leq 2.5$)



1 ANCHORAGE TO CONCRETE SLAB-ON-GRADE WITH 4" MIN. THICKNESS



2 ANCHORAGE TO CONCRETE SLAB ON METAL DECK



A STEEL BENT PLATE DETAIL



SUN Structural Engineering, Inc.
Consulting Structural Engineers



GERMFREE LABORATORIES
HORIZONTAL LAMINAR AIRFLOW WORKSTATION
BZ3, BZ4 AND BZ6

DATE: 09/14/2017

SHEET NO. 8

OF 8