

OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL	OFFICE USE ONLY									
OF MANUFACTURER'S CERTIFICATION (OPM)	APPLICATION #: OPM-0310-13									
OSHPD Preapproval of Manufacturer's Certification (OPM)										
Type: ☐ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number:										
Manufacturer Information										
Manufacturer: Samsung										
Manufacturer's Technical Representative: Ninad Gujar										
Mailing Address: 14 Electronics Ave., Danvers, MA. 01923										
Telephone: On File Email: Don File										
Product Information	ON									
Product Name: X-Ray Universal Arm System										
Product Type: Other Mechanical and Electrical Components 13	Z.									
Product Model Number: X-Ray Universal Arm System	lin H									
General Description: Digital Radiography system component										
DATE: 08/24/2016										
Ţ.	\(\frac{\chi}{\chi}\)									
Applicant Information	\$ °									
Applicant Company Name: EASE Co.	COA									
Contact Person: Jonathan Roberson, S.E.										
Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709										
Telephone: (909) 606-7622 Email: J.Rob I hereby agree to reimburse the Office of Statewide Health Faccordance with the California Administrative Code, 2013.	erson@EASECo.com Planning and Development review fees in									
Signature of Applicant:	Date: 2/9/16									
Title: Principal Engineer Company Name: EASE	Co.									

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

os Dpd

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY OSH-FD-700 (REV 1/24/13)

Page 1 of 2



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations								
Company Name: EASE Co.								
Name: Jonathan Roberson, S.E. California License Number: S4197								
Mailing Address: 5877 Pine Ave. Suite 210, Chino Hills, CA. 91709								
Telephone: 909-606-7622 Email: J.Roberson@EASECo.com								
OSHPD Special Seismic Certification Preapproval (OSP)								
 □ Special Seismic Certification is preapproved under OSP-(Separate application for OSP is required) □ Special Seismic Certification is not preapproved 								
Certification Method(s)								
☐ Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-10 ☐ Other* (Please Specify):								
*Use of test criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing. Analysis								
OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2013 ONLY								
Signature:								

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"

Page 2 of 2



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

5877 Pine Ave, Ste. 210 Chino Hills, CA. 91709 Phn: (909) 606-7622

Office of Statewide Health Planning and Development
PREAPPROVAL OF MANUFACTURER'S CERTIFICATION
OPM-0310-13

THIS PREAPPROVAL CONFORMS TO THE 2013 CALIFORNIA BUILDING CODE

MANUFACTURER:

SAMSUNG

Sheet: 1 of 7

EQUIPMENT NAME:

GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM

Date: 5/31/16

GENERAL NOTES

- 1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2013 CBC. THE DEMANDS (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2013 CBC
- 2. THIS DOCUMENT MAY ONLY BE USED WITH THE EXPRESS WRITTEN CONSENT OF THE MANUFACTURER LISTED ABOVE FOR THE SPECIFIC PROJECT SITE AND INSTALLATION LOCATION. THIS DOCUMENT IS INVALID WITHOUT SUCH CONSENT.
- 3. THIS PREAPPROVAL CONFORMS TO THE 2013 CALIFORNIA BUILDING CODE WHERE SDS IS NOT GREATER THAN 2.00.
- 4. FORCES PER ASCE 7-10 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE SDS ≥ 2.00, ap = 1.0, lp = 1.5, Rp = 1.5, z/h = 0 AT CONCRETE SLAB & z/h < 1 AT CONCRETE SLAB ON METAL DECK. SEE FOLLOWING SHEETS FOR Ω₀
- 5. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
- ALL DESIGN FORCES SHOWN ON THE DRAWINGS ARE FACTORED LOADS THAT SHALL BE USED FOR STRENGTH DESIGN.
- 7. CONCRETE SLAB ON METAL DECK DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION IN THE BUILDING. (i.e. z/h < 1)
- 8. CONCRETE SLAB ON GRADE DETAIL VALID FOR DEMANDS SHOWN AT ANY ELEVATION BELOW GRADE. (i.e. z/h = 0)

9. RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD OF THE BUILDING

- A. PROVIDE SUPPORTING STRUCTURE TO SUPPORT WEIGHTS AND FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
- B. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2013 CBC AND WITH THE DETAILS, MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PREAPPROVAL DOCUMENTS.
- C. VERIFY THAT PROJECT SPECIFIC VALUES OF SDS & z/h RESULT IN SEISMIC FORCES (Eh, Ev) THAT DO NOT EXCEED THE VALUES ON THE DETAILS.
- D. VERIFY THAT THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ANCHORED MEETS THE REQUIREMENTS OF THE APPLICABLE ICC ESR.
- E. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS (SEE TYPICAL DETAIL ON SHEET 2).
- F. VERIFY THAT ALL NEW OR EXISTING ANCHORS ARE AN ADEQUATE DISTANCE FROM THE UNIT ATTACHMENTS AND CHECK FOR INTERACTION WHERE OTHER ANCHORS ARE WITHIN 18" OR 6hef FROM THIS UNIT'S ANCHORS.



EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

11-1528

www.EquipmentAnchorage.com

SAMSUNG

DES. J. ROBERSON

2

GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM

DATE 5/31/16

JOB NO.


_f 7 sheets

10. EXPANSION ANCHORS:

A. ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

Anchor Diameter	Concrete Type	Min. f'c (psi)	Anchor Type	ICC Report No.	Min. Embed.	Min. Spacing	Min. Edge Dist.	Min. Conc. Thickness	Torque Test	Direct Tension
3/8"	Sand Light Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	2"	N/A	N/A	See Sheet 7 of 7	25 FT-LB	1173 lb
5/8"	Normal Weight	3000	Hilti Kwik Bolt TZ	ESR-1917	4"	6"	24"	6"	60 FT-LB	3405 lb

- B. THIS PREAPPROVAL ALLOWS FOR UP TO A MAXIMUM OF 2 ADJACENT CONCRETE SLAB EDGES, 24" AWAY MINIMUM (i.e. CORNER). SEE ADJACENT DETAIL FOR ADDITIONAL MINIMUM ALLOWABLE CONCRETE EDGE DISTANCES.
- C. TESTING OF EXPANSION ANCHORS PER 2013 CBC, 1913A.7:
 TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO OSHPD
 - (i) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION,
 DIRECT PULL TENSION TEST OR TORQUE TEST AT LEAST 50% OF
 THE ANCHORS.
 - (ii) ACCEPTANCE CRITERIA:
 - DIRECT TENSION TEST: THE ANCHOR SHOULD HAVE NO 2 0 1 6
 OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY
 TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER
 BECOMES LOOSE.
 - TORQUE TEST: THE APPLICABLE TORQUE MUST BE ACHIEVED WITHIN THE FOLLOWING LIMITS: WEDGE TYPE: 1/2 TURN OF THE NUT
 - (iii) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.
- D. AVOID DAMAGING EXISTING STEEL REINFORCING IN CONCRETE SLAB WHEN INSTALLING CONCRETE EXPANSION ANCHORS.
- E. PROVIDE FOR FULL THREAD ENGAGEMENT OF NUT & WASHER.
- 11. BOLTS THROUGH CONCRETE ON METAL DECK
 - A. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE SNUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
 - B. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE.
 - C. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



TYPICAL CONCRETE EDGE DETAIL

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

DES. J. ROBERSON

www.EquipmentAnchorage.com

SAMSUNG

A/GU60A-65

GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM


DATE 5/31/16

SHEET 3

SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB



NOTES:

1. FORCES ARE DETERMINED PER 2013 CALIFORNIA BUILDING CODE AND ASCE 7-10

STRENGTH DESIGN IS USED. (SDS = 2.00, a_p = 1.0, I_p = 1.5, R_p = 1.5, Ω_o = 1.5, z/h = 0)

HORIZONTAL FORCE (En) = 0.90 Wp UILDING

HORIZONTAL FORCE (Emh) = 1.35 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (Ev) = 0.40 Wp

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER

LOADS THAT MAY BE PRESENT.

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

5/31/16

www.EquipmentAnchorage.com

SAMSUNG


DES. J. ROBERSON


SHEET 4

GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM јов но. 11**-1528**

DATE

OF 7 SHEETS





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

www.EquipmentAnchorage.com

SAMSUNG

JOB NO. 11-1528

DES. J. ROBERSON

5

SHEET


GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM

DATE 5/31/16

7 SHEETS

SEISMIC SUPPORTS & ATTACHMENTS

CONCRETE SLAB ON METAL DECK



NOTES:

1. FORCES ARE DETERMINED PER 2013 CALIFORNIA BUILDING CODE AND ASCE 7-10.

STRENGTH DESIGN IS USED. (SDS = 2.00, ap = 1.0 lp = 1.5 Rp = 1.5, Ω_0 = 1.5, z/h \leq 1)

HORIZONTAL FORCE (En) = 2.40 Wp

HORIZONTAL FORCE (Emh) = 3.60 Wp (FOR CONCRETE ANCHORAGE)

VERTICAL FORCE (E_V) = 0.40 W_p

2. CENTER OF GRAVITY (C.G.) AND WEIGHT ARE THE GOVERNING PARAMETERS FOR DESIGN. THIS PREAPPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.

3. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE

DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN IN COMBINATION WITH ALL OTHER

LOADS THAT MAY BE PRESENT.

EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

11-1528

www.EquipmentAnchorage.com

SAMSUNG

DES. J. ROBERSON

6


SHEET


GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM

DATE 5/31/16

JOB NO.

OF 7 SHEETS





EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

11-1528

www.EquipmentAnchorage.com

SAMSUNG

DES. J. ROBERSON


SHEET 7

GU60A/GU60A-65 X-RAY UNIVERSAL ARM SYSTEM

DATE 5/31/16

JOB NO.

OF 7 SHEETS




MIN STEEL DECK REQUIREMENTS AND STRUT DETAIL

BY: WITHROUGH STRUTEFLANGE.

CANNOT BE PROVIDED AT TOP SIDE

OF STRUT, PROVIDE TAPPED HOLE



L3.5 X 3.5 X 5/16" X 1'-5" MIN

(A36) AT EACH ANCHOR