


ATTACHMENT 2: TEST SPECIMEN SUMMARY

ATTACHMENT PAGE | 6 OF 6

UUT-2002-8 Con Box - Type B w/ Optional Display Console Kit						
MANUFACTURER: Canon Medical Systems Corporation						
IDENTIFICATION: Model No.: KCCN-020C / 1B S/N.: 1BA1977001						
DESCRIPTION: Component of the Aquilion Lightning 40/80/160 CT System Includes <u>Optional Display Console Kit</u> Model No.: CGS-72B / 1B; Serial No.: 1BA1977001 Requires Seismic Kit P/N: ANC-CA-002COX						
MOUNTING: <u>Rigid Base (Floor) Mounted:</u> (4) – 3/8" dia. SAE J429 Grade 5 Hex Head Bolts to aluminum interface plate.						
PROPERTIES:						
DIMENSIONS (in.)			Weight (lb.)	LOWEST RESONANT FREQUENCY (Hz.)		
Width	Depth	Height		X-Axis	Y-Axis	Z-Axis
23.2	37.6	53.7	652	> 50	19.05	12.36
ICC-ES AC156 SHAKE TABLE TEST PARAMETERS						CODE: 2019 CBC
S _{DS} (G)	z/h	I _P	A _{FLX-H} (G)	A _{RIG-H} (G)	A _{FLX-V} (G)	A _{RIG-V} (G)
2.0 2.5	1 0	1.5	3.20	2.40	1.68	0.68
Unit satisfied AC156 requirements for structural integrity and manufacturer requirements for functionality after AC156 test.						